

# 10. (October) 2019

---

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Anja

FRIDAY, 11/10/2019

---

## Anja

**Cytotoxicity (LDH) Assay** 48h after transfection (according to Cytotoxicity (LDH) Assay protocol)

- 3 technical replicates were taken from each well
- Triton-X-100-PI-DMEM mix volume adjusted to volume in 24-well plates (scale-up)
- done with background controls including the only full DMEM and Zeocin+full DMEM mix wells

## Joshi

**HEK 293T harvesting for long-term cell-monitoring experiment** from 24-well plate for qPCR analysis (48h after transfection)

- 600µl harvested and centrifuged (according to VLPs and exosome harvesting)
- 175µl of Total Exosome Isolation reagents added to 350µl supernatant (top part after centrifugation)
- incubation in the fridge ON

**EXO HiBiT Assay for HEK293T long-term cell-monitoring experiment** from 24-well plate (48h after transfection, according to VLPs

/ Exosomes harvesting & HiBiT assay from 24-well plates)

- 40µl used for supernatant lysed (SNlysed)
- 10µl used for supernatant unlysed (SNunlysed)

**Cell culture: Cytotoxicity assay medium change** on 72h plate (samples including or not including Zeocin™ antibiotic in the media)

- long-term cell-monitoring plate: 600µl harvested and saved for HiBiT and qPCR analysis, 600µl fresh media with/out antibiotic added
- cytotoxicity assay 48 & 72h plates: 600µl removed, 600µl fresh media with/out antibiotic added

**Cell culture: Splitting and seeding**

- splitted 1:6
- seeded 1x 24-well plate with 150.000 cells/well into 800µl full DMEM