

07. (July) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Theresa Keil

MONDAY, 8/7/2019

DNA digestion of samples from 04.07.2019

Table3: Protocol for VLP samples ...

	Reagent	Volume (µl)
1	RNA	0.6
2	Reaction Buffer	10
3	DNase	1
4	H2O	88.4

Table4: Protocol for DNase treatm...

	Reagent	Volume (µl)
1	RNA	11.6
2	Reaction Buffer	10
3	DNase	1
4	H2O	77.4

- samples were incubated for 10min at 37°C & then splitted into half
- one half EDTA treatment: Add 1 µl of 0.5 M EDTA (to a final concentration of 5 mM), Heat inactivate at 75°C for 10 min
other half direct Phenol-/Chloroform treatment