

08. (August) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Theresa Keil

FRIDAY, 9/8/2019

exosome precipitation: continued with protocol from invitrogen form step 4

Prepared samples for SDS-PAGE and WB for monday

mixed 25 µl of 1a, 2a, 3a, 4a, 5aa, with 25 µl laemmli (2x), incubated 15 min at 95 °C

DNase treatment: 10 µl DNase I + 1/10 volume buffer, 1 h , 37 °C, 200 rpm

Pooled the triplicates and made duplicated of 170 µl

Lysis added 1 volume 170 µl 1% triton in 1x PBS

Spiking of E1a-5a with 1 µl XPA (5 ng/µl)

RNA-Extraction with Trizol like protocol used 400 µl of Trizol

Table5				
	A	B	C	D
1	sample	concentration ng/µl	260/280	260/230
2	E1a	102.8	1.56	0.34
3	E2a	371.1	1.74	0.66
4	E3a	433.7	1.93	0.81
5	E4a	724.8	1.47	0.37
6	E5a	450.2	1.87	0.83
7	E1b	253.4	1.68	0.47
8	E2b	351.3	1.80	0.63
9	E3b	459.0	1.94	0.92
10	E4b	440.0	1.83	0.81
11	E5b	228.4	1.56	0.48

1 % agarose gel was performed with B16 and B17 from 08.08 qPCR
the bands of fluc were cut off and gel digest was performed