

Protocol for cell culture and subculture

- 1. Culture in a humidified incubator at 37 °C, 5% CO₂.**
- 2. Check the color medium, healthy usually leaves medium slightly orange; observe the cells under phase microscope.**
- 3. If the cells are not very confluent, remove the growth medium and replace with fresh complete medium.**
- 4. If the cells are at least 70 %~ 80% confluence, subculture them.**
- 5. Subculturing cells (60 mm×15 mm cell culture plate):**
 - (1) Aspirate the growth medium from the cells.and wash the cells with PBS.
 - (2) Aspirate the PBS, then add 600μL trypsin EDTA solution.
 - (3) Incubation for several minutes at 37 °C and check the culture with inverted microscope to be sure that the cells are rounded up which indicates they are detached from the surface.
 - (4) Add 2 mL fresh complete medium,then pipet all the adherent cells into cell suspension and pipet up and down.
 - (5) Transfer the cell suspension into a 15 mL centrifuge tube and centrifuge at 1200 rpm for 5 minutes.
 - (6) Aspirate the filtrate ,then add 1mL PBS and resuspend the cell.
 - (7) Place 5 ml of fresh medium to each new culture .Add equal volume of cell suspension to each of the fresh plates that have been appropriately labelled.
 - (8) Culture in a humidified incubator at 37 °C, 5% CO₂.