


# HiFi cloning

WOENSDAG 11-9-2019

Hifi Assembly protocol

We linearized the vector AraC3 by PCR trying temperatures :65/68 and 72 with following result



At 72 °C PCR results in strong band of the right size and very few unspecific product. We PCR purified and ran the sample again with a positive result and a concentration of 38.8ng/μl.

We then proceeded and prepared the hifi assembly for bb1.2 and 3 as well as for a negative control (no enzyme) and the positive control from the kit.

## AraC3+bb1.2

Optimal ng amounts according to NEB are:

Max total 0,2pmol for 2-3 inserts

25,10 fmol Vector (3224bp) = 50ng = (38.8ng/μl) = 1,29μl

50.2 fmol bb1.2\_Amp (1380) = 46ng = (147ng/μl) = 0,32μl

Water: 8,39μl

Total: 0,0753 pmol

## AraC3+bb3

Optimal ng amounts according to NEB are:

Max total 0,2pmol for 2-3 inserts

25,10 fmol Vector (3224bp) = 50ng = (38.8ng/μl) = 1,29μl

50.2 fmol bb3\_Amp (1335) = 41.41ng = (76ng/μl) = 0,55μl

Water: 8,16

Total: 0,0753 pmol

Hifi mix was transformed into competent cells and plated overnight, starting at 10 PM

ZATERDAG 14-9-2019

Redoing the ligation of AraC3 with BB3

This time I reduced the amount of template for vector linearization to 0,1ng to reduce background transformation

PCR product will be shown here:

after purification I will follow the same recipe for ligation :

## AraC3+bb3

Optimal ng amounts according to NEB are:

Max total 0,2pmol for 2-3 inserts

25,10 fmol Vector (3224bp)= 50ng = (27ng/μl)=1,85 μl

50.2 fmol bb3\_Amp (1335) = 41.41ng =(76ng/μl) = 0,55μl

Water:

Total: 0,0753 pmol

10μl Hifi Master mix

Negative control:

reaction mix with water instead of enzyme

Of those ligations i will transform 1/2/3/5 μl both for experiment and negative controll, cause i think the amount of DNA used for transformation will also influence the background. This way i kind find the amoutn of DNA for transformation that does lead to thelowest background while still resulting in colonies

Negative con

Table39

	A	B	C	D	E
1		Reaction		negative	
2	Insert	0.55		0.55	
3	Vector	1.85		1.85	
4	Water	7.6		17.6	
5	Mastermix	10μl		0μl	
6					
7					
8					

Transformation

	A	B	C	D	E	F
1		Volume				positive
2	KCM	2	2	2	2	2
3	Water	7	6	5	3	3
4	DNA	1	2	3	5	5(1ng/μl)araC3

#### MAANDAG 16-9-2019

AraC3\_BB1.2

25,10 fmol Vector (3224bp)= 50ng = (27ng/μl)=1,85 μl

50.2 fmol bb1.2\_Amp (1380) = 42.80ng(+5 %) =(147ng/μl) = 0,31μl

Water: 7,84μl

Total: 0,0753 pmol

10μl Hifi Master mix

**Arac3+OEbb1+OEbb2**

25,10 fmol Vector (3224bp)= 50ng = (27ng/μl)=1,85 μl  
 50,20 fmol OEbb1 (1261bp)= 39.11(+2%)ng =( 142,85 ng/ul)= 0,28μl  
 50,20 fmol OEbb2 (2019) = 62,62(+2%) ng = (169 ng/μl) =0,38μl  
 Total : 0,1256 pmol  
 Water: 7,49  
 10μl HIFI mix

negative:

25,10 fmol Vector (3224bp)= 50ng = (27ng/μl)=1,85 μl  
 50.2 fmol bb1.2\_Amp (1380) = 42.80ng(+5 %) =(147ng/μl) = 0,31μl  
 Water: 17,84μl  
 Total: 0,0753 pmol  
 0 μl Hifi Master mix

**VRIJDAG 20-9-2019**

Transformation of AraC3\_1.2 aswell as AraC3\_1+2 led to low background.  
 colonies were scr

**Arac3+OE1.2+OEbb3**

25,10 fmol Vector (3224bp)= 50ng =(27ng/μl)= 1,85μl  
 50,20 =30%=(65)fmol OE1.2 (1351bp)= 54.72ng =( 95ng/ul)= 0,576μl  
 50,20 =50%= (100,4)fmol OEbb2 (1311) = 81,34 ng = (44ng/μl) =1,84μl  
 Total 190,4 fmol  
 Water 5,73  
 MM 10μl  
 Negative control:  
 25,10 fmol Vector (3224bp)= 50ng =(27ng/μl)= 1,85μl  
 50,20 =30%=(65)fmol OE1.2 (1351bp)= 54.72ng =( 95ng/ul)= 0,576μl  
 50,20 =50%= (100,4)fmol OEbb2 (1311) = 81,34 ng = (44ng/μl) =1,84μl  
 Total 190,4 fmol  
 Water 15,73  
 MM 0μl

Since the gel shows unspecific bands i guesed the real concentration to be 30% and 50% lower for OE1.2 and OE3 respectively, and therefore calculated 30 and 50% higher amount of moles. I am still below the maximum of 200fmol for hifi cloning. The controll as allways will be the reaction mix without enzyme. A positive transofmrtaion controll will be AraC3

Sequencing of

**AraC3\_1+2: Primers:**

VR  
 Seq2R  
 Seq1F  
 Seq1R  
 VF2

12μl per reaction = 60μl at 40-100ng/μl

AraC3\_1.2: Primers

VR

Seq1R

VF2

12µl per reaction: 36µl at 40-100ng/µl

Primers to send:

Seq2R\_HIFI: 10µM 50µl

Seq1F\_HIFI: 10µM 50µl

Seq1R\_HIFI :

10µM 50µl

---

## MAANDAG 23-9-2019

### Arac3+OE1.2+OEbb3

Colony PCR

Primers:

Seq3R

Seq1F

Amplicon : 1223 bp

Elongation Time 1min 30 sec

Annealing temp: 95,3°C

Sequencing results for **Arac3+OEbb1+OEbb2** came back #18 has almost no mutations and no stopcodons, which should render it functional. However some gaps need to be sequenced again to make sure.

I made two liquid cultures, 15 ml TSB 15µl Chloramphenicol, 15µl Ecoli containing SP, and added 300µl of a 50x 10%(20%?) Arabinose solution to one of the flasks.

Sequencing of **Arac3\_1.2** showed a mutation in the C protein, creating a stopcodon and therefore rendering it nonfunctional, we have to screen more colonies, since the cloning worked other successful clones should have plasmids without fatal mutations.

---

## DINSDAG 24-9-2019

Use SeqP1R to sequence the missing gap of RFP in **Arac3+OEbb1+OEbb2 (Done)**

### AraC3\_1.2


We need to screen the AraC3\_1.2 colonies for more positives since #6 has a mutation.

Screening will be done with P1F and VR which should result in a Amplicon of 770 bp(1min). #6 can be used as positive control.

We need to screen more Model plasmid since none of the 14 screened colonies has the right insert.


We took colonies from the AraC\_1.2 Plate 6-22 and MP plate 15-27. AraC\_1.2 and MP PCR was done in two machines, one at 1min elongation and the other at 2min 10 sec elongation time.





Difficulty is that we don't have a positive control for the PCR reaction, therefore we use P1F and VR, like we do for 1.2 screening. This amplifies the region between BB1 and VR, containing BB3. Positive colonies should have an amplicon of 2028 bp (2 Min 10 sec).

I also tested the ligation of the MP and AraC3\_BB3, they both seemed to have worked, due to the disappearance of the insert bands.




#### DONDERDAG 26-9-2019

To do list:

Screen AraC3\_bb3 plates for positives.

Send positive AraC3\_BB1.2 for sequencing-



Do PCR for OE1.2 and OE3

Redo the ligation for MP

Take SP and delete bb1 and bb2 to create AraC3\_bb1 and AraC3\_bb2 respectively.

### PCR for amplified1.2 and amplified BB3


TM 65°C

Elongation time 45sec

Max\_primer\_BB1/BB1.2

Tobi\_BB3

Q5 Biobrick amplification reaction mix							
	A	B	C	D	E	F	G
1	Recepie		Biobrick	negative control	STEP	TEMP	TIME
2	Q5 High-Fidelity 2X Master Mix	12.5 µl	12,5µl	12,5µl	Initial Denaturation	98°C	30s
3	10 µM Forward Primer	1.25 µl			25 Cycles	98°C 66°C 72°C	10 seconds 30 seconds 30 seconds/kb= 1 min
4	10 µM Reverse Primer	1.25 µl	MM 1,25 ml	MM 1,25µl	Final Extension	72°C	2 minutes
5	Template DNA	10ng/µl	1µl	0µl	Hold	4°C	
6	Nuclease-Free Water	to 25 µl	10,25µl	11,25µl			



One of each BB was used for PCR purification and subsequent OE PCR

### PCR for OE1.2 and OE3

OE\_1.2\_1.2+3 67°C


OE\_3\_1.2+3 69°C

Elongation time 45 sec


One of each was used in Gelpurification(hope it works this time, i am using washbuffer and elution buffer from the PCR kit, since they seemd to work. Wahbuffer of the PG kit didnt contain any enthanol(didnt smell like it))

This does not work, for some reason OE1.2 is shorter then OE 3, I think this is because the PCR purified BB1.2 and 3 contain unspecific bands leading to some kind of shit.

I will gel purify the second BB1.2 and 3 tomorrow. The gel purified BB1.2 and 3 will then be used as template for OE PCR. If this doesnt work fuck it, why am I doing this anyway.



 Gel purified OE1.2 and 3



## ZATERDAG 28-9-2019

### Arac3+OE1.2+OEbb3

25,10 fmol Vector (3224bp)= 50ng =(27ng/μl)= 1,85μl

50,20 fmol OE1.2 (1351bp)= 41.91ng =( 17,4ng/ul)= 2,41μl

50,20 fmol OEbb2 (1311) = 40,67 ng = (12,3ng/μl) =3,31μl

Total 125,5 fmol

Water 2,43

MM 10μl

Negative control:

25,10 fmol Vector (3224bp)= 50ng =(27ng/μl)= 1,85μl

50,20 fmol OE1.2 (1351bp)= 41.91ng =( 17,4ng/ul)= 2,41μl

50,20 fmol OEbb2 (1311) = 40,67 ng = (12,3ng/μl) =3,31μl

Total 125,5 fmol

Water 12,43

MM 0μl

Screening of the MP

PRimers:

PIF

VR

Amplicon: 2028

Anealing temp 59.3

## MAANDAG 30-9-2019

Sequencing

AraC3\_BB1.2\_20 3 sequences =36μl

AraC3\_BB1.2\_21 3 sequences=36μl

AraC3\_BB3\_4 3 sequences=36μl

AraC3\_BB3\_7 3 sequences=36μl

MP\_8 4 sequences=48μl

MP\_13 4 sequences=48μl

New primer to send: P3R 50μl of 10μM solution

All the other primers are at microsynth

DONDERDAG 3-10-2019

---

Plan: screen colonies BB3 - send to sequencing

THURSDAY - site directed mutagenesis - PCR 3 temperatures -

### Deltetion Primers

<b>Delete BB1 R</b>	<b>GTTGAACAGTACGAACGTGCCGAGG</b>	Snapgene: 64°
<b>Delete BB1F_II</b>	<b>AGCTAGCACTGTACCTAGGACTGAGCTAGC</b>	NEB: 72° Snapgene: 65°
Delete BB2 R_II	GCAAATAATCAATGTTGGCCGGCTTGACG	Snapgene: 65°
Delete BB2_F_II	ACTACGCTGACGCTTCTTAATAACTGTAACAGAGC	NEB: 72° Snapgene: 64°

Following the protocol from NEB, but using our regular q5 mastermix, since we only have 10 reactions which will be saved until the construct with the correct RBS arrives.

<https://international.neb.com/protocols/2013/01/26/q5-site-directed-mutagenesis-kit-quick-protocol-e0554>

Temperatures tried: 65°, 68°, 72°

**Cycling conditions:** Q5 regular protocol, but elongation time 2:50 min

Stela used the pcr products from 72°C and went on with the ligation and transformation protocol.

Deltion was succesfull, now we need to screen the colonies for correct clones.

Screening:

### Machine 1: elongation time: 1min

AraC3\_BB1

primers:

seq1F

VR

Amplicon: 782bp

Anealing temp: 59.3

### Machine 2: elongation time 2:30

AraC3\_BB2:

primers:

Seq2R

VF2

Amplicon: 2329

Anealing temp 59.3

Screening of the MP


PRimers:

Seq 1F

VR

Amplicon: 2028

Anealing temp 59.3



Three colonies of each were taken.

#### MAANDAG 7-10-2019

Sequencing of colonies

AraC3\_BB1:

Primers:

VR

Seq1R-HIFI

AraC3\_BB2

Primers:

VR

Seq2R-HIFI

MP

Primers:

Seq3F-II need to send

Seq1R-II

Ara3\_BB1.2

Primers

VR

Seq1R

Ara3\_BB3

Primers

VR

P3R

AraC3:

Primers

VF2

#### DONDERDAG 10-10-2019

Adding RBS site to BB1 and BB2

BB1:

Primers:

Annealing temp 72°C

Elongation time: 2:30

BB2:

Primers:

Annealing temp 72°C

Elongation time: 2:50

Put both in the same machine at 72°C and elongation time 3 minutes

Elongation time with Q5 is 30sec/KB



PCR only worked with BB2, so I proceeded with BB2. cloning was done, control was lysed bacteria, positive was arac3 plasmid

## VRIJDAG 11-10-2019

We send some colonies of BB2 for sequencing bei ENS

Today stella did mutagenesis of BB1, with successful PCR. She did cloning and Saturday we can take the colonies.

## WOENSDAG 16-10-2019

thermo fisher construct is here

I dilute in 50µl of water to a concentration of 100ng/µl

I will use: Ampicilin plate

2µl DNA

2µl KCM

6µl DNA

Positive control: Chloramphenicol plate

AraC3:

2µl DNA

2µlKCM  
6µl water

negative control Ampicilin plate  
2µlKCM  
8µl water

Working timeline:  
Do Colony PCR with  
Seq1F  
Seq2R  
at 59.3 °C  
Amplicon 1451bp ( 1min 45 sec)

grow over night:  
the positive selected colonies  
RFP control  
ask for gfp control  
at 8 am dilute to 0.5 OD and induce with 0.2 Arabinose until 2 PM  
Meassure at FACS