

08. (August) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

TUESDAY, 6/8/2019

Thess

cell culture: transfection for VLP Purification

- 550 µL medium were exchanged before transfection
- all conditions were seeded as duplicates

transfection scheme Purification assay 30/07/2019

	condition	V8	V27	V28
1	W1	250 ng	250 ng	-
2	W2	167 ng	333 ng	-
3	W3	333 ng	167 ng	-
4	W4	250 ng	-	250 ng
5	W5	167 ng	-	333 ng
6	W6	333 ng	-	167 ng

- Transfection:

Transfectionmix 24-well plate 30/0...

	A	B
1	DNA per well	500 ng
2	P3000 Reagent per well	1 µL
3	Lipofectamine 3000 reagent per well	0.75 µL
4	OptiMEM per well	2 x 25 µL

- Transfect cells according to the following table. Use the indicated volume of DNA and P3000™ Reagent with each of the two volumes of Lipofectamine™ 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.

Timeline	Steps	Procedure Details (Two Reaction Optimization)			
Day 0	1. Seed cells to be 70-80% confluent at transfection	Control well	24-well	6-well	
	Adherent cells	4.4 x 10 ⁵	0.5-2 x 10 ⁵	0.25-1 x 10 ⁶	
Day 1	2. Dilute Lipofectamine™ 3000 Reagent in Opti-MEM® Medium (2 tubes) - Mix well	Opti-MEM® Medium	5 µL x 2	25 µL x 2	125 µL x 2
	Lipofectamine™ 3000 Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL	
Day 1	3. Prepare master mix of DNA by diluting DNA in Opti-MEM® Medium, then add P3000® Reagent - Mix well	Opti-MEM® Medium	10 µL	50 µL	250 µL
	DNA (0.5-5 µg/µL)	0.2 µg	1 µg	5 µg	
Day 1	4. Add Diluted DNA to each tube of Diluted Lipofectamine™ 3000 Reagent (1:1 ratio)	P3000® Reagent (2 µL/µg DNA)	0.4 µL	2 µL	10 µL
	Diluted DNA (with P3000® Reagent)	5 µL	25 µL	125 µL	
Day 1	5. Incubate	Diluted Lipofectamine™ 3000 Reagent	5 µL	25 µL	125 µL
	Incubate for 10-15 minutes at room temperature.	Opti-MEM® (per well)	5 µL/well	25 µL/well	125 µL/well
Day 2-4	6. Add DNA-lipid complex to cells	DNA-lipid complex	10 µL	50 µL	250 µL
	DNA amount	100 ng	500 ng	2500 ng	
Day 2-4	7. Visualize/analyze transfected cells	P3000® Reagent	0.2 µL	1 µL	5 µL
	Lipofectamine™ 3000 Reagent used	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL	

Anja

cell culture: transfection for hArc

- 550 µL medium were exchanged before transfection
- all conditions were seeded as duplicates

transfection scheme hArc assay 06/08/2019			
	V8	V25	V30
1	-	-	500 ng
2	500 ng	-	-
3	-	500 ng	-

- Transfection:

Transfectionmix 24-well plate 06/0...		
	A	B
1	DNA per well	500 ng
2	P3000 Reagent per well	1 µL
3	Lipofectamine 3000 reagent per well	0.75 µL
4	OptiMEM per well	2 x 25 µL

- Transfect cells according to the following table. Use the indicated volume of DNA and P3000™ Reagent with each of the two volumes of Lipofectamine™ 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.

Timeline	Steps	Procedure Details (Two Reaction Optimization)			
Day 0	1 Seed cells to be 70–80% confluent at transfection	Component	6-well	24-well	96-well
		Adherent cells	4×10^5	0.5×10^5	0.25×10^5
		Opti-MEM [®] Medium	5 μ L \times 2	25 μ L \times 2	125 μ L \times 2
Day 0–2	2 Dilute Lipofectamine [®] 3000 Reagent in Opti-MEM [®] Medium (2 tubes)– Mix well	Lipofectamine [®] 3000 Reagent	0.15 and 0.3 μ L	0.75 and 1.5 μ L	3.75 and 7.5 μ L
		Opti-MEM [®] Medium	10 μ L	50 μ L	250 μ L
		DNA (0.5–5 μ g/ μ L)	0.2 μ g	1 μ g	5 μ g
Day 1	3 Prepare master mix of DNA by diluting DNA in Opti-MEM [®] Medium, then add P3000 [®] Reagent– Mix well	P3000 [®] Reagent (2 μ L/ μ g DNA)	0.4 μ L	2 μ L	10 μ L
		Diluted DNA (with P3000 [®] Reagent)	5 μ L	25 μ L	125 μ L
		Diluted Lipofectamine [®] 3000 Reagent	5 μ L	25 μ L	125 μ L
Day 1–4	4 Add Diluted DNA to each tube of Diluted Lipofectamine [®] 3000 Reagent (1:1 ratio)				
	5 Incubate	Incubate for 10–15 minutes at room temperature.			
	6 Add DNA-lipid complex to cells	Component (per well)	6-well	24-well	96-well
	DNA-lipid complex	10 μ L	50 μ L	250 μ L	
	DNA amount	100 ng	500 ng	2500 ng	
	P3000 [®] Reagent	0.2 μ L	1 μ L	5 μ L	
	Lipofectamine [®] 3000 Reagent used	0.15 and 0.3 μ L	0.75 and 1.5 μ L	3.75 and 7.5 μ L	
Day 2–4	7 Visualize/analyze transfected cells	Incubate cells for 2–4 days at 37°C. Then, analyze transfected cells.			