

#1 Preperation of solutions and LB-media

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 17/6/2019

Experimenter: Vamsi, Karthik and Johanna

Aim: To prepare solutions for the LB-media

Protocol: Lab Manual Synthetic Biology, LB medium page 98

600 ml of LB medium was prepared with:

NaCl 0,17 M which is 5,96 grams

1% (w/v) Bacto Tryptone which is 6 grams

0,5% yeast extract- 3 grams

ddH₂O- 600 ml

5M NaOH- 100ml (prepared from another group)

Autoclave

Result: 600 ml LB media

MONDAY, 8/7/2019

Experimenter: Johanna and Gabriel

Aim: To prepare an LB-media in order to use it for further lab-use

Protocol: Lab Manual, Synthetic Biology, LB medium p. 98

600ml of LB-medium was prepared with:

0,17 M (5,96 g) NaCl

1% (w/v) (6 g) Bacto Tryptone

0,5% (3g) yeast extract

600 ml ddH₂O

5M 100 µl NaOH

Autoclave

Result: 600 ml LB-media

#4 Making of 0,9% NaCl

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 17/6/2019

Experimenter: Vamsi, Karthik and Johanna

Protocol: Lab Manual, Synthetic Biology, Making of 0,9% NaCl page 91

Aim: making the 0,9% of NaCl solution

0,9% NaCl-10ml

0,09 grams NaCl

9,91ml ddH₂O

Result: 10 ml of 0,9% NaCl Solution

#3 Making of the 1 M CaCl₂

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 17/6/2019

Experimenter: Vamsi, Karthik and Johanna

Aim: Making of the CaCl₂ 1M solution

Protocol: Lab Manual Synthetic Biology, Preparation of Solution CaCl₂ page 93

ddH₂O to the 10 ml mark

1,1098 gram CaCl₂

Result: 1 M CaCl₂

1M CaCl₂ - 10 ml

#2 Preparation of LB medium + Agar + Kanamycin for 500 ml

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 17/6/2019

Experimenter: Vamsi, Karthik and Johanna

Aim: to prepare LB medium with agar to be able to make agar plates for growth of bacteria.

Protocol: Lab Manual Synthetic Biology, LB Agar Plates and addition of Antibiotics

NaCl 0,17 M- 4,967 grams

1% Bacto Tryptone- 5 grams

0,5% yeast extract- 2,5 grams

ddH₂O- 500 ml

5 M NaOH- 83,33 microlitre

Agar- 7,5 gram

Autoclave

Result:

Error: instead of agar we added agarose which made the LB medium agar unusable, therefore we are not continuing on making the agar plates.

#5Preparation of competent cells (DH5 α) using CaCl₂

Project: iGEM uppsala 2019

Authors: Vamsi Krishna Polavarapu

TUESDAY, 18/6/2019

Experimenter: Vamsi, Karthik and Johanna

Protocol: Lab Manual, Synthetic Biology, Page .no 111

Diluted the overnight culture of DH5 α 500 μ l to 50ml L.B medium and kept in 37 °C room

Measured the OD₆₀₀ till it reach 0.4.

Centrifuged the 50ml of the inoculum at 3500rpm for 5min at 4°C and decanted the supernatent.

Pellet was suspended in 100 μ l of icecold CaCl₂ and mixed well .Then added 15ml of icecold CaCl₂ ,left for 30mins incubation.

Again centrifuged at 3500rpm for 5min at 4°C and decanted the supernatent .

Added 2ml of 0.1 Micecold CaCl₂ mixed with 20% of glycerol and incubated it for 45 mins in Ice.

Prepared aliquotes of 50 μ l in 1.5 eppendorf tube.By putting them in Liquid No₂ ,snap freezed it.

Competent cell vials are stored in -80°C.

#6 Transformation of CaCl₂ competent E.Coli cells

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 19/6/2019

Experimenter: Johanna, Vamsi, Karthik

Aim: to transform the competent cells prepared from protocol #5 and to also plate them on agar plates.

Sample: Competent cells DH5 α 18.06.19 #5

Protocol: Lab Manual, Synthetic Biology, Transformation of CaCl₂-Competent E.Coli cells p. 113-115

Material:

5 ng PSB-1C3

50 μ l deionized water

5 μ l PSB-1C3 plasmid diluted 1:10 from stock solution that was 5 ng diluted with 50 μ l deionized water.

50 μ l competent cells snap frozen

50 μ l competent cells stored over night in cold room


950 μ l LB media


Agar plates

Stock solution of PSB-1C3 was diluted with 1 μ l stock solution and 9 μ l deionized water. This was then used with the competent cells.

Result: Agar plate with competent cells.

PSB-1C3_NOT_SNAP_Frozen.jpg





INSERT PHOTO OF COLONIES

DNA-plasmid PSB-1C3

Agar plate incubating in 37 degree celsius room.

WEDNESDAY, 3/7/2019

Experimenter: Johanna and Sofia (group 2)

Aim: To transform the two plasmids PSB1A3 and PSB1K3 into competent cells of DH5 α E.coli strains.

Protocol: Lab Manual, Synthetic biology, protocol #6 page 113

1 μ l PSB1A3

1 μ l PSB1K3

50 μ l DH5 α

3x 950 μ l SOB media

1 μ l ddH₂O

Sample: Escherichia Coli DH5 α from 13/9-2018

BBa- PSB1A3

BBa- PSB1K3

Result: 6 agar plates, 3 with Amp and 3 with Cm (wrong penicillin, so the plates did not grow).

The Amp grew to much, the plan is to re-do the plates with new ampicillin.

Stored: 37°C room for growth

MONDAY, 8/7/2019

Experimenter: Sofia and Johanna

Aim: To transform the wanted DNA into competent cells

Protocol: Lab manual synthetic biology, Protocol #6

Sample: PSB1A3

DH5αs

Result: Agar plates with transformed

Stored:

#7 Overnight Culture of K5-00000 and PSB-1C3 with chloramphenicol

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

THURSDAY, 20/6/2019

Experimenter: Johanna, Louise

Aim: To make overnight cultures with DH5 α bacteria with the plasmids K5-00000 and PSB-1C3.

Sample: DNA plasmid PSB-1C3 from 18.06.2019 protocol #6 from Victors labgroup

DNA plasmid K5-00000 protocol # 6 from Victors labgroup (19.06.2019)

LB medium

Chloramphenicol 1 μ l/ml

Protocol: Lab Manual, Synthetic Biology, Overnight Cultures with Antibiotics page 102-103

Different from the protocol: we mixed 12 ml LB media with 12 μ l chloramphenicol and then mixed it gently. This was then aliquoted into two 5 ml solutions for the K5-00000 colonies and 2 ml for the PSB-1C3 colonies.

The colonies were picked up with the tip of a pipet and then put into the falcontube for the corresponding colonie, the pipet was left over night with the solution.

Results: 3 tubes with colonies with LB medium + chloramphenicol + colonies from agar plates and was left for an overnight growth.

K5-00000 falcon tube stored in the 37 degree celsius room

PSB-1C3 falcon tube stored in the 37 degree celsius room.

8 Ligation of 3A-Assembly

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 24/6/2019

Experimenter: Johanna and Louise

Aim: To ligate the plasmids with the prevailing promoter (KB0-8006) with the ribosome binding site, coding sequence (lignin peroxidase) and the backbone with AmpR.

Protocol: Lab Manual, Synthetic Biology 3A Assembly p.107

Material:

2 μ l of each (the three) digestion mixtures

11 μ l ddH₂O

2 μ l 10x reaction buffer (for T4 DNA ligase)

1 μ l T4 DNA ligase

Sample: 3A Assembly

Result: Eppendorf tube with the samples.

Stored: used directly for transformation see protocol #10.

#11 GenElute Plasmid- Mini prep kit

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 24/6/2019

Experimenter: Johanna and Louise

Aim: To concentrate DNA from the samples prepared with the 3A assembly from protocol #8 (from Viktors group)

Protocol: GenElute Plasmid mini prep kit protocol.

Sample: K5-00000 in LB media O/N culture with 1:100 chloramphenicol

Result: 2 eppendorf tubes with concentrated DNA. The tube with sample #1 had 89,5 ng/μl and the tube with sample #2 had 113,5 ng/μl amount of DNA lignin peroxidase (K5-00000).

Stored: in -20°C freezer in two eppendorf tubes.

#10 Transformation of 3A assembly into DH5 α

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 24/6/2019

Experimenter: Johanna and Louise

Aim: To transform the ligated 3A assembly into DH5 α bacteria.

Protocol: Protocol #6

Material:

5 μ l ligation mixture (from protocol #8)

50 μ l DH5 α

950 μ l LB medium

Sample: 3A assembly from protocol #8 and from group 5s 3A assembly

Results: 3 Agar plates with AmpR, one with ddH₂O, one with non centrifuged bacteria and centrifuged ones.

Stored: in 37°C room for growth.

#9 Measurement of concentration from O/N culture from protocol #7

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 24/6/2019

Experimenter: Johanna and Louise

Aim: To measure the concentration of the O/N culture from protocol #7

Protocol: The nano-drop was cleaned with deionized water. The elution buffer was used as a blank for the samples. 2 μ l of the samples with K5-00000 was put on the measurer and then closed, the measuring was started and the values noted.

Material:

2 μ l ddH₂O

2 μ l Elute buffer

2 μ l Sample of K5-00000

Sample: K5-00000 from overnight culture from protocol #7

Results: The samples were tested for the concentration and sample #1 prepared by Johanna showed 89,5 ng/ μ l, sample #2 prepared by Louise showed 113,5 ng/ μ l.

Stored: Sample #1 and sample #2 was saved in a -20 °C freezer. The samples were marked with Lignation Peroxidase and the corresponding number.

#12 Make bacteria colonies on agar plate for PSB1A3 with K608006 from O/N 25.6.19

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 26/6/2019

Experimenter: Johanna and Vanja

Aim: To investigate whether the overnight culture had the relevant bacteria

Protocol: Lab Manual, Synthetic biology, protocol #6 Transformation of CaCl₂ competent E.coli cells from p.114 step 11-12

1 Agar plate

1 Coloni

Sample: PSB1A3, K608006 +Amp, colony #4 from 25.06.19

Result: Two agar plates with 100 µl each of the O/N cultures spread on Agar plates

Stored: 37°C room for growth.

#13 Overnight Culture with Ampicillin

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 26/6/2019

Experimenter: Johanna and Vanja

Aim: To make O/N culture from colonies on agarplate from 25/6-19

Protocol: Lab manual, Synthetic biology, protocol #2 Overnight Culture with antibiotic p.102.

Gloves were not used when the first coloni was taken, gloves was used for the second coloni.

5 µl ampicillin

1 coloni

5 ml LB media

Sample: PSB1A3, K608006 + Amp (Lignin peroxidase)

Result: 5 falcon tubes with one coloni each and: ampicillin, LB media.

Stored: 37°C room on shake tray, left over night.

#14 3A Assembly, Digestion of promotor and RBS

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 26/6/2019

Experimenter: Johanna and Vanja

Aim: To digest the promotor and restriction enzymes in order to have one the plasmids needed for the assembly.

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly

14µl K608006

29µl dd H₂O

5µl 10x reaction buffer

1µl EcoR1

1µl Spe1

Sample: K608006 + EcoR1 + Spe1

Result: One eppendorf tube with a plasmid with the promotor and RBS

Stored: in -20°C fridge (in lab)

#15 Running the gel

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 26/6/2019

Experimenter: Johanna and Vanja

Aim: To verify that the digestion from notebook #14 was succesful

Protocol: Lab Manual, Synthetic Biology, Protocol #3 Biobrik 3A assembly and gel analysis and Protocol #4 Agarose Gelectrophores

The gel was running with 120 V

well #1: 1µl ladder + 1µl loading dye + 4µl ddH₂O

well #2: 10 µl uncut plasmid + 10 µl ddH₂O + 4µl loading dye


well #3: 20 µl sample from group 5 (--.06.19) + 4µl loading dye

well #4: 20 µl sample from 26.06.19 + 4 µl loading dye

Sample: 3A assembly 26.06.19 from note#14 and 3A assembly from group 5 (--.06-19)

Result: A gel with clear lines of the samples. (insert photo)

Stored: Gel discarded



#17 Transformation of 3A Assembly

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

THURSDAY, 27/6/2019

Experimenter: Johanna and Vanja

Aim: to transform the earlier ligated 3A Assembly

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation

5 µl ligation mixture

50 µl competent cells, DH5α

950 µl LB medium

Sample: PSB1A3 + K68006 + Amp + Competent cells from iGEM 2018

Result: 3 Agar plates with the prevailing solutions

Stored: 37°C room for growth

Error: The incubation for 1-1,5 h was made with LB media and ligation mixture from entry #14, the error was detected when the agar plates were being plated. Therefore the incubation on ice lasted not for 30 minutes but for approximately 3,5 h before the cells were heat shocked and the incubation on 37°C and the protocol was continued

#18 Thermo scientific GenJet plasmid miniprep kit, from O/N culture 26.06.19

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

THURSDAY, 27/6/2019

Experimenter: Johanna and Vanja

Aim: To concentrate the DNA from the overnight culture

Protocol: Thermo Scientific GenJet plasmid miniprep kit

Sample: O/N PSB1A3 from 26.06.19

Result: 4 eppendorf tubes with the different colonies for O/N cultures.

Stored: Discarded

Nanodrop:

Experimenter:

Aim: To examine the concentration of the DNA in the solution

Protocol: Blank with elution buffer, and then measure 2 μ l of sample

Sample: PSB1A3

Result:

#6 had -3,8 ng/ μ l DNA

#8 had 7,5 ng/ μ l DNA

#9 had 27,5 ng/ μ l DNA

#10 had 6,9 ng/ μ l DNA

Stored: Discarded, since there were to little DNA in the samples.

#16 Ligation for 3A assembly from entry #14

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

THURSDAY, 27/6/2019

Experimenter: Johanna and Vanja

Aim: To ligate the three digestion mixture into one plasmid

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick assembly ligation

2 µl of each digestion mixture + 11 µl ddH₂O

2 µl 10x reaction buffer for T4 DNA ligase

1 µl T4 DNA ligase

Sample: Sample #2 from 20.06.19 (group 5) and sample #3 from 20.06.19(group 5) and sample #1 from entry #14

Result: A sample with, hopefully, the three digestion mixtures ligated together

Stored: Used directly for transformation

Enzymatic Assay- Piercy Quantitative Peroxide Assay Kit [23280]

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

MONDAY, 1/7/2019

Experimenter: Johanna and Vamsi

Aim: To investigate whether the kit that had arrived would work with our enzymes and to test the kit

Protocol: Piercy Quantitative Peroxide Assay Kit, the protocol for the aqueous solution.

A stock solution of H_2O_2 was prepared. H_2O_2 was diluted from a 35% solution to a 1 M solution. It was done by mixing 9,7 μ l of the 35% solution with 99,99 μ l of dd H_2O .

The 1 M solution was diluted to a:

1 mM solution: 0,9 ml H_2O_2 and 0,1 ml H_2O

0,9 M solution: 0,8 ml H_2O_2 and 0,2 ml H_2O

0,8 M solution: 0,7 ml H_2O_2 and 0,3 ml H_2O

0,7 M solution: 0,6 ml H_2O_2 and 0,4 ml H_2O

0,5 M solution: 0,5 ml H_2O_2 and 0,5 ml H_2O

0,4 M solution: 0,4 ml H_2O_2 and 0,6 ml H_2O

0,3 M solution: 0,3 ml H_2O_2 and 0,7 ml H_2O

0,2 M solution: 0,2 ml H_2O_2 and 0,8 ml H_2O


0,1 M solution: 0,1 ml H_2O_2 and 0,9 ml H_2O

1 μ M solution: 0,001 ml H_2O_2 and 9,999 ml H_2O

Result:

The table are from a spectrophotometry with the absorbance measured at 560 nm

	A	B	C	D
1	Samples micro molar	Spec #1 [15 min incubation]	Spec #2 [20+ min incubation]	Spec #3 [30 min incubation]
2	1000	0,078	0,078	0,009
3	700	0,067	0,071	0,002
4	400	0,059	0,062	-0,011
5	200	0,054	0,058	-0,016
6	10	0,048	0,053	-0,018



Stored: The not used concentrations were left in the fridge along with 1M stock solution.

TUESDAY, 2/7/2019

Experimenter: Johanna, Gaberial and Vamsi

Aim: To perform the enzymatic assay by using pierce quantitative peroxide assay kit

Procedure: Followed the protocol provided in the kit (Thermo fisher - protocol 23280)

From the stock solution 35% of H₂O₂ we prepared 1M and 1mM stock solutions.

We later diluted the 1mM H₂O₂ stock solutions to 900μM, 800μM, 700μM, 500μM, 400μM, 300μM, 200μM, 100μM, 10μM and 1μM concentrations.

By mixing 120ul of reagent A in 12 ml of reagent B on ice we prepared working reagent for 12 samples.


We mixed 100ul of H₂O₂ sample and 1ml of working reagent. We did it for all 11 samples.

After 15 minutes of incubation in 37 °C, we measured the OD at 590nm.

Table2

	A	B	C
1	Samples	Time 15 mins	Time 25 mins
2	1000µM	0.025	0.009
3	900µM	0.022	0.004
4	800µM	0.019	0.003
5	700µM	0.017	0.001
6	500µM	0.012	-0.005
7	400µM	0.009	-0.008
8	300µM	0.004	-0.011
9	200µM	0.002	-0.015
10	100µM	0.000	-0.016
11	10µM	-0.002	-0.017
12	1µM	-0.003	-0.019

A graph from the 11 samples.



Enzymatic Assay with Platerreader Tecan

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

TUESDAY, 9/7/2019

Experimeter: Johanna and Gabriel

Aim: To investigate whether the kit that was ordered worked on platereaders as well.

Protocol: Piercy Quantitative Peroxide Assay Kit

Concentrations:

1000 μM



700 μM

400 μM

200 μM

10 μM

Result:



	A	B	C
1	Sample	Concentration [μ M]	Absorbance
2	1	1000	0,1688
3	2	700	0,1628
4	3	400	0,157
5	4	200	0,1296
6		10	0,142

Biobrick Assembly

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 10/7/2019

Experimenter: Johanna and Gabriel

Aim: To get the desired DNA from the iGEM plates and directory.

Protocol:

10 µl ddH₂O

2-3 pg/µl DNA from plates

Punch hole on the aluminum cover, discard first pipettip. Take the ddH₂O and pipet up and down in the well with DNA. Let sit for 5 min.

Take the sample from the well and put it in eppendorf tubes.

Sample:

asPink (BBa_K1033933)- plate 4, well 21C

amilCp -blur/purple chromoprotein (BBa_K103930)- plate 6, well 9L

Promotor (BBa_J23119)- plate 2, well 18P

Promotor and RBS (BBa_K081005)- plat 3, well 5E

Result: 4 tubes of the collected DNA, asPink (BBa_K1033933), amilCp-blue/purple chromoprotein (BBa_K103930), Promotor for ribosome standby site (BBa_J23119), promotor and RBS (BBa_K081005)

Stored: in 4°C fridge

THURSDAY, 11/7/2019

Transformation of the plasmids from 10/7-19

Experimenter: Johanna, Gabriel and Vamsi

Aim: To transform the appropriate plasmids into E.Coli DH5α

Protocol: Lab Manual, Synthetic Biology, protocol #6

instead of adding 5 µl of plasmid we used 1 µl of the plasmid from the samples aquired yesterday from the iGEM registry.

Sample:

300 µl E.Coli, DH5α

1 µl asPink (BBa_K1033933)

1 µl amilCp -blur/purple chromoprotein (BBa_K103930)

1 µl Promotor (BBa_J23119)

1 µl Promotor and RBS (BBa_K081005)

Result: 10 agar plates with chlorampicillin, 5 were with the samples and a negative control that was not centrifuged and 5 were the samples and negative control but centrifuged.

Stored: in 37°C room for growth

Making of LB plates with chlorampicillin

Experimenter:

Aim: To make agar plated with the correct antibiotic

Protocol: Lab manual, Synthetic Biology protocol #1 page 100, LB agar plates and addition of antibiotics.

Result: Agar plates with chlorampicillin

Stored: stacked together in the 4°C fridge

Preparation of competent *E.coli* (BL21)

Experimenter: Gabriel, Johanna

Aim: Start making competent cells.

Protocol: Lab manual, Synthetic Biology protocol #5 page 111. One colony of BL21 was taken from plate and inoculated in 5 mL SOB.

Incubated at 37° with shaking.

Result: Overnight culture was started.

Stored: in 37°C room for growth.

FRIDAY, 12/7/2019

Prepare 0,1 M CaCl₂ and 0,1 CaCl₂ with 20% glycerol

Experimenter: Johanna

Aim: To prepare the correct amount for the different solutions in order to use them for competent cells

Protocol: Lab manual, Synthetic Biology protocol # 1 page 92-93

Result:

	A	B	C	D	E
1	Sample	1 M CaCl ₂	50% glycerol	ddH ₂ O	Total Volume
2	0,1 M CaCl ₂	4 ml	-	36 ml	40 ml
3	0,1 M CaCl ₂ + 20% glycerol	0,6 ml	2,4 ml	3 ml	6 ml

Making competent cells of *E.Coli* BL21

Experimenter: Johanna, Vamsi and Gabriel

Aim: To make competent cells with BL21 to be able to use it for future transformation

Protocol: Lab Manual, Synthetic biology protocol#5, page 111

Sample: BL21 from may from the overnight culture

Result:

Stored: In -80° freezer

Prepare ON culture of transformed cells

Experimenter: Johanna, Vamsi and Gabriel

Protocol: Start overnight culture by inoculating colonies into 5mL of SOB media with shaking at 37°C

Sample: Plates with transformed cells of:

1 µl asPink (BBa_K1033933)

1 µl amilCp -blur/purple chromoprotein (BBa_K103930)

1 µl Promotor (BBa_J23119)

1 µl Promotor and RBS (BBa_K081005)

Result: 4 falcon tubes tubes with 5mL of SOB containing one of the transformed cells.

Stored: In 37°C room for growth.

SATURDAY, 13/7/2019

Plasmid purification - GeneJET plasmid Miniprep

Experimenter: Gabriel, Johanna

Aim: To concentrate plasmid DNA from transformed cells

Protocol: GeneJET plasmid Miniprep kit X2

Sample: 4 O/N culture in LB media of cells with one of following plasmids each:

asPink (BBa_K1033933)

amilCp -blur/purple chromoprotein (BBa_K103930)

Promotor (BBa_J23119)

Promotor and RBS (BBa_K081005)

Result: 2 eppendorfs with purified plasmid of each plasmid. 8 tubes in total.

Gabriels samples:

46,7 ng/µl asPink (BBa_K1033933)

42,6 ng/µl amilCp -blur/purple chromoprotein (BBa_K103930)

30,1 ng/µl Promotor (BBa_J23119)

38,5 ng/µl Promotor and RBS (BBa_K081005)

Johannas samples:

53,7 ng/µl asPink (BBa_K1033933)

39,1 ng/µl amilCp -blur/purple chromoprotein (BBa_K103930)

35,9 ng/µl Promotor (BBa_J23119)

37,8 ng/µl Promotor and RBS (BBa_K081005)

Stored: in -20°C freezer in 4x2 eppendorf tubes.

TUESDAY, 16/7/2019

Begin 3A assembly

Digestion

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly.

Chromoproteins digested with E and X

Promotors digested with E and S

0,7 agarose gel was used

Sample:

Gabriels samples:

46,7 ng/µl asPink (BBa_K1033933)

42,6 ng/µl amilCp -blur/purple chromoprotein (BBa_K103930)

30,1 ng/µl Promotor (BBa_J23119)

38,5 ng/µl Promotor and RBS (BBa_K081005)

Result:

3 eppendorfs with 30 µl of 100 ng/µl (300 ng total) of digested plasmids:

asPink (BBa_K1033933)

amilCp -blur/purple chromoprotein (BBa_K103930)

Promotor (BBa_J23119)

Stored: In -20° freezer.

Digestion 2

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly.

Promotor and RBS (BBa_K081005) digested with S+P and E+S

Codon optimized AsPink digested with E+S

Sample:

38,5 ng/µl Promotor and RBS (BBa_K081005)

Codon optimized AsPink (optimized BBa_K1033933)

Result:

Successful digestion according to gel.

3 eppendorfs

(1) Promotor and RBS (BBa_K081005) digested with S+P

(2) Promotor and RBS (BBa_K081005) digested with E+S

(3) Codon optimized AsPink digested with X+P

Stored: In -20° freezer.

Ligation

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly

AmilCP(Digested with E+X) + Promotor(no RBS)(Digested with E+S)

AsPink (non-optimized)(Digested with E+X) + Promotor-RBS(digested with E+S)(2)

AsPink (optimized)(Digested with X+P) + Promotor-RBS(digested with S+P)(1)

Sample:

Digested plasmids in eppendorf:

AmilCP (Digested with E+X)

AsPink (non-optimized) (Digested with E+X)

AsPink (optimized) (Digested with X+P)

Promotor(no RBS) (Digested with E+S)

Promotor-RBS(digested with E+S) (2)

Promotor-RBS(digested with S+P) (1)

Result: 3 eppendorfs with 20µl of ligated products.

Stored: In -20° freezer

WEDNESDAY, 17/7/2019

Transformation using ligated products

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation of CaCl₂-competent *E.coli* cells (X2)

Sample: Ligated products in eppendorfs:

- AmilCP(Digested with E+X) + Promotor(no RBS)(Digested with E+S)
- AsPink (non-optimized)(Digested with E+X) + Promotor-RBS(digested with E+S)(2)
- AsPink (optimized)(Digested with X+P) + Promotor-RBS(digested with S+P)(1)
- Competent DH5-alpha *E.coli*

Result: 7 LB plates (Each ligation X2 + 1 negative control)

Stored: In 37° room

THURSDAY, 18/7/2019

Digestion 3

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly

Digest "asPink (BBa_K1033933)" and "amilCp -blur/purple chromoprotein (BBa_K103930)" with X+P.

Sample: Purified plasmids

Johannas samples:

53,7 ng/µl asPink (BBa_K1033933)

39,1 ng/µl amilCp -blur/purple chromoprotein (BBa_K103930)

Result:


Digested: "asPink (BBa_K1033933)" and "amilCp (BBa_K103930)". In eppendorfs with concentration of 10 ng/µl.

Stored: In -20° freezer

GEL

6 samples+ ladder.

Ladder, asPink(control+digest), amilCp(control+digest), Two different PSB1A3 with RFP (digested).



Phosphorylation of primers(amilCP)

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #10 Inverse PCR mutagenesis

Sample:

FWD primer

REV primer

Result: Phosphorylated primers in eppendorfs with primer concentration of 5 μ M. 20 μ l total volume.

Stored: In -20° freezer

FRIDAY, 19/7/2019

Ligation

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly

AmilCP(Digested with X+P) + Promotor(no RBS)(Digested with E+S) + PSB1A3 with RFP (digested with E+P)

AsPink (non-optimized)(Digested with X+P) + Promotor-RBS(digested with E+S)(2) + PSB1A3 with RFP (digested with E+P)

AsPink (optimized)(Digested with X+P) + Promotor-RBS(digested with E+S)(1) + PSB1A3 with RFP (digested with E+P)

Sample:

Digested plasmids in eppendorf:

AmilCP (Digested with X+P)

AsPink (non-optimized) (Digested with X+P)

AsPink (optimized) (Digested with X+P)

Promotor(no RBS) (Digested with E+S)
Promotor-RBS(digested with E+S) (2)
PSB1A3 with RFP (digested with E+P)

Result: 3 eppendorfs with 20µl of ligated products

Stored: -20° freezer

Transformation

Experimenter: Gabriel, Katherin

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation of CaCl₂-competent *E. coli* cells

Number of plates prepared of which sample. Total of 8 plates:

3x AsPink (optimized)
2x AsPink (non-optimized)
2 AmilCP
1x H₂O (control)

Sample: Ligated products in eppendorfs:

1. AmilCP(Digested with X+P) + Promotor(no RBS)(Digested with E+S) + PSB1A3 with RFP (digested with E+P)
2. AsPink (non-optimized)(Digested with X+P) + Promotor-RBS(digested with E+S)(2) + PSB1A3 with RFP (digested with E+P)
3. AsPink (optimized)(Digested with X+P) + Promotor-RBS(digested with E+S)(1) + PSB1A3 with RFP (digested with E+P)

Result: 8 plated amp plates

Stored: 37°room

WEDNESDAY, 24/7/2019

Transformation

Experimenter: Gabriel, Katherin

Aim: Insert plasmids into competent *E. coli*

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation of CaCl₂-competent *E. coli* cells

Number of plates prepared of which sample. Total of 3 plates:

Sample: pSB1C3 with GFP, pSB1C3 with asPink, competent DH5α *E. coli*

Result: No growth on plates

Making competent cells of *E. coli* DH5α

Experimenter: Gabriel, Katherin

Aim: To make competent cells with DH3to be able to use it for future transformation

Protocol: Lab Manual, Synthetic biology protocol#5, page 111

Sample: plate with DH5α from fridge

Result: ON prepared

Stored: 37°room

THURSDAY, 25/7/2019

Transformation

Experimenter: Gabriel, Katherin, Qian

Aim: Insert plasmids into competent *E. coli*

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation of CaCl₂-competent *E. coli* cells

Sample: pSB1C3 with GFP, pSB1C3 with asPink, competent DH5α *E. coli*

Result: 6 plates, 2x GFP, 2x asPink, 1x H₂O

Stored: 37° room

Digestion with DpnI

Experimenter: Gabriel, Katherin

Aim: Remove plasmid template, from PCR product.

Protocol: Lab Manual, Synthetic biology page 129

Sample: Purified product of amilCP PCR

Result: Digestion performed on PCR product sample

Stored: -20° freezer

Making competent cells of *E. coli* DH5α

Experimenter: Gabriel, Katherin

Aim: To make competent cells with DH3to be able to use it for future transformation

Protocol: Lab Manual, Synthetic biology protocol#5, page 111

Sample: ON culture

Result: Eppendorfs with 100μl of competent cells

Stored: -80° freezer

FRIDAY, 26/7/2019

Transformation

Experimenter: Gabriel, Katherin

Aim: Insert plasmid pSB1C3 with cjBlue(BBa_k1073026) into competent DH5α *E. coli*

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation of CaCl₂-competent *E. coli* cells

Sample: pSB1C3 with cjBlue, 2 different batches of competent DH5α *E. coli* (*Gabriel and Irina*)

Result: 4 plates: 2x cjBlue(Gabriel and Irinas *E. coli*), 2x H₂O (Gabriel and Irinas *E. coli*)

Stored: 37° room

MONDAY, 29/7/2019

Ligation

Transformation

Experimenter: Gabriel, Katherin, Johanna, Jonas

Aim: Insert plasmid pSB1C3 with cjBlue(BBa_k1073026) into competent DH5α *E. coli*

Protocol: Lab Manual, Synthetic Biology, protocol #6 Transformation of CaCl₂-competent *E. coli* cells

Sample: pSB1C3 with cjBlue, 2 different batches of competent DH5α *E. coli* (*Gabriel and Irina*)

Result: 4 plates: 2x cjBlue(Gabriel and Irinas *E. coli*), 2x H₂O (Gabriel and Irinas *E. coli*)

Stored: 37° room

Ligation

Experimenter: Gabriel and Johanna

Aim: To ligate the PCR sequences in order to make them circular (plasmid)

Protocol: Lab Manual, Synthetic Biology, protocol #10 page 130-131

one tube had both p1 and p2 in it and where therefore doubled in volume since the concentration would otherwise be too high. Tubes with just p1 and p2 were also prepared.

	A	B	C	D
1	Chemical	Dpnl p1	Dpnl p2	Dpnl p1 and p2 mix
2	Dpnl p1	2 µl	-	2 µl
3	Dpnl p2	-	2 µl	2 µl
4	T4 ligation buffer	5 µl	5 µl	10 µl
5	T4 ligation	1 µl	1 µl	2 µl
6	ddH ₂ O	42 µl	42 µl	84 µl

Sample: Dpnl p1 and p2 from group 5, the samples are the same sequence but labeled different.

Result: One tube with Dpnl p1, one tube with Dpnl p2 and one tube with the mix of Dpnl p1 and p2

Concentration of Dpnl p1 was 24,0 ng/µl

Concentration of Dpnl p2 was 24,9 ng/µl

Stored: -20° freezer in the group 4 box.

O/N Culture of cjBlue from 26/7-19

Experimenter: Johanna

Aim: To make an O/N culture to make a miniprep the next day in order to concentrate the DNA

Protocol: Lab Manual Synthetic Biology, protocol #2 page 102

5 ml LB media

5 µl Chloroampicillin

1 colony from plate

Sample: cjBlue from transformation 26/7-19

Result: Two falcon tubes with one colony each

Stored: 37°C room overnight (the room is at 35°C)

TUESDAY, 30/7/2019

Miniprep of the transformed cjBlue from 29/7-19

Experimenter: Johanna

Aim: To concentrate the DNA of cjBlue

Protocol: Thermo Scientific GenJET Plasmid miniprep Kit

Sample: cjBlue from transformation 29/7-19

Result: 4 eppendorf tubes 2 different colonies.

cjBlue #1: 29,9 ng/μl

cjBlue #1 40,9 ng/μl

cjBlue #2 50,1 ng/μl

cjBlue #2 46,0 ng/μl

Stored: -20°C freezer

Transformation

Experimenter: Johanna and Gabriel

Aim: To transform cjBlue into *E.coli* BL21

Protocol: Lab Manual, Synthetic Biology, protocol #6

Sample: cjBlue from miniprep made 30/7-19

Result: 4 agar plates that were spread with cjBlue

Stored: 37°C

Making of Agar plates

Experimenter: Johanna and Gabriel

Aim: To make agar plates with chloramphenicol for future use

Protocol: Lab manual, Synthetic Biology protocol #1 page 100, LB agar plates and addition of antibiotics.

Result: 23 plates of agar with chloramphenicol

Stored: 4°C

Overnight Culture

Experimenter: Johanna

Aim: To make an overnight culture in order to do a miniprep the following day.

Protocol: Lab Manual Synthetic Biology, protocol #2 page 102

Sample: cjBlue, DpnI p1 and DpnI p1+p2

Result: 3 falcon tubes with O/N cultures

Stored: 37°C

WEDNESDAY, 31/7/2019

Miniprep

Experimenter: Gabriel

Aim: To concentrate the DNA of cjBlue, P1 and P1+P2 mix

Protocol: GeneJET plasmid Miniprep kit. Cells were harvested from 1 ml of ON culture. Concentration measured with NanoDrop.

Sample: Overnight cultures of transformed cells with cjBlue, DpnI p1 and DpnI p1+p2

Result:

P1+P2 mix 53,6 ng/μl

cjBlue: 37,6 ng/μl

P1: 47,3 ng/μl

Stored: -20° freezer in group 4 box

Miniprep

Experimenter: Gabriel, Johanna

Aim: To concentrate the DNA of cjBlue, P1 and P1+P2 mix, harvesting cells from 2ml instead of 1ml as used earlier this day.

Protocol: GeneJET plasmid Miniprep kit. Cells were harvested from 2 ml of ON culture. Concentration measured with NanoDrop.

Sample: Overnight cultures of transformed cells with cjBlue, Dpnl p1 and Dpnl p1+p2

Result:

P1: 124 ng/μl

P1+P2: 143,4 ng/μl

cjBlue: 91,1 ng/μl

Stored: -20° freezer in group 4 box

Preparation of sample (P1) for sequencing

Experimenter: Gabriel, Johanna

Aim: Prepare sample for sequencing

Protocol: Followed protocol on package of Mix2Seq from Eurofins.

Primers used: VR, VF2

Sample: Concentrated P1 from Miniprep (124 ng/μl)

Result: One tube ready for sequencing

Stored: Sent for sequencing

Restreak of transformed E. coli

Experimenter: Johanna

Aim: Too visually observe the rate at which cells express cjBlue.

Protocol: Picked 1 colony

Sample: cjBlue colonies from transformation performed on friday 26/7.

Result: 1 plate with 1 streaked colony

Stored: 37°C room

Overnight Culture

Experimenter: Johanna

Aim: To make an overnight culture in order to observe rate of color expression.

Protocol: Lab Manual Synthetic Biology, protocol #2 page 102

Sample: cjBlue colonies from transformation performed on friday 26/7.

Result: 1 falcon tube with O/N cultures

Stored: 37°C room

THURSDAY, 1/8/2019

Digestion of P1(amiICP with RSS) with P and X

Experimenter: Gabriel, Johanna


Aim: Digest P1 for ligation with PSB1C3 and Promotor (BBa_J23119).

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly.

Sample: Eppendorf with P1 with concentration of 124 ng/μl.

Result: Succesfull digestion according to gel.

From left: Digested P1, undigested plasmid, ladder.



Stored: In -20° freezer.

Ligation

Experimenter: Gabriel, Johanna

Aim: To ligate digested P1, PSB1C3 and Promotor (BBa_J23119).

Protocol: Lab Manual, Synthetic Biology, protocol #3 Biobrick, 3A assembly.

Sample: digested P1, PSB1C3 and Promotor (BBa_J23119).

Result: Eppendorf with ligation mixture.

Stored: In -20° freezer.

Transformation

Experimenter: Gabriel, Johanna

Aim: To transform ligation mixture of P1, PSB1C3 and Promotor (BBa_J23119) into DH5 α

Also transform eforRed (BBa_K864402) into DH5 α

Protocol: Lab Manual, Synthetic Biology, protocol #6

Sample: eforRed, ligation mixture of P1, PSB1C3 and Promotor. Competent DH5 α , 3 amp plates, 3 chloramphenicol plates.

Result: 6 plates: 2 ligation mixture, 2 eforRed, 2 controls (one for each type of plate).

Stored: 37°C room

FRIDAY, 2/8/2019

Transformation

Experimenter: Gabriel and Johanna

Aim: To transform eforRed, amilCp and asPink into DH5 α

Protocol: Lab manual, Synthetic Biology protocol #6

Sample: eforRed, amilCp and asPink on chloramphenicol plates

Result: 3 agar plates with the corresponding plasmids

Stored: 37°C room

Sequencing

Experimenter: Gabriel and Johanna

Aim: To sequence the plasmid

Protocol: Mix2Seq Kit

Sample: DpnI p1 (amilCp)

Result: One tube of 15 µl DNA plasmid and 2µl VF2 primer.

A sequence was returned that showed that our overhang was successful in the insertion. One nucleotide was missing for the Xba1 restriction site which means that we have to order new primers and redo a PCR in order to submit the biobrick to iGEM.

MONDAY, 5/8/2019

Restreak

Experimenter: Johanna and Gabriel

Aim: To restreak the colonies that should have colour, since they are a high number of colonies on the plate.

Protocol: Agar plates with chloramphenicol

one colonie was taken from each agar plate with colonies and then restreaked onto the new agar plates.

Sample: amilCp and asPink

Result: two agar plates with restreaking of colonies

Overnight culture

Experimenter: Johanna and Gabriel

Aim: to make overnight cultures

Protocol: Lab Manual Synthetic Biology, protocol #2 page 102

Sample: asPink and amilCp from transformation

Result: two overnight cultures with SOB media. one with asPink and the other with amilCp

Stored: 37°C room

TUESDAY, 6/8/2019

Miniprep. Thermo Scientific GeneJET plasmid miniprep Kit

Experimenter: Johanna

Aim: to concentrate DNA

Protocol: Thermo Scientific GeneJET plasmid miniprep kit

Sample: O/N culture from 5/8-19

Result: two eppendorf tubes with concentrated DNA.

asPink: 65,5 ng/µl

amilCp: 84,4 ng/µl

Stored: -20°C freezer

WEDNESDAY, 7/8/2019

Ligation

Experimenter: Gabriel and Johanna

Aim: to ligate different genes that are desired to be in one plasmid

Protocol: Lab manual, Synthetic Biology, protocol #3 page 107

Sample:

amilCp

asPink

Promotor

Promotor + RBS

Result: 2 eppendorf tubes with ligation mixtures of asPink and amilCp

Stored: -20°C freezer

Transformation

Experimenter:

Aim: To transformed the ligation mixtures that was

Protocol: lab manual, Synthetic Biology, protocol #6

Sample: Ligation mixtures of asPink and amilCp

Result: 2 agar plates with ampicillin, one with amilCp and one with asPink (one negative as well).. No plate had any colonies, we suspect the plates

Stored:37°C room

PCR

Project: iGEM uppsala 2019

Authors: Johanna Cederblad


TUESDAY, 30/7/2019

Experimenter: Johanna

Aim: To multiply DNA

Protocol: Thermo Scientific DreamTaq DNA polymerase

A temperature gradient was used for the samples (50-60°C)



Sample: pPicz

Result: 8 tubes into the PCR machine

Stored: Left overnight in the PCR machine and put in a gel the day after

WEDNESDAY, 31/7/2019

Agarose gel to detect the DNA from PCR

Experimenter: Johanna

Aim: To detect whether the PCR was succesful and if the temperature gradient was useful

Protocol: Lab manual, Synthetic Biology protocol #4 page 108-112

Sample: the pPICz from the PCR done 30/7-19

Result: Nothing showed on the gel, due to probable mix-up with the tubes

Stored: Discarded

Degrading of Lignin- Proof of concept

Project: iGEM uppsala 2019

Authors: Johanna Cederblad

WEDNESDAY, 7/8/2019

Fractionating lignin

Drag and drop a file from the sidebar to attach or

Upload a File

Experimenter: Johanna and Gabriel

Aim: to collect fractions of the lignin in order to degrade it further in an easier way.

Protocol:

Methanol fraction:

2 gram lignin and 20 ml methanol was stirred for 2 h (500 varv/min) and then filtered with level 1 filterpaper, the supernatant was saved.

2 gram lignin and 20 ml methanol was stirred for 2 h (500 varv/min) and then centrifuged for 5 min in 4000 rpm. the supernatant was saved.

Aceton fraction:

the pellet from the filter fraction was mixed with 20 ml aceton and then stirred for 2h (500 varv/min) and then filtered through a level 3 filter paper.

the pellet from the centrifuged samples was mixed with 20 ml aceton and then stirred for 2h (500 varv/min) and then centrifuged for 10 min in 4000 rpm. supernatant and pellet was saved seperatly.

Sample: Lineo Prime W

Result:

Methanol:

One falcon tube with filtered supernatant

one falcon tube with centrifuged supernatant

Acetone:

one falcon tube with filtered supernatant

one falcon tube with centrifuged supernatant

one falcon tube with the pellet from the centrifuged sample.

two eppendorf tubes with the scraped pellet from the filterpaper.

Stored: room temperature

THURSDAY, 8/8/2019

Native-PAGE

Experimenter: Johanna and Gabriel

Aim: To detect whether the lignin is detectable in a gel.

Protocol:

Sample buffer(2X):

49µl 1,5 M tris HCl, pH 8,8

50µl 25% glycerol

1µl 1% bromophenol Blue

All was mixed into eppendorf tube.

Running buffer(10X):

275 ml ddH₂O

125 ml Tris HCl pH8,8

100 ml glycerol (1,92M)


All was mixed and diluted to 1X when need to run the gel.

For running the gel 10µl sample was mixed with 10µl sample buffer

Sample: Filtered MeOH, Filtered Acetone, Centrifuge MeOH, Centrifuged Acetone, Original Lignin, Lignin pellet from Acetone

Result:

 Native-PAGE.png



SDS-PAGE 20% separation gel, 4% stacking gel

Experimenter: Johanna and Gabriel

Aim: To see if lignin is visible on an SDS-PAGE

Protocol:

Separation gel:

0,59 ml ddH₂O

10/20/2019

Degrading of Lignin- Proof of concept · Benchling

6,7 ml acrylamid/bis (30%)
2,5 ml Tris-HCl (1,5 M, pH 8,8)
100 µl SDS (10%)
10 µl TEMED
100 µl APS


Stacking gel:

6,1 ml ddH₂O
1,3 ml Acrylamide/bis (30%)
Tris-HCl (0,5M pH6,8)

Sample:

Result:

 Screenshot_#1_20%_SDS-PAGE.png



MONDAY, 12/8/2019

Fractionating of lignin in high pH

Experimenter: Johanna and Gabriel

Aim: To make the kraft lignin (Lineo Prime) soluble with the help of high pH


Protocol: 2 gram of Lineo Prime was mixed with 20 ml of 0,5M NaOH and stirred over night
a lump of pellet from the acetone fraction was mixed with 20 ml of 0,5 M NaOH

Sample: Lineo Prime W

Result: Two 50 ml falcon tubes with soluble lignin.

Stored: On bench in room temperature

Lignin_NaOH.png



TUESDAY, 13/8/2019

Degrade Lignin with HRP in Phosphate buffer

Experimenter: Johanna and Gabriel

Aim: To degrade lignin with HRP and H₂O₂ in phosphate buffer

Protocol: Lignin blandades i olika lösning (se tabell) and were tested with a nanodrop.

Table2


	A	B	C	D	E	F	G
1	Sample	Lignin	H ₂ O	Phosphate buffer	H ₂ O ₂	Heat-inactivated enzyme	ICA HRP
2	A	100 µl	9,9 ml	-	-	-	-
3	B	100 µl		9,9 ml	-	-	-
4	C	100 µl		9,74 ml	-	0,16 ml	-
5	D	100 µl		9,74 ml	-	-	0,16 ml
6	E	100 µl		9,41 ml	0,33 ml	-	0,16 ml
7	F	100 µl		9,41 ml	0,33 ml	0,16 ml	-
8	G	100 µl		9, 57 ml	0,33 ml	-	-

Sample: Sulfonated lignin

Result: Pellet in the samples that had the ICA HRP, both in the denaturated HRP and the regular one. The denaturated lignin might not be completely inactive.

samples were taken every hour. These were loaded on an SDS-PAGE later.

Screenshot_Lignin_samples_phosphate_buffer.png



Run of SDS-PAGE of lignin samples 15%

Experimenter: Johanna and Gabriel

Aim: To see if any pattern of lignin could be detected on the gel.


Protocol: The samples were mixed with running dye and then loaded into the wells of the gel. The gel was started at 80 V and then increased to 130 V after the samples had run down the gel.

Sample: The samples from the degrading of lignin.

Result: a gel se picture

[Ladder: E0:E1:E2:E4:kill_it_0:kill_it_1:A:B:C:D:F:G]

Screenshot_#2_15%_SDS-PAGE.png



WEDNESDAY, 14/8/2019

SDS-PAGE #3 of lignin samples 15%

Experimenter: Johanna and Gabriel

Aim: To detect difference in the gel and whether the lignin has been degraded.


Protocol: The samples were mixed with running dye and then loaded into the wells of the gel. The gel was started at 80 V and then increased to 130 V after the samples had run down the gel.

Sample: from lignin degradation earlier

Result:

[ladder: centrifuged sonicated acetone: centrifuged sonicated MeOH: Filter MeOH:Filter Acetone: centrifuge MeOH: bubble: centrifuge acetone: original Lineo prime: Original sulfonated]

Screenshot_#3_15%_SDS-PAGE.png



Sonication of lignin

Experimenter: Johanna

Aim: To degrade the lignin with sonication

Protocol: 50% amplitude, 5 seconds puls 5 seconds rest, for 30 min

Sample: sulfonated lignin

Result: a falcon tube with sonicated lignin

THURSDAY, 15/8/2019

Degrade lignin with HRP in Tris-HCl buffer

Experimenter: Johanna and Gabriel

Aim: To degrade lignin with HRP

Protocol: lignin was mixed with different solutions in order to see if any degradation would occur.

Table3

	A	B	C	D	E	F	G	H
1	Sample	Tris-HCl	ICA HRP	Sigma HRP	Heat-inactivated ICA	Heat-inactivated Sigma	H2O2	Lignin
2	A	9,74 ml	-	160 µl	-	-	-	100 µl
3	B	9,74 ml	160 µl	-	-	-	-	100 µl
4	C	9,41 ml	-	160 µl	-	-	330 µl	100 µl
5	D	9,41 ml	160 µl	-	-	-	330 µl	100 µl
6	E	9,41 ml	-	-	-	160 µl	330 µl	100 µl
7	F	9,41 ml	-	-	160 µl	-	330 µl	100 µl
8	G	9,90 ml	-	-	-	-	-	100 µl

Sample: sulfonated lignin

Result: Pellet and sign of change in the ICA HRP, otherwise some decrease in curves.

Prepare 0,1 M tris-buffer

Experimenter: Johanna

Aim: To get the correct M of the buffer

Protocol: diluted a 1 M tris-buffer to a 0,1 M tris-buffer

Sample: 1 M tris- buffer from group 2

Result: a 0,1 M tris-buffer

SDS-PAGE 18%

Experimenter: Johan Erkers

Aim: To make SDS-plate for analysis of lignin samples.

Protocol:

The following protocoll was used:

☰ SDS Page Gel (10 and 15%)


The amount of sample used was:

Table1		
	A	B
1	Material	Quantities for 2 Gels
2	water	1.334ml
3	Akrylamide	6.066ml
4	1,5 M Tris pH 8.8	2,5ml
5	SDS 10%	100 ul
6	APS 10%	100 ul
7	TEMED	10ul


Sample: See protocoll

Result: One plate was succesfully made

📎 Screenshot_#4_18%_SDS-PAGE.png



Screenshot_#4_TLC_SDS-PAGE.png

**Platereader with lignin samples**

Experimenter: Johanna and Gabriel

Aim: To examine several different samples at the same time

Protocol:

Sample:

Result:

SUNDAY, 18/8/2019


Degradation of lignin with HRP with continuous addition of hydrogen peroxide

Experimenter: Johanna, Gabriel

Aim: To degrade lignin with HRP and measure changes in nanodrop by full spektrum analysis

Protocol: The reactions were performed in room temperature and the final volume was 10ml. Lignin (1 mg/ml) was solved in deionized water and stirred. HRP was added to a concentration of 0,03mg/ml or 0,003mg/ml, and hydrogen peroxide (90 μ M) was added every 120 minutes. These two samples were stirred for 4 hours*****

IMG_4549.jpg



Sample: Lignosulfanate lignin, Horse radish peroxidase (thermo fischer), hydrogen peroxide, deionized water.

Result: Peaks that showed a decrease of lignin on the nano-drop

MONDAY, 19/8/2019

SDS-PAGE #5

Experimenter: Johanna and Gabriel

Aim: To see a change of the different treated lignin running on the gel.

Protocol: 10 µl Sample were mixed with 10 µl loading dye and then loaded into the wells of the SDS-PAGE and then run from 80V and later increased to 120V

Sample: [1mg/ml lignin: 0,03_0h_HRP: 0,03_2h_HRP: 0,03_4h_HRP: sonicated sulfonated lignin 20mg/ml: sonicated MeOH lignin 20mg/ml: soniated acetone lignin 20mg/ml: sulfonated lignin 20mg/ml: Pellet_kill_it_4h: Kill_it_1_ : kill_it_2: kill_it_3:]

Result:

Degrading lignin with HRP

Experimenter: Johanna and Gabriel

Aim: to degrade lignin with enzyme

Protocol: : The reactions were performed in room temperature and the final volume was 10ml. Lignin (1 mg/ml) was solved in deionized water and stirred. HRP was added to a concentration of 0,03mg/ml or 0,003mg/ml, and hydrogen peroxide (90 µM) was added every 120 minutes. These two samples were stirred for 4 hours

Sample: Sulfonated lignin

Result:

Making of 50mM sodium acetate

Experimenter: Gabriel and Vamsi

Aim:

Protocol:

Sample:

Result:

TUESDAY, 20/8/2019

Degrading lignin with HRP

Experimenter: Johanna and Gabriel

Aim: to degrade lignin with HRP

Protocol: : The reactions were performed in room temperature and the final volume was 10ml. Lignin (1 mg/ml) was solved in deionized water and stirred. HRP was added to a concentration of 0,03mg/ml or 0,003mg/ml, and hydrogen peroxide (90 µM) was added every 120 minutes. These two samples were stirred for 4 hours

Sample: Sulfonated lignin

Result:

Bradford test

Experimenter: Johanna and Gabriel

Aim: To investigate what concentration the MnP from Sigma has

Protocol: BSA(Bovine Serum Albumin) was added in different concentration and checked with a spectrophotometer. The MnP was then tested with the same wavelength (595 nm), and later adjusted to the curve from the BSA

Sample: MnP from Sigma

Result: the enzyme concentration was around 0,005 mg/ml

run of SDS-PAGE #6

Experimenter: Johanna and Gabriel


Aim: to see a change of the degrading of lignin on SDS-PAGE

Protocol: 10 μ l Sample were mixed with 10 μ l loading dye and then loaded into the wells of the SDS-PAGE and then run from 80V and later increased to 120V


Sample: [20mg/ml sulfonated lignin: 20mg/ml sonicated lignin: 50 mg/ml sulfonated lignin: 50 mg/ml boiled sulfonated lignin: 20 mg/ml boiled sulfonated lignin: 20 mg/ml boiled sonicated lignin]

Result:

 Screenshot_#6_SDS-PAGE.png



 Screenshot_#6_SDS-PAGE_TLC.png



WEDNESDAY, 21/8/2019

Degrading lignin with HRP with continuous addition of hydrogen peroxide

Experimenter: Johanna and Johan

Aim: To degrade lignin with HRP and measure changes in nanodrop by full spektrum analysis

Protocol: The reactions were performed in room temperature and the final volume was 10ml. Lignin (1 mg/ml) was solved in deionized water and stirred. HRP was added to a concentration of 0,03mg/ml or 0,003mg/ml, and hydrogen peroxide 33 was added every 120 minutes. These two samples were stirred for 4 hours*****

Sample: Lignosulfanate lignin, Horse radish peroxidase (thermo fischer), hydrogen peroxide, deionized water.

Result:

THURSDAY, 22/8/2019

Degrading lignin with HRP

Experimenter: Johanna and Johan

Aim: To degrade lignin with HRP and H₂O₂ and to measure the changes

Protocol:

	A	B	C	D	E
1	Sample	Kraft Lignin	HRP	H2O2	ddH2O
2	A	0,02 g	-	33 ul	9,67 ml
3	B	0,02 g	250 ul	33 ul	9,717 ml
4	C	0,02 g	-	-	10 ml
5	D	0,02 g	250 ul	-	9,75 ml
6					

Sample: Kraft Lignin and sulfonated lignin

Result:

SDS-PAGE #7 from lignin degrading samples

Experimenter: Johanna and Johan


Aim: to see a change in the degraded lignin, on the gel.

Protocol: 10 µl Sample were mixed with 10 µl loading dye and then loaded into the wells of the SDS-PAGE and then run from 80V and later increased to 120V

Sample: [B_0h: B_2h: B_4h: B_16h: C : sulfonated lignin 20 mg/ml]

Result:

 Screenshot_#7_SDS-PAGE_TLC.png



Making 18% SDS_PAGE with 4% stacking gel

Experimenter: Johanna

Aim: To make a gel in order to run the lignin on it

Protocol:

Sample:

THURSDAY, 29/8/2019

Degrading Lignin with HRP and higher amounts of lignin

Experimenter: Gabriel and Johanna

Aim: Add HRP and hydrogen peroxide to break down lignin.

Protocol: The reactions were performed in room temperature and the final volume was 10ml. Lignin (Lineo prime, wet) (0,02 mg/ml) was solved in deionized water and stirred. HRP was added to a concentration of 0,03mg/ml or 0,003mg/ml, and hydrogen peroxide (90 μ M)***** was added at time 0. Samples were taken and put on ice at minute: 1, 5, 10, 15, 20, 25, 30, 40, 50, 60 and 90.

Sample:

Results:

Making 18% SDS_PAGE with 4% stacking gel

Experimenter: Gabriel

Aim: To make a gel in order to run the lignin on it

Protocol:

4% stacking gel			
	A	B	C
1	What	Quantity	
2	Water	6,1 ml	
3	Akrylamide/Bis 30%	1.3 ml	
4	Tris HCl 0,5M pH 6.8	2. 5 ml	
5	SDS 10%	100 ul	
6	APS 10%	100 ul	
7	Temed	10 ul	

18% resolving gel			
	A	B	C
1	What	Quantity	
2	Water	1,334	
3	Akrylamide/Bis 30%	6,066	
4	Tris HCl 1,5M pH 8,8	2. 5 ml	
5	SDS 10%	100 ul	
6	APS 10%	100 ul	
7	Temed	10 ul	

Sample:

Results:

Sonication of yeast pellet

Experimenter: Johanna, Gabriel

Aim: To lyse cells before testing samples on enzymatic assay.

Protocol: Cells were sonicated for 30 minutes.

Sample: Pellet of "X-33 HRP-2A-AAO pick XV VI (induced + uninduced)".

Results: Two falcon tubes with sonicated pellets suspension.

Gel with lignin

Experimenter: Johanna, Gabriel

Aim: To observe potential breakdown of lignin on SDS gel.

Protocol: 1, 5, 10, 15, 20, 25, 30, 40, 50, 60 and 90.

Sample:

Results:

Enzymatic assay on sonicated pellet of (AAO, induced+uninduced)

Experimenter: Johanna, Gabriel

Aim: To perform enzymatic assay with sonicated pellet of (X-33 AAO, induced+uninduced)*** using pierce quantitative peroxide assay kit.

Procedure: Followed the protocol provided in the kit (Thermo fisher - protocol 23280)

Induced and uninduced samples were tested, and each sample was tested with both reagent A which was made in the lab and reagent A that came with the kit. A total of four samples.

1mL of working reagent was added to 1,5mL cuvettes, the sample was added and mixed. After 15 minutes of incubation at room temperature, the OD was measured at 560nm.

Sample: Induced and uninduced pellet of "X-33 AAO"

Results: Sedimentation occurred in the induced samples when added to the working reagent

IMG_4792.JPG




Figure 1. Samples in working reagent after 15 minutes at room temperature. Order from left: reagent A from kit (Induced, uninduced, blank), self made reagent A (induced, uninduced, blank).

Table 1		A	B
1	Table 1. Absorbance measurements of induced and uninduced samples at 560nm		
2	Sample	abs (560nm)	
3	Original: induced	1,12	
4	Original: uninduced	1,356	
5	Self made: induced	0,773	
6	Self made: uninduced	1,828	

Enzymatic assay of Peroxidase using ABTS as substrate

Project: iGEM uppsala 2019

Authors: Gabriel Torres Muñoz

FRIDAY, 9/8/2019

Experimenter: Gabriel, Johanna

Aim: To test if the assay works for HRP by using HRP extracted from horseradish. Obtain the $\Delta A_{405 \text{ nm}}/\text{minute}$ for the maximum linear rate for both the Test and Blank to calculate Units/mg solid.

Protocol: Enzymatic Assay of Peroxidase (EC 1.11.1.7) 2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid) as a Substrate (CR SOP-DEK-ENZ42). Two different dilutions were tested: 0,005 and 0,0005.

Sample: HRP derived from horseradish, substrates needed according to instructions


Result: When using 0,005, the increase in absorbtion was beyond the maximum readings of the machines (>3). The maximum linear rate of the 0,0005 reaction was used to calculate the units/mg solid:***



image.png

$$\text{Units/mg solid} = \frac{(\Delta A_{405 \text{ nm}}/\text{min}(\text{test}) - \Delta A_{405 \text{ nm}}/\text{min}(\text{blank})) * 3.05 * \text{DF}}{36.8 * 0.05}$$

units/mg solid =

Table1				
	A	B	C	D
1	Table 1. Absorbance values obtained from two different dilutions of HRP and test with Bovine serum albumin			
2	Time (s)	Test (abs)	Dilute 0,005 (abs)	Dilute 0,0005 (abs)
3	0	-0,193	0,298	-0,19
4	10	-0,214	0,501	-0,18
5	20	-0,222	0,732	-0,171
6	30	-0,222	0,966	-0,173
7	40	-0,22	1,211	-0,166
8	50	-0,221	1,459	-0,151
9	60	-0,221	1,738	-0,134
10	70	-0,219	2,06	-0,117
11	80	-0,219	2,45	-0,098
12	90	-0,219	>3	-0,081
13	100	-0,219	>3	-0,064
14	110	-0,219	>3	-0,047
15	120	-0,218	>3	-0,028
16	130	-0,218	>3	-0,011
17	140	-0,218	>3	0,007
18	150	-0,218	>3	0,023
19	160	-0,217	>3	0,042
20	170	-0,217	>3	0,058
21	180	-0,217	>3	0,077





MONDAY, 12/8/2019

Experimenter: Gabriel, Johanna




Aim: To test again if the assay works for HRP by using HRP extracted from horseradish. Obtain the $\Delta A_{405\text{ nm}}/\text{minute}$ for the maximum linear rate for both the Test and Blank to calculate Units/mg solid.

Protocol: Enzymatic Assay of Peroxidase (EC 1.11.1.7) 2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid) as a Substrate (CR SOP-DEK-ENZ42). A dilution of 0,0005 was used.

Sample: HRP derived from horseradish. Substrates used according to instructions.

Result: Maximum linear rate for 0,0005 was 0,007, and for the test with ovine serum albumin it was 0,0001. In this test the absorbance values were positive as expected, compared to earlier test.

The maximum linear rate of the 0,0005 reaction was used to calculate the units/mg solid ***



TUESDAY, 13/8/2019

Experimenter: Gabriel, Johanna

Aim: To observe the activity of the HRP obtained from Thermo Fischer

10/20/2019

Enzymatic assay of Peroxidase using ABTS as substrate · Benchling

assay works for HRP by using HRP extracted from horseradish. Obtain the $\Delta A_{405\text{ nm}}/\text{minute}$ for the maximum linear rate for both the Test and Blank to calculate Units/mg solid.

Protocol:

Sample:

Results: