



DYS SEE

PCR - Q5.



Protocols

PCR protocol for NEB's Q5 High-Fidelity 2X Master Mix

Reaction setup according to NEB

Component	25 μ l Reaction	50 μ l Reaction	Final Concentration
Q5 High-Fidelity 2X Master Mix	12.5 μ l	25 μ l	1X
10 μ M Forward Primer	1.25 μ l	2.5 μ l	0.5 μ M
10 μ M Reverse Primer	1.25 μ l	2.5 μ l	0.5 μ M
Template DNA	Variable	variable	< 1,000 ng
ddH ₂ O	to 25 μ l	to 50 μ l	

- Prepare all reactions on ice.
- Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.
- Transfer PCR tubes to a PCR machine and begin thermocycling.

Thermocycling Conditions for a Routine PCR

STEP	TEMP	TIME
Initial Denaturation	98°C	30 seconds
35 Cycles	98°C	10 seconds
	*50–72°C	15 seconds
	72°C	10 seconds/kb
Final Extension	72°C	2 minutes
Hold	4–10°C	



Protocols