

# 10. (October) 2019

---

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Johanna Wallner

THURSDAY, 17/10/2019

---

## Johanna + Sarah

### Purification

- samples:
  - new Coiled-Coil L7Ae
  - new Coiled-Coil MCP
  - old Coiled-Coil MCP
  - Gag-L7Ae Fusionconstruct (purification with 10 CV highest salt concentration at the end)
- **Preparation**
  - Prepare Stock Solutions, filter them!
  - Prepare Buffer A, Buffer B, Equilibration Buffer, Wash Buffer
  - Prepare Elution Buffer (20 x 10 mL -> label 20 x 15 mL Falcons)
  - Prepare XY x FT Eppis, 10 x Wash Eppis, 20 x Elution Eppis
  - Equilibrate the heparin resin with 5 CV **Equilibration Buffer**
  - Prepare the sample: harvest the supernatant from the 6-well plate. Centrifuge 10 minutes at 2000 g.
  - Filtration: Filter the harvested supernatant through a 0.8 µm syringe filter
  - Load the sample on the resin -> collect Flow-through = XY mL = XY fractions
- **Wash**
  - Wash with 10 CV Wash Buffer = 10 fractions
- **Elution**
  - Elute with a salt gradient (120 - 20000 mM NaCl) with 20 CV = 20 fractions

Salt Gradient					
	NaCl [mM]	Buffer A [ $\mu$ L]	Buffer B [ $\mu$ L]	Buffer A [mL]	Buffer B [mL]
1	150	925	75	9.25	0.75
2	200	900	100	9	1
3	300	850	150	8.5	1.5
4	400	800	200	8	2
5	500	750	250	7.5	2.5
6	600	700	300	7	3
7	700	650	350	6.5	3.5
8	800	600	400	6	4
9	900	550	450	5.5	4.5
10	1000	500	500	5	5
11	1100	450	550	4.5	5.5
12	1200	400	600	4	6
13	1300	350	650	3.5	6.5
14	1400	300	700	3	7
15	1500	250	750	2.5	7.5
16	1600	200	800	2	8
17	1700	150	850	1.5	8.5
18	1800	100	900	1	9
19	1900	50	950	0.5	9.5
20	2000	0	1000	0	10

10 CV at the end of the gradient with **100 % Buffer B**

#### Regeneration

10 CV 1 M NaOH

10 CV 30 % 2-Propanol

- total elution was about 50-60 % from total load
- after purification the samples were concentrated by using 30 kDa centricons -> the concentrated samples are used for Western Blot and DLS