

BG-11 medium preparation Protocol and Growth Conditions

Introduction

To grow *Synechococcus elongatus* sp. *UTEX 2973*, BG-11 was used. This protocol originates from The University of Texas (2009). BG-11 Medium Recipe [on line]. Available at <https://utex.org/products/bg-11-medium> (15/05/2019)

Materials

- › Stock 1 (1L)
 - › Na₂MG EDTA 0.1g/liter
 - › Ferric ammonium citrate 0.6g
 - › Citric acid.1H₂O 0.6g
 - › CaCl₂.2H₂O 3.6g
- › Stock 2 (1L)
 - › MgSO₄.H₂O 7.5g
- › Stock 3 (1L)
 - › K₂HPO₄ 3.05g
- › Stock 4:Trace metal solution (1L)
 - › H₃BO₃ 2.86g
 - › MnCl₂.4H₂O 1.81g
 - › ZnSO₄ . 7H₂O 0.222g
 - › CuSO₄ . 5H₂O 0.079g
 - › COCl₂ . 6H₂O 0.050g
 - › Na₂MoO₄ . 2H₂O 0.391g
 - › or MoO₄ (85%) 0.018g

Procedure

BG-11 Medium

1. Prepare stock 1-4: in 900 ml dH₂O add components in the order specified on a magnetic stirrer. Bring the total volume to 1L with dH₂O. Autoclave at 15 psi for 30 minutes at 121°C.
2. Combine the following solutions and adjust pH to 7,5 (Use 1.0M HCl or NaOH).

Stock	Per liter of medium
Stock 1	10 ml
Stock 2	10 ml
Stock 3	10 ml
Na ₂ CO ₃	0,02 g
Stock 4	1,0 ml

NaNO₃ 1.5 g (For 100x BG-11 solution don't add NaNO₃)

For **solid media**: add 1.5% agar to the liquid BG11 media.

3. Aliquot into flasks (50 ml/125 ml flask) with cotton stoppers on top and autoclave at 15 psi for 30 minutes at 121°C.
4. After autoclaving and cooling, the pH might have changed, so check and adjust.

Growth Conditions

5. The day of arrival, the culture was first put under light (50 $\mu\text{mol}/\text{m}^2\cdot\text{s}$) for 20 hours, and further incubated in dark for 24 hours to avoid bleaching. Afterwards, cultures were inoculated in 30 mL of BG-11 medium and grown in conditions as described below.
6. *S. elongatus* sp. *UTEX 2973* was grown in filtertop erlenmeyers (0.22 μm) in 28°C at 100 rpm and a constant irradiation of 50 $\mu\text{mol}/\text{m}^2\cdot\text{s}$.