

Flow Cytometry Protocol

E.coli containing the backbone(AraC3), C protein characterization construct (Lysogenic module: [BBa_K3024004](#)), Inducible Cox characterization construct (Lytic module: [BBa_K3024011](#)), Switch Plasmid (http://parts.igem.org/Part:BBa_K3024015), fluorescence reporters GFP and RFP were grown overnight in 37°C in TSB and induced for 5 hours with Arabinose (0.2%).

Step 1: Preparation

50 µl bacteria were taken out of each flask after 5 hours of induction with Arabinose (0.2%)

The samples were centrifuged at 13 000 rpm for 3 minutes.

The supernatant of each sample was discarded and the cells were resuspended in 750 µl Phosphate-Buffered Saline (PBS).

Step 2: Analysis

The Flow Cytometry machine (Gallios Flow Cytometer from Beckman Coulter) was set up with gating systems for detection of RFP and GFP in *E.coli* and cleaned according to protocol.

The data was collected and analyzed in the Flow Cytometry analysis software Kaluza.