

Name: Chiara, Kennex, Rehmat, Amirah, Shakera, Xuecheng Ye Yilin Lu Zexi Guo Jiayi Lan

Yujie Huang

Date: 8/21/19

Goal:

1. Transform DinolII-GFP into *S. Microadriaticum* and *O. Marina*
2. Finish DinolII-RFP ethanol precipitation
3. New overnights of mcherry from glycerol stocks
4. Measure the ABS600 and fluorescence of the overnight culture Interlab

Name: Chiara

Date: 8/21/19

Goal:

1. Transform DinolII-GFP into S. Microadriaticum and O. Marina

Protocol:

1. Cultured *Symbiodinium Microadriaticum* cells in ASP-8A with filtered seawater medium under natural sunlight for about 2 months and *O. Marina* cells in f/2 media for about 2 months with no antibiotics

- Cultures:
  - S. Microadriaticum: ASP 8a-SW media 75 mL labeled 8/21/19; [cell]=  
 $2.52 \times 10^6$  cells/mL
  - ii. O. Marina: f/2 , FSW media 75 mL labeled 6/28/19; [cell]=  $3.8 \times 10^5$  cells/mL

2. Harvested the cells by centrifugation at 800 g for 5 min at 4°C.

- a. S. Microadriaticum: Pelleted 5 mL=  $12.6 \times 10^6$  cells
- b. O. Marina: Pelleted 1 mL=  $3.8 \times 10^5$  cells

3. Used 500 µL of 0.1M EDTA to resuspend the cell pellet

4. Centrifuged at 800 g for 2 min at 4°C.

5. Washed cells with 10% Glycerol 3 times, centrifuged at 800 g for 2 min in 4°C.

6. Re-measured the cell concentration of each sample

- a. [S. Microadriaticum]=  $5.0 \times 10^6$  cells/mL
- b. [O. Marina]= 0 cells":s\[PTRSGO]

i. O. Marina seem to have either died and bursted or were lost in a step.  
Therefore, we discontinued the protocol for this sample

7. Resuspended the pellet in 1 mL of 10% Glycerol1

8. Incubated ~100µl of cells with 17 µl (~1µg) of DinolIII-GFP DNA (60 ng/µL 8/16/19 ethanol precipitation sample) on ice for 5 min.

9. Put cells into a 0.2 cm cuvette, and electroporated using SHS (2.0 kV, 1 pulse), SC2 (1.5 kV, 1 pulse), or DIC (1.0 kV, 2 pulses, 1.0 msec) program with Bio-Rad MicroPulser 165-2100.

10. Added 1mL of ASP-8A SW medium to the 0.2 cm cuvette, mixed well, and transferred to a 15 mL tube

11. Left the cultures to grow under the programmed lights

12. Observed the cells under a normal microscope in 1-3 days and according to the need.

Name: Shakera

Date: 8/21/19

Goal:

1. New overnights of mcherry from glycerol stocks

Protocol:

1. Scraped ice from mcherry glycerol stocks and placed in a solution consisting of 7  $\mu$ L LB and 7  $\mu$ L ampicillin
2. Let incubate overnight at 37 degrees C

Name: Xuecheng Ye Yilin Lu Zexi Guo Jiayi Lan Yujie Huang

Goal:

1. Measure the ABS600 and fluorescence of the overnight culture

Results:

<u>LB</u>	<u>positive control</u>	<u>control</u>	<u>device 1</u>	<u>device 5</u>
<u>colony1</u>	<u>0.149</u>	<u>0.165</u>	<u>0.157</u>	<u>0.170</u>
	<u>0.158</u>	<u>0.172</u>	<u>0.163</u>	<u>0.175</u>
	<u>0.158</u>	<u>0.174</u>	<u>0.170</u>	<u>0.173</u>
	<u>0.158</u>	<u>0.179</u>	<u>0.163</u>	<u>0.171</u>
<u>colony2</u>	<u>0.198</u>	<u>0.197</u>	<u>0.207</u>	<u>0.206</u>
	<u>0.210</u>	<u>0.210</u>	<u>0.213</u>	<u>0.209</u>
	<u>0.210</u>	<u>0.207</u>	<u>0.208</u>	<u>0.212</u>
	<u>0.214</u>	<u>0.207</u>	<u>0.206</u>	<u>0.209</u>

Table 1: The ABS600 of the overnight culture in LB medium

	<u>positive conyrol</u>	<u>negetive control</u>	<u>device1</u>	<u>device 5</u>	<u>YM mediuum</u>

	<u>0.111</u>	<u>0.109</u>	<u>0.149</u>	<u>0.141</u>	<u>0.041</u>
<u>colony</u> <u>1</u>	<u>0.099</u>	<u>0.111</u>	<u>0.149</u>	<u>0.142</u>	<u>0.045</u>
	<u>0.108</u>	<u>0.114</u>	<u>0.142</u>	<u>0.145</u>	<u>0.043</u>
	<u>0.099</u>	<u>0.115</u>	<u>0.171</u>	<u>0.143</u>	<u>0.045</u>
<u>colony</u> <u>2</u>	<u>0.128</u>	<u>0.145</u>	<u>0.138</u>	<u>0.143</u>	<u>0.045</u>
	<u>0.130</u>	<u>0.138</u>	<u>0.132</u>	<u>0.146</u>	<u>0.047</u>
	<u>0.134</u>	<u>0.143</u>	<u>0.131</u>	<u>0.141</u>	<u>0.046</u>
	<u>0.133</u>	<u>0.142</u>	<u>0.130</u>	<u>0.155</u>	<u>0.052</u>

Table2: the Abs600 of overnight culture in YM medium

	<u>positiv e</u>	<u>negetiv e</u>	<u>device 1</u>	<u>device 5</u>	<u>YPD medium</u>
<u>colony</u> <u>1</u>	<u>0.206</u>	<u>0.192</u>	<u>0.155</u>	<u>0.200</u>	<u>0.053</u>
	<u>0.207</u>	<u>0.184</u>	<u>0.159</u>	<u>0.202</u>	<u>0.052</u>
	<u>0.200</u>	<u>0.193</u>	<u>0.147</u>	<u>0.198</u>	<u>0.052</u>
	<u>0.197</u>	<u>0.181</u>	<u>0.142</u>	<u>0.190</u>	<u>0.053</u>
<u>colony</u> <u>2</u>	<u>0.206</u>	<u>0.214</u>	<u>0.202</u>	<u>0.189</u>	<u>0.051</u>
	<u>0.207</u>	<u>0.214</u>	<u>0.200</u>	<u>0.194</u>	<u>0.051</u>
	<u>0.200</u>	<u>0.203</u>	<u>0.198</u>	<u>0.186</u>	<u>0.081</u>
	<u>0.197</u>	<u>0.204</u>	<u>0.185</u>	<u>0.189</u>	<u>0.052</u>

Table 3: the Abs600 of overnight culture in YPD medium

—	<u>positiv e</u>	<u>negetiv e</u>	<u>device 1</u>	<u>device 5</u>	<u>TB medium</u>

	<u>0.274</u>	<u>0.294</u>	<u>0.238</u>	<u>0.302</u>	<u>0.059</u>
<u>colony</u> <u>1</u>	<u>0.294</u>	<u>0.304</u>	<u>0.246</u>	<u>0.292</u>	<u>0.059</u>
	<u>0.302</u>	<u>0.305</u>	<u>0.248</u>	<u>0.309</u>	<u>0.058</u>
	<u>0.277</u>	<u>0.286</u>	<u>0.241</u>	<u>0.302</u>	<u>0.056</u>
	<u>0.246</u>	<u>0.252</u>	<u>0.225</u>	<u>0.242</u>	<u>0.058</u>
<u>colony</u> <u>2</u>	<u>0.242</u>	<u>0.266</u>	<u>0.236</u>	<u>0.246</u>	<u>0.059</u>
	<u>0.241</u>	<u>0.257</u>	<u>0.232</u>	<u>0.258</u>	<u>0.058</u>
	<u>0.252</u>	<u>0.256</u>	<u>0.223</u>	<u>0.241</u>	<u>0.059</u>

Table4 : the Abs600 of overnight culture in TB medium

—	<u>positive</u>	<u>negetive</u>	<u>device1</u>	<u>device5</u>
<u>colony</u> <u>1</u>	<u>23122</u>	<u>17274</u>	<u>18319</u>	<u>18383</u>
	<u>23330</u>	<u>17980</u>	<u>17996</u>	<u>19002</u>
	<u>23416</u>	<u>17618</u>	<u>18460</u>	<u>19288</u>
	<u>22707</u>	<u>17638</u>	<u>18186</u>	<u>19015</u>
<u>colony</u> <u>2</u>	<u>27916</u>	<u>17312</u>	<u>19609</u>	<u>18248</u>
	<u>28484</u>	<u>18485</u>	<u>18452</u>	<u>18774</u>
	<u>29335</u>	<u>17821</u>	<u>18420</u>	<u>18806</u>
	<u>29018</u>	<u>18313</u>	<u>18564</u>	<u>18542</u>

Table5 : the fluorescence of overnight culture in LB medium

—	<u>positiv</u> <u>e</u>	<u>negetiv</u> <u>e</u>	<u>device</u> <u>1</u>	<u>device</u> <u>5</u>	<u>YM</u> <u>medium</u>
<u>colony</u>	<u>6881</u>	<u>6551</u>	<u>10220</u>	<u>6544</u>	<u>6380</u>

<u>1</u>	<u>6992</u>	<u>6459</u>	<u>10057</u>	<u>6478</u>	<u>6527</u>
	<u>6818</u>	<u>6930</u>	<u>9351</u>	<u>6678</u>	<u>6520</u>
	<u>6893</u>	<u>6528</u>	<u>10984</u>	<u>6499</u>	<u>6638</u>
<u>colony</u> <u>2</u>	<u>6551</u>	<u>6433</u>	<u>9372</u>	<u>10138</u>	<u>6619</u>
	<u>6738</u>	<u>6194</u>	<u>9263</u>	<u>10907</u>	<u>6636</u>
	<u>6395</u>	<u>6367</u>	<u>8891</u>	<u>10744</u>	<u>6525</u>
	<u>6823</u>	<u>6638</u>	<u>8909</u>	<u>11096</u>	<u>6850</u>

Table6 : the fluorescence of overnight culture in YM medium

<u>—</u>	<u>positiv e</u>	<u>negetiv e</u>	<u>device 1</u>	<u>device 5</u>	<u>YPD medium</u>
<u>colony</u> <u>1</u>	<u>19512</u>	<u>19361</u>	<u>19857</u>	<u>19447</u>	<u>19807</u>
	<u>19522</u>	<u>18431</u>	<u>20573</u>	<u>19712</u>	<u>19358</u>
	<u>19696</u>	<u>19178</u>	<u>19152</u>	<u>19675</u>	<u>18729</u>
	<u>19051</u>	<u>18257</u>	<u>19330</u>	<u>18473</u>	<u>18371</u>
<u>colony</u> <u>2</u>	<u>19317</u>	<u>19041</u>	<u>20068</u>	<u>19763</u>	<u>19247</u>
	<u>19514</u>	<u>19347</u>	<u>20403</u>	<u>19822</u>	<u>19762</u>
	<u>19588</u>	<u>18751</u>	<u>20459</u>	<u>19214</u>	<u>38173</u>
	<u>18614</u>	<u>18552</u>	<u>19638</u>	<u>18759</u>	<u>19743</u>

Table7 : the fluorescence of overnight culture in YPD medium

<u>—</u>	<u>positiv e</u>	<u>negetiv e</u>	<u>device 1</u>	<u>device 5</u>	<u>TB medium</u>
<u>colony</u> <u>1</u>	<u>27036</u>	<u>25805</u>	<u>39170</u>	<u>26497</u>	<u>24979</u>

	<u>27089</u>	<u>26100</u>	<u>39941</u>	<u>24860</u>	<u>24569</u>
	<u>27210</u>	<u>25849</u>	<u>39884</u>	<u>25485</u>	<u>23808</u>
	<u>25944</u>	<u>24820</u>	<u>39088</u>	<u>24538</u>	<u>23824</u>
<u>colony</u> <u>2</u>	<u>26819</u>	<u>25522</u>	<u>40757</u>	<u>37385</u>	<u>25328</u>
	<u>26155</u>	<u>25852</u>	<u>40950</u>	<u>37091</u>	<u>24987</u>
	<u>26068</u>	<u>25564</u>	<u>40793</u>	<u>37177</u>	<u>24399</u>
	<u>26130</u>	<u>25114</u>	<u>39834</u>	<u>35650</u>	<u>24819</u>

Table8 : the fluorescence of overnight culture in TB medium

Conclusion:

It seems device 5 colony 1 didn't grow well in any kind of medium. In comparison, the TB medium is better for the culture.