07. (July) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

FRIDAY, 26/7/2019

Miniprep: V31 & V32 (eluted in 20µL MilliQ)

-> V31

#1:80.5ng/μL#8:112.1ng/μL#12:101.0ng/μL

-> V32

#1:79.4ng/μL#4:77.0ng/μL#11:77.7ng/μL

• sequencing each with 15µL sample

Cell culture :

- splitting of johannas mother plate (passage 16), 600µL for inoculation
- splitting of Alejandros mother plate (passage 16), 800µL for inoculation

Cell culture:

- Medium exchange in transfected wells for HiBit assay (plates 1 & 2, conditions 1 to 10)
 - -> 600 of the 800µL were exchanged after ~48h

Maxi Prep inoculation:

• V15#12 in 150µL LB + Amp

Fluorescence microscopy assay:

- plate 3 was used
- The medium over the cells was removed and replaced with 500μL Medium containing a 1μL/mL solution of Alexa Fluor 546 biocytin (1:1000 of stock solution)
- the cells were incubated at 37°C for 35min
- the medium with the fluophore was removed and cells were washed with 37°C warm PBS (600μL)
- cells were imaged at Evos microscope
 - -> V27 appears to allocate on the cellmembrane

file://tmp/tmpt273U7.html