

# 08. (August) 2019

---

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Anja

MONDAY, 26/8/2019

---

**Joshi**

**Cell culture: HEK293T splitting and seeding**

- 6 wells, 760.000 cells/well, 2.5ml full DMEM/well
- 400ul splitted

**Cell culture: Harvesting transfected flask**

- 13ml taken, centrifuged on 5.5rpm, 10min
- 11ml taken + 5.5ml Total Exosome Isolation Reagent
- Mix & incubate ON

**MidiPrep** inoculation (E16, 150ml LB+Amp)

**Anja**

**Cell culture: HEK293T seeding**

- For cytotoxicity assay (96-well plate, 25.000 cells/well, 150 ul full DMEM/well), H row only DMEM for background controls