

## **SDS-PAGE**

**Connor Trotter – 01/08/2019**

### **Materials**

Protein or cell sample

1.5 mL Eppendorf tubes

5x SDS-loading buffer (63  $\mu$ L 1M Tris-HCL pH 6.9; 240  $\mu$ L 100% glycerol, 687  $\mu$ L ddH<sub>2</sub>O, 10  $\mu$ L bromophenol blue, 500  $\mu$ L 10% SDS solution)

SDS-PAGE gel

Gel percentage differs between target protein size

Buffer dam required if only 1 gel

SDS-PAGE gel tank

1x SDS-Running buffer

Coomassie Brilliant Blue

Destain buffer (20% v/v ethanol, 10% v/v acetic acid)

### **Methods**

1. Resuspend cell pellets in 50  $\mu$ L 5x SDS-loading buffer or add 1:1 ratio for protein samples (eg from purification) in a 1.5 mL Eppendorf tube
2. Boil samples at 100 °C for 10 minutes, opening every 2/3 minutes and vortexing to avoid pressure build up
3. Load the SDS-PAGE gel into the SDS-PAGE tank
  - a. Add fresh 1x SDS-running buffer to the inner gasket; recycled 1x SDS-running buffer can be used for the external basket if visibly free from contamination
4. Load 5-7  $\mu$ L of protein ladder to the first well
  - a. Volume depends on size of wells
5. Load 10-15  $\mu$ L of sample to each well
6. Plug in the gel and run at 180 V for until blue dye front is at the bottom of the gel
  - a. Typically >40 minutes
7. Drain the tank and obtain the gel
8. Stain gel in Coomassie Brilliant Blue for 2/3 hours

9. Destain the gel in destain buffer overnight
10. Image the gel using appropriate gel dock equipmen