



# Lethbridge HS iGEM 2019

## Batch Affinity Purification of Cas13a

Protocol obtained from the Kothe lab. Buffer composition compared to [Gootenberg et al., 2017](#)

### **Buffers**

Buffer A:

20 mM Tris-HCl (pH 8.0)

400 mM KCl

5% v/v glycerol

30 mM Imidazole

1 mM  $\beta$ ME

0.5 mM PMSF

Buffer B:

20 mM Tris-HCl (pH 8.0)

400 mM KCl

5% v/v glycerol

500 mM Imidazole

1 mM  $\beta$ ME

S200 Buffer:

10 mM HEPES (pH 7.5)

1 M NaCl

5 mM  $\text{MgCl}_2$

2 mM DTT



# Lethbridge HS iGEM 2019

## **Preparation of Ni Sepharose Slurry**

- 1) Gently shake resin bottle to create a homogenous medium and remove 5 mL slurry using a serological pipette. Transfer slurry to a 50 mL falcon tube, label with Ni-Sepharose, date, initials, and the protein
- 2) Spin the Ni-NTA sepharose column at 500 xg for 2 min to sediment the resin. Remove the supernatant.
- 3) Wash the resin with 3 volumes of sterile water. Gently mix to create a slurry. Spin at 500 xg for 5 min.
- 4) Wash the resin 6 times with Buffer A, leave 3 mL of Buffer A on the resin to obtain a 50% slurry.
- 5) Divide resin into more falcon tubes if lysate volume is greater than 50 mL.

## **Cell opening**

- 1) Add PMSF and  $\beta$ ME to buffer A
- 2) Resuspend the frozen cell pellet in buffer A (~5 mL/g of cell) in a small beaker. Thaw cells by stirring mixture slowly on ice (30-60 min).
- 3) Add lysozyme (1 mg/mL) and incubate the cell suspension on ice for 30 min
- 4) Add sodium deoxycholate (12.5 mg/g of cell). Continue stirring on ice for 30 min
- 5) Open cells using a sonicator while sample is on ice: at least 10 mins for 1 min, short pause in-between, shake beaker, intensity level 6, duty cycle at 60%
- 6) Centrifuge cell lysate for 30 mins at 30000 xg at 4 °C in a JA-25.50 rotor.
- 7) Take an  $A_{280}$  reading of cell lysate to monitor cell opening
- 8) Remove 50  $\mu$ L of cell lysate for SDS-PAGE analysis
- 9) Store pellet for later analysis

## **Purification by Ni-Sepharose slurry**



# Lethbridge HS iGEM 2019

- 1) Apply clarified cell lysate to washed Ni-sepharose evenly and gently mix thoroughly. Incubate for 60 mins inverting periodically to bind protein to resin on ice.

NOTE: all further steps are performed using cold buffers (4 °C).

- 2) Spin slurry at 500 xg for 5 min, remove the supernatant and store at 4 °C. Take a 50 µL sample for SDS PAGE analysis and store at 4 °C
- 3) Wash the resin 6 times with Buffer A, spinning at 500 xg for 5 min each time. Pool the washes, store at 4 °C. Take a 50 µL sample for SDS PAGE analysis and store at 4 °C.
- 4) Add βME to Buffer B.
- 5) Elute the protein 6 times using 90% resin volume of Buffer B. Incubate with buffer B for 5 min on ice. Spin for 5 min at 500 xg. Save elutions in 15 mL falcon tubes, store at 4 °C, and take a 50 µL sample for SDS PAGE analysis and store at 4 °C.
- 6) Run an SDS-PAGE of all samples at 180 V for ~3 hours.

## **Regeneration of Ni-sepharose**

- 1) Wash resin 5 more times with 20 mL buffer B, centrifuge for 2 min at 500 xg. Check the last wash on a gel to check for protein presence.
- 2) Wash resin 3 times with 40 mL of d<sub>2</sub>H<sub>2</sub>O.
- 3) Wash resin once with 40 mL 20% ethanol. Leave 1.5 mL 20% ethanol on resin to obtain a 50% slurry for storage at 4 °C.

## **Concentration**

- 1) Rinse Vivaspin MWCO 30000 with S200 buffer (2 mL) and centrifuge at 4000 xg for 10 min
- 2) Pool elutions together.
- 3) Remove buffer from Vivaspin and add pooled elutions.
- 4) Centrifuge at 4000 xg for 5 min
- 5) Keep centrifuging: check for speed of concentrating. Concentrate to 5 mL. Watch for precipitate!



# Lethbridge HS iGEM 2019

- 6) Remove filtrate from bottom of Vivaspinn, and pool with other filtrate and store at 4 °C.
- 7) Store concentrated sample at 4 °C if not proceeding directly to size exclusion chromatography.