

PureLink™ RNA Mini Kit



Required Materials

- 96–100% ethanol, 2–mercaptoethanol, 70% ethanol (in RNase-Free Water), 1.5 mL RNase-free microcentrifuge tubes
- Homogenizer, RNase-free syringe (1 mL) with 18–21-gauge needle or, Rotor-stator homogenizer Microcentrifuge capable of centrifuging $12,000 \times g$
- PBS (for samples with $>10^7$ cells)
- 15 mL RNase-free tubes (for samples with $>10^7$ cells), RNase-free pipet tips

Buffer Preparation

When using Wash Buffer II for the first time, add 60 mL 96–100% ethanol or 300 mL 96–100% ethanol. Mark the label to indicate that ethanol is already added.

Prepare fresh Lysis Buffer containing 1% 2-mercaptoethanol. Add 10 μ L 2-mercaptoethanol for every 1 mL Lysis Buffer.

Cell Number	Lysis Buffer Required for Each Sample
$\leq 1 \times 10^6$	0.3 ml (0.6 mL if using a rotor-stator for lysis/homogenization)
1×10^6 – 5×10^6	0.6 mL
5×10^6 – 5×10^7	0.6 mL per 5×10^6 cells (e.g., use 1.2 mL for 1×10^7 cells)

Lysis and Homogenization

$\leq 5 \times 10^6$ Suspension Cells

1. Transfer the cells to an RNase-free tube and centrifuge at $2,000 \times g$ for 5 min at 4°C to pellet. Discard the growth medium.
2. Add 0.3 or 0.6 mL Lysis Buffer with 2-mercaptoethanol to the sample (see table above for volume).
3. Vortex until the cell pellet is dispersed and the cells appear lysed.
4. Proceed to Homogenization below.

$\leq 5 \times 10^6$ Monolayer Cells

1. Remove the growth medium from the cells, then add 0.3 or 0.6 mL Lysis Buffer with 2-mercaptoethanol (see table above for volume).
2. Vortex until the cell pellet is dispersed and the cells appear lysed.
3. Proceed to Homogenization below.

5×10^6 - 5×10^7 Suspension Cells

1. Transfer cells to a 15-mL tube and centrifuge at $2,000 \times g$ for 5 min at 4°C . Discard the supernatant.
2. Add 0.6 mL Lysis Buffer with 2-mercaptoethanol (see table above for volume).
3. Vortex until the cell pellet is dispersed and the cells appear lysed.
4. Homogenize at room temperature with a rotor-stator homogenizer (see Homogenization below).

Frozen Cell Pellets

1. Transfer cells to a 15-mL tube and add 0.6 mL Lysis Buffer with 2-mercaptoethanol (see table above for volume).
2. Vortex until the cell pellet is dispersed and the cells appear lysed.
3. Homogenize at room temperature with a rotor-stator homogenizer (see Homogenization below).

Homogenization

1. Proceed with one of the following homogenization options at room temperature:
 - Transfer the lysate into a clean homogenization tube, and perform manual homogenization. Centrifuge the homogenate at $12,000 \times g$ for 2 minutes.
 - Pass the lysate 5–10 times through an 18- to 21-gauge syringe needle.
 - Transfer the lysate into a clean tube, and homogenize using a rotor-stator homogenizer at maximum speed for ≥ 45 s. Centrifuge the homogenate at $26,000 \times g$ for 5 minutes, then transfer the supernatant to a clean RNase-free tube.
2. Proceed to RNA Purification.

RNA Purification

Binding, Washing, and Elution of RNA

1. Add one volume 70% ethanol to each volume of cell homogenate.
2. Vortex to mix thoroughly and to disperse any visible precipitate that may form after adding ethanol.
3. Transfer up to 700 μ L of the sample (including any remaining precipitate) to the spin cartridge (with the collection tube).
4. Centrifuge at $12,000 \times g$ for 15 seconds at room temperature. Discard the flow-through, and reinsert the spin cartridge into the same collection tube.
5. Repeat Steps 3–4 until the entire sample has been processed.
6. Add 700 μ L Wash Buffer I to the spin cartridge.
7. Centrifuge at $12,000 \times g$ for 15 seconds at room temperature. Discard the flow-through and the collection tube. Place the spin cartridge into a new collection tube.
8. Add 500 μ L Wash Buffer II with ethanol to the spin cartridge.
9. Centrifuge at $12,000 \times g$ for 15 seconds at room temperature. Discard the flow-through.
10. Repeat Steps 8–9 once.
11. Centrifuge the spin cartridge at $12,000 \times g$ for 1–2 minutes to dry the membrane with bound RNA. Discard the collection tube and insert the spin cartridge into a recovery tube.
12. Add 30–100 μ L RNase-free water to the center of the spin cartridge.
13. Incubate at room temperature for 1 minute.
14. Centrifuge the spin cartridge for 2 minutes at $\geq 12,000 \times g$ at room temperature to elute the RNA from the membrane into the recovery tube. Note: If the expected RNA yield is $>100 \mu$ g, perform 3 sequential elutions of 100 μ L each. Collect the eluates in a single tube.
15. Store your purified RNA or proceed to downstream application.

Storage

Store the purified RNA on ice for immediate use.

For long-term storage, keep the purified RNA at -80°C. Perform DNase I treatment after purification (refer to the PureLink® RNA Mini Kit manual) to assure highly pure RNA without genomic DNA contamination.

Source: ThermoFisher Scientific