

# 10. (October) 2019

Project: iGEM\_Munich2019 Shared Project

Authors: Johanna Wallner

SATURDAY, 5/10/2019

Alejandro  
Cell Culture: Transfection

- 96-well plate
  - transfect triplicates

Transfectionmix 96-well plate 05/1...		
	A	B
1	DNA per well	100 ng
2	P3000 Reagent per well	0.2 µL
3	Lipofectamine 3000 reagent per well	0.15 µL
4	OptiMEM per well	2 x 5 µL

- Transfect cells according to the following table. Use the indicated volume of DNA and P3000™ Reagent with each of the two volumes of Lipofectamine™ 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.

Timeline		Procedure Details (Two Reaction Optimizations)			
Day 0	Steps	Component	96-well	24-well	6-well
1	Seed cells to be 70-90% confluent at transfection	Adherent cells	1-4 x 10 <sup>4</sup>	0.5-2 x 10 <sup>5</sup>	0.25-1 x 10 <sup>6</sup>
		Opti-MEM™ Medium	5 µL x 2	25 µL x 2	125 µL x 2
		Lipofectamine™ 3000 Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
2	Dilute Lipofectamine™ 3000 Reagent in Opti-MEM™ Medium (2 tubes) - Mix well	Opti-MEM™ Medium	10 µL	50 µL	250 µL
		DNA (0.5-5 µg/µL)	0.2 µg	1 µg	5 µg
		P3000™ Reagent (2 µL/µg DNA)	0.4 µL	2 µL	10 µL
3	Prepare master mix of DNA by diluting DNA in Opti-MEM™ Medium, then add P3000™ Reagent - Mix well	Diluted DNA (with P3000™ Reagent)	5 µL	25 µL	125 µL
		Diluted Lipofectamine™ 3000 Reagent	5 µL	25 µL	125 µL
		Incubate for 10-15 minutes at room temperature.			
4	Add Diluted DNA to each tube of Diluted Lipofectamine™ 3000 Reagent (1:1 ratio)	Component (per well)	96-well	24-well	6-well
		DNA-lipid complex	10 µL	50 µL	250 µL
		DNA amount	100 ng	500 ng	2500 ng
5	Incubate	P3000™ Reagent	0.2 µL	1 µL	5 µL
		Lipofectamine™ 3000 Reagent used	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
		Incubate cells for 2-4 days at 37°C. Then, analyze transfected cells.			
6	Add DNA-lipid complex to cells				
7	Visualize/analyze transfected cells				

- Transfection scheme

Transfection scheme - 05/10/19 in ng per well - 96-well plate- HiBit									
	condition	wells	V8	V41	V42.1	V42.2	V15	V34	
1	1	6	-	-	-	-	-	-	1
2	2	6	60 ng	-	-	-	40 ng	-	-
3	3	6	-	60 ng	-	-	40 ng	-	-
4	4	6	-	-	-	-	40 ng	60 ng	-
5	5	6	-	-	60 ng	-	40 ng	-	-
6	6	6	-	-	-	60 ng	40 ng	-	-

Joshi

Cell Culture: Transfection


- min6 VLP transfection for qPCR
- 550 µL medium were exchanged before transfection








transfection scheme qPCR min6 05/10/2019 - in ng per well						
	condition	V30	V15	V42.1	V42.2	V14
1	Mock	500 ng	-	-	-	-
2	Only V15	300 ng	200 ng	-	-	-
3	Fusion 1 o.k.	-	200 ng	300 ng	-	-
4	Fusion 1 x	-	-	300 ng	-	200 ng
5	Fusion 2 o.k.	-	200 ng	-	300 ng	-
6	Fusion 2 x	-	-	-	300 ng	200 ng

- Transfection:

Transfectionmix 24-well plate 05/1...		
	A	B
1	DNA per well	500 ng
2	P3000 Reagent per well	1 µL
3	Lipofectamine 3000 reagent per well	0.75 µL
4	OptiMEM per well	2 x 25 µL

- Transfect cells according to the following table. Use the indicated volume of DNA and P3000™ Reagent with each of the two volumes of Lipofectamine™ 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.



Timing		Steps
1		Seed cells to be 70-90% confluent at transfection
2		<p>Dilute Lipofectamine<sup>™</sup> 3000 Reagent to target volume</p> <p>Prepare master mix of DNA by diluting DNA in Opti-MEM<sup>®</sup> Medium (2 tubes - Mix well)</p>
3		<p>Prepare master mix of DNA by diluting DNA in Opti-MEM<sup>®</sup> Medium (2 tubes - Mix well)</p> <p>Prepare master mix of P3000<sup>™</sup> Reagent, then add P3000<sup>™</sup> Reagent - Mix well</p>
4		Add Diluted DNA to each tube of Diluted Lipofectamine <sup>™</sup> 3000 Reagent (1:1 ratio)
5		Incubate
6		Add DNA-lipid complex to cells
7		Visualize/analyze transfected cells

Procedure Details (Two Reaction Optimization)				
Component	10-wells	24-wells	6-wells	
Adherent cells	1 × 10 <sup>5</sup>	0.5 × 10 <sup>5</sup>	0.25 × 1 × 10 <sup>5</sup>	
Opti-MEM <sup>®</sup> Medium	5 µL × 2	25 µL × 2	125 µL × 2	
Lipofectamine <sup>™</sup> 3000 Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL	
Opti-MEM <sup>®</sup> Medium	10 µL	50 µL	250 µL	
DNA (0.5-5 µg/µL)	0.2 µg	1 µg	5 µg	
P3000 <sup>™</sup> Reagent (2 µL/kg DNA)	0.4 µL	2 µL	10 µL	
Diluted DNA (with P3000 <sup>™</sup> Reagent)	5 µL	25 µL	125 µL	
Diluted Lipofectamine <sup>™</sup> 3000 Reagent	5 µL	25 µL	125 µL	
Incubate for 10-15 minutes at room temperature.				
Component	10-wells	24-wells	6-wells	
DNA-lipid complex	10 µL	50 µL	250 µL	
DNA amount	100 ng	500 ng	2500 ng	
P3000 <sup>™</sup> Reagent	0.2 µL	1 µL	5 µL	
Lipofectamine <sup>™</sup> 3000 Reagent used	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL	
Incubate cells for 2-4 days at 37°C. Then, analyze transfected cells.				