

Protocol: Purification via Heat-Soluble Precipitation

Purpose: To get purified IDPs from *E.coli* cell lysis from induction of protein expression via IPTG by heat shocking and verified the plasmid through SDS-PAGE

Material:

- Cell pellets and supernatants from induction of protein expression
- 25 mM Tris/HCl
- Water bath
- Centrifuge
- 2X SDS-PAGE sample buffer
- 1X SDS running buffer (see recipe)
- SDS stain solution (see recipe)
- SDS destain solution (see recipe)
- Mini protean TGX gels
- Power supply
- Electrophoresis Chambers

Methods:

1. Resuspend cell pellets in 500 μ L of 25 mM Tris/HCl and take 500 μ L aliquot from supernatants.
2. Heat shock in water bath at 80°C for 15 minutes.
3. Centrifuge at 16,000 rcf for 40 minutes at room temperature.
4. Separate supernatant and pellets in different eppendorf tubes. (Vortex them and make sure well dissolved)
5. Mix 6 μ L of supernatants and pellets with 6 μ L of 2X SDS-PAGE sample buffer (Vortex and spin down)
6. Heat shock in water bath at 80°C for 5 minutes.
7. Load 10 μ L of the sample in gel and run the gel at 200V for around half an hour.
8. Take out the gel as soon as finish running. Gently open it and wash it with buffer in the wash box.
9. Stain the gel with 25 ml of stain solution for 30 to 60 minutes with shaking.
10. Collect the stain solution for reuse. Destain gel with 25ml of destain solution with shaking overnight.
11. Next morning, destain twice for one hour each time.

Recipe:

- SDS running buffer: Tris-Glycine Electrophoresis Buffer (10X pH8.3)

Reagents	Volume	Final Concentration
Tris Base	30.3 g	250mM
Glycine	144 g	1.9M
SDS	10 g	1%
Water	To 1 L	

- SDS stain solution: “end of the day” stain (0.25% coomassie blue)

Reagents	Final Concentration
Coomassie Blue	0.25%
Acetic Acid	10%
Methanol	30%

- SDS destain solution

Reagents	Final Concentration
Acetic Acid	10%
Methanol	30%