

TSS PREPARATION

ADAPTED FROM THE BONNET TEAM PROTOCOL REPOSITORY

MATERIAL

- PEG 8000
- $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$
- DMSO
- LB medium
- Syringe
- 0.2 μm syringe filter
- 15mL Falcon tube

PROTOCOL

1. Mix 1.5g PEG 8000 and 25g of $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ in 625 μL of DMSO and 15mL of LB.
2. Filter sterilize it under the hood in the Falcon tube.