

# Vibrio transformation journal

---

## MONDAY, 7/1/2019

---

different methods of electroporation and chemical transformation have been tested.  
the working protocol is a chemical transformation protocol

## SUNDAY, 8/4/2019

---

**Plan:** check Transformation protocol in Vibrio, clear scheme on planned strains

3 Transformation samples: 1) no DNA, 2) plasmid with RFP coding device and AmpR, 3) plasmid with lacI-rbs-amiICP and KanR  
procedure of transformation in Vibrio: 3.5 ul of overnight culture in LBv2 diluted with 350 ul of fresh LBv2 and 30 ng of DNA  
incubated statically for 4 hrs at 30 °C; add 1 ml of LBv2 and incubate for 1 more hour at 30 °C shaking; plated on plates to grow  
overnight at 37 °C (no DNA on plate w/o antibiotics, Kan, Amp; plasmid 2 on plate w/o antibiotics & Amp; plasmid 17 on plate w/o  
antibiotics & Kan)

## MONDAY, 8/5/2019

---

transformation did not work, the transformation protocol used yesterday needs a gene (from *V. cholerae*) (for details checkout  
the link above)

found new protocol in <https://www-nature-com.proxy-ub.rug.nl/articles/nmeth.3970>; nuclease activity seems to be particularly  
high in Vibrio -> discuss ko of dns for chemical competency, try electroporation

## TUESDAY, 8/6/2019

---

made electrocompetent cells: inoculated 100 ml of BHlv2 with 1 ml of o/n culture, grown at 37 °C to reach OD 0.5, pelleted at 4°C,  
20 min, 6500 rpm, washed in electroporation buffer (680 mM sucrose in 7 mM KPi pH 7) 3 times (decant, resuspend, centrifuge 4  
°C, 20 min, 6750 rpm), resuspended to final OD 16 and aliquoted 400 µl into chilled tubes, frozen at -80 °C

## WEDNESDAY, 8/7/2019

---

Electroporation trial: Electrocompetent cells were thawed on ice, combined with 3 µl of plasmid with Kanamycin resistance (3  
from second isolated batch), transferred into 4 mm electroporation cuvette (chilled) and electroporated at 700 V, 25 µF, 200  
ohm, recovered with 500 µl recovery medium (BHlv2 supplemented with 680 mM sucrose, filtered), cultured for 1 hr at 37 °C, 50  
µl plated on small plate and grown o/n at 37 °C

unfortunately the cuvettes leaked in both cases, tried to get out what was there and washed them out with BHI once and  
combined for culturing in one tube. Medium not turbid after one hr of growth, still plated.

## THURSDAY, 8/8/2019

---

from yesterday's electroporation: growth on standard plate, no growth on plate with kanamycin, left in the incubator for another  
24 hr

Electroporation trial 2: same reagents as yesterday, 9 µl of plasmid with Amp resistance (plasmid two from second batch of  
isolation), 1 mm cuvette, 1: 600 V, 2: 800 V, cuvettes didn't leak this time (transferred immediately into culture tube)

## FRIDAY, 8/9/2019

---

one colony on plate with Kanamycin (I've got hopes, don't let me down Vibrio), growth on plates with and without antibiotic, also  
in the negative control. FIGURE OUT AMPICILLIN RESISTANCE

## TUESDAY, 8/20/2019

---

start of culture for chemically competent cells

**WEDNESDAY, 8/28/2019**

---

preparation of chemically competent cells and transformation

**THURSDAY, 8/29/2019**

---

colony PCR to confirm the plasmid in Vibrio after chemical transformation

wooop woop the transformation worked!

picked, and grew clones 2, 3 & 4

**TUESDAY, 9/3/2019**

---

1. isolate plasmid & try PCR again
2. try to get RFP by inducing with IPTG