

# 3A assembly side project

---

**Project:** iGEM uppsala 2019

**Authors:** JINWENYUAN

THURSDAY, 13/6/2019

---

MONDAY, 17/6/2019

---

## Prepare solutions for future usage (Johan)

### Aim:

prepare solutions as follows:

0.9% NaCl (10ml)

1 M CaCl<sub>2</sub> (10ml)

LB media (600ml)


LB Agar Plates (600ml) and Addition of Antibiotics

### Materials:

NaCl, CaCl<sub>2</sub>, Yeast extract, Bacto™ tryptone, agar, double distilled H<sub>2</sub>O and a 5M NaOH solution.

### Method:

A NaCl solution was prepared according to the protocol:

 0.9% NaCl (10ml)

The amount of NaCl that was used was 0.093g

A CaCl<sub>2</sub> solution was prepared according to the protocol but 1.11g of CaCl<sub>2</sub> was used

 1 M CaCl<sub>2</sub> (10ml)

Liquid-Broth media was produced according to the protocol:

 LB media (600 mL)

But the amount of chemicals used can be seen in the table below: NaCl used was 5.9619g and the amount of tryptone was 6g.

Amount of chemicals used.			
	A	B	C
1	NaCl	5.9619g	
2	BactoTMtryptone 1% (w/v)	6.0078g	
3	yeast extract 0.5% (w/v)	3.0009g	
4	ddH2O	600ml	
5	5 M NaOH	100ul	

A LB agar media was prepared by following the protocol:

☰ LB Agar Plates (600 ml LA) and Addition of Antibiotics

Table1			
	A	B	C
1	NaCl	5.9592g	
2	BactoTMtryptone 1% (w/v)	6.0027g	
3	yeast extract 0.5% (w/v)	3.0069g	
4	ddH2O	600ml	
5	5 M NaOH	100ul	
6	Agar	9.0016g	

## Transformation of pSB1C3\_K608006 (Johan & Manish)

### Aim:

Since other groups prepared: pSB1C3\_K592009 ( amciCP ) for the CDS part, our group only need to transform the plasmid: pSB1C3\_K608006 for the Promoter+RBS part, from the well 5E on plate 1 in 2019 IGEM distribution kit into DH5α competent cells for the amplification.

### Materials:

DNA Plasmid: pSB1C3\_K608006 (conc. 200~300 pg/uL, 1uL)

DH5α Competent cells (50 uL)

LB media

LA+Cm

### Procedure:

Following steps were done according to the protocol: Transformation

1. Dilution of DNA sample was done according to the instruction of 2019 IGEM distribution kit. (10  $\mu$ L ddH<sub>2</sub>O was added into the well 5E on Plate1 for resuspension of the plasmid:pSB1C3\_K608006 )
2. Competent cells was thawed on ice for more than 25 mins.
3. 1  $\mu$ L of plasmid resuspension solution and 1 $\mu$ L of ddH<sub>2</sub>O (as negative control) were added above to 50 $\mu$ L of competent cells.
4. Then the mixtures were Incubate for 30 minutes on ice
5. 950 $\mu$ L of LB Media (pre-heated to 37°C) were added to the mixtures.
6. Mixtures then were incubated for 80 minutes at 37°C, with occasional gentle mixing by inversion of the tubes.
7. For the mixture containing plasmid: pSB1C3\_K608006, it was mixed gently and plated 100 $\mu$ L (=1/10th) on an agar plate containing kanamycin.
8. For the rest mixures, they were spined down at 4000rpm for 5min, at RM (Room Temperature).
9. All but 100 $\mu$ L of the supernatant of the remaining mixtures were discarded and the pellets were resuspended in the remaining 100 $\mu$ L.
10. The remaining suspension were spread on an agar plate.

**To-to-list for following days :**

- Day 1:** Tranformation
- Day 2:** Calculate transformation efficiency (colonies/ $\mu$ g) of the competent cells. In the evening, re-streak appropriate colonies
- Day 3:** Set overnight cultures for preparation of plasmid miniprep.
- Day 4:** Make plasmid miniprep and glycerol cell stocks.

TUESDAY, 18/6/2019

---

**Continued from 2019/06/17:****Conclusion and discussion:**

We failed the transformation yesterday, no colonies on plates of all groups. It may be due to the bad transformation efficiency of our competent cells since the competent cells we used yesterday were from last year. So, we plan to make our competent cells again.

**Preperation of competent E.coli(DH5a) cells (Johan, Manish & Lilli)****Aim:**

Due to the failed transformation yesterday new competent cells are neccery to perform a new transformation of the plasmid. The new copotent cells are hoped to bring a succesfull transformation.

**Protocol:**

The copotent cells are produced by the following protocol:

 Preperation of E.coli cells for transformation

WEDNESDAY, 19/6/2019

---

**Prepare solutions for future usage (Johan & Manish)****Aim:**

prepare solutions as follows:  
2xLB media (600mL)

LB media with low NaCl concentration (600mL)

## Materials:

NaCl , Yeast extract , Bacto™tryptone, agar, double distilled H<sub>2</sub>O and a 5M NaOH solution.

## Method

Two Liquid-Broth media was produced according to the protocol:

LB media (600 mL)

But the NaOH was added after autoclaving and was then autoclaved again.

A Liquid-Broth media with low NaCl concentration was produced according to the protocol:

LB media (600 mL)

But the amount of NaCl that was added was 2.98 g.

THURSDAY, 20/6/2019

## 3A Assembly and gel analysis (Lilli)


Aim: to build a J23110+RBS+amciCP(CDS)+vector(pSB1K3) construct

## Materials:

ddH<sub>2</sub>O, 2 DNA plasmids: psBIC3 and psBIK3/psBIAS, 10x reaction buffer, 3 restriction endonucleases: EcoRI, XbaI, and SpeI.

## Method:

Proceeded according to protocol:





MONDAY, 1/7/2019

---

## 3A assembly on Jinwens samples (Lilli)

### Digestion

#### Aim:

To build a Lac+Bba\_K50003 construct (we are not gonna do the backbone now, another group)

#### Materials:

94.5 ng/μl Gene (Bba\_K50003), 83,1 ng/μl Promotor (Lac), ddH<sub>2</sub>O, 5 μl Reaction Buffer, 1 μl EcoR1, Spe1, and Xba1.

#### Methode:

Proceeded according to protocol 3.

We also did control samples of the gene and promotor with a total of 20 μl.

### Gel analysis

#### Aim:

See if the digestion was succesful.

#### Materials:

Our digestion samples (Promotor, Control promoter, Gene, Control gene) , Purple Dye, Loading dye, GeneRuler. (The gel was already pre-done)

#### Methode:

Methode proceeded according to protocol 3.

#### Result:



Ladder, Promotor, Control Promotor, Gene, Control Gene

TUESDAY, 2/7/2019

---

## Digestion of backbone pSB1C3 (Lilli)

### Aim:

Part of yesterdays 3A assembly since we didnt have a backbone. We want to cut the backbone so its complemetary to our gene and promotor samples.

### Materials:

25 ng/ $\mu$ l Backbone (pSB1C3), ddH<sub>2</sub>O, Solution Buffer, EcoR1 and Pst1

### Method:

Method proceeded according to protocol 3.

Made a 500 ng sample of backbone with enzymes and a control sample of 200 ng backbone.

## Gel analysis of digested backbone

Proceeded like yesterdays gel analysis.

### Results:

Not sucesfull. No result.

## Ligation

### Aim:

Since our backbone digestion not was sucesful we are gonna use Vanjas groups backbone (PSB1A). This we are gonna ligate with our Dlgested gene and promotor from yesterday.

Bba\_K50003 + Lac + PSB1A3

### Materials:

2 µl of Bba\_K50003, Lac and PSB1A3, 11 µl ddH<sub>2</sub>O, 2 µl reaction buffer, 1 µl T4 DNA ligase

### Method:

Method proceeded according to protocol 3.

During the incubation in room temperature i accidentally put the sample in ice for about 5 min. So it was incubated in room temperature for about 15 min then put on ice, then incubated at 80 degrees.

## Transformation

### Aim:

Putting our, hopefully ligated, plasmids in to bacteria.

### Materials:

Everything in protocol 6.

### Method:

Proceeded according to protcol 6.

---

THURSDAY, 4/7/2019

## Transformation (Lilli)

### Aim:

We are changing our pSBIK3 backbone because we did not digest it properly according to our gel analysis. We are instead using the ampicilin resistance pSBIA3 bacbone. Transforming a plasmid with the compondmnts Lac+ GLOX+ pSBIA3. This so we can put our plasmid in to backterias.

### Materials/Method:

See protocol 6.

---

MONDAY, 8/7/2019

## 3A assembly(Johan and Manish)

**Aim:**

The aim of this 3A assembly is to produce a Lac+GLOX+pSBIA3 plasmid.

## Digestion (Johan and Manish)

**Aim:**

To prepare the different genes for ligation.

**Materials:**

See the following protocols:

 3A assembly-Digestion and ligation

The concentration of the samples used and the volume used to reach a amount of 500 ng.

	<b>Components</b>	<b>Concentration s</b>	<b>Volume</b>
1	GLOX	83.1 ng/mL	6 $\mu$ L
2	Lac	90.4 ng/mL	5.5 $\mu$ L
3	pSBIA3	25 ng/mL	20 $\mu$ L

**Method:**

See the following protocols:

 3A assembly-Digestion and ligation

## Ligation (Johan and Manish)

**Aim:**

To prepare the construct for transformation.

**Materials:**

See the following protocols:

 3A assembly-Digestion and ligation

**Method:**

See the following protocols:

 3A assembly-Digestion and ligation

## Transformation(Johan and Manish)

### Aim:

To introduce the construct into bacteria.

### Materials:

See the following protocols:

 3A assembly-Digestion and ligation

### Method:

See the following protocols:

 3A assembly-Digestion and ligation

With the exception that the heat shock was carried out in 80°C instead of 42°C

WEDNESDAY, 10/7/2019

---

No new growth was detected on agar plates with the transformed cells.

THURSDAY, 11/7/2019

---

## 3A assembly(Johan and Manish)

## Ligation (Johan and Manish)

### Aim:

To prepare the construct for transformation.

### Materials:

See the following protocols:

 3A assembly-Digestion and ligation

### Method:

See the following protocols:

The amount of GLOX used was 9  $\mu$ L (90 ng) of GLOX instead of 2  $\mu$ L (20 ng).

3A assembly-Digestion and ligation

## Transformation(Johan and Manish)

### Aim:

To introduce the construct into bacteria.

### Materials:

See the following protocols:

3A assembly-Digestion and ligation

### Method:


See the following protocols:

3A assembly-Digestion and ligation

The temperature for heat shock was 80°C instead of 42°C

## 0.7% Agarose gel analysis (Johan and Manish)

### Results:



**Conclusion:**

There is suspected that there is a problem with the enzymes EcoRI

MONDAY, 15/7/2019

---

**3A assembly(Johan and Manish)****Aim:**

The aim of this 3A assembly is to produce a Lac+GLOX+pSBIA3 plasmid.

**Digestion (Johan and Manish)****Aim:**

To prepare the different genes for ligation.

**Materials:**

- (1) 8-tube strip, or (3) 0.6ml thin-walled tubes
- BioBrick Part in BioBrick plasmid (Purified DNA, > 16ng/ul)
- dH2O
- NEB Buffer 2
- BSA
- Restriction Enzymes: EcoRI, SpeI, XbaI, PstI

	<b>Components</b>	<b>Concentration s</b>	<b>Volume</b>
1	GLOX	83.1 ng/mL	3 $\mu$ L
2	Lac	90.4 ng/mL	3 $\mu$ L
3	pSBIA3	25 ng/mL	3 $\mu$ L

**Method:**

1. Add 250ng of DNA to be digested, and adjust with dH2O for a total volume of 16ul.
2. Add 2.5ul of NEBuffer 2.
3. Add 0.5ul of BSA.
4. Add 0.5ul of EcoRI.
5. Add 0.5ul of PstI.
6. There should be a total volume of 20ul. Mix well and spin down briefly.
7. Incubate the restriction digest at 37C for 30min, and then 80C for 20min to heat kill the enzymes. *We incubate in a thermal cycler with a heated lid*
8. Run a portion of the digest on a gel (8ul, 100ng), to check that both plasmid backbone and part length are accurate.

**Validity analysis(Johan and Manish)****Aim:**

The aim of this was to analyse the validity of the restriction enzymes EcoRI, XbaI, SmaI and PstI

### Material:

Agarose  
1xTBE buffer

### Method:

A 0.7% w/w agarose gel was made and the digest was put in the wells and electrophoresis was run at ~90 V for two hours.

## Ligation (Johan and Manish)

### Aim:

To prepare the construct for transformation.

### Materials and Method:

- Add 2ul of digested plasmid backbone (25 ng)
- Add equimolar amount of EcoRI-HF SmaI digested fragment (< 3 ul)
- Add equimolar amount of XbaI PstI digested fragment (< 3 ul)
- Add 1 ul T4 DNA ligase buffer. **Note:** Do not use quick ligase
- Add 0.5 ul T4 DNA ligase
- Add water to 10 ul
- Ligate 16C/30 min, heat kill 80C/20 min

## Transformation (Johan and Manish)

### Aim:

To introduce the construct into bacteria.

### Materials:

See the following protocols:

 3A assembly-Digestion and ligation

### Method:

See the following protocols:

 3A assembly-Digestion and ligation

TUESDAY, 16/7/2019

---

Two cultures had appeared on the plates of the transformed cells one red and one transparent

Overnight cultures were run by putting cell colonies in LB and leaving overnight for 14 hours in 37°C




WEDNESDAY, 17/7/2019

Miniprep was run on the overnight cultures in accordance with the instructions for the miniprep kit. Analytical digestion and gel electrophoresis was done in accordance with the protocol:

#### Analytical Digestion and Agarose Gel Electrophoresis

iGEM\_trans\_before\_gibbs.jpg



THURSDAY, 18/7/2019

Cells from prior this week was digested with the digestion part of the protocol:

#### 3A assembly-Digestion and ligation

Colony PCR was performed on the digested cell in accordance with the protocol:

#### PCR for synthesizing DNA fragments

Gel electrophoresis was performed on the digest and bacteria in accordance with:

 3A assembly-Digestion and ligation

TUESDAY, 23/7/2019

---

## 3A Assembly (Lilli)

# Main project

---

**Project:** iGEM uppsala 2019

**Authors:** JINWENYUAN

THURSDAY, 20/6/2019

---

FRIDAY, 21/6/2019

---

## Transformation of Shuttle vector (Lilli)

To amplify our shuttle vector pPICZaB through transformation.

### Materials:

150 µl Competent cells, 1 ng/µl pPICz α-B dilution and ddH<sub>2</sub>O.

### Methode:

The methode was proceeded according to this protocol:

 Transformation of Competent E.coli cells

However in step 2 our competent cells were thawed for 30 min on ice instead o 15 min.

In step 3 we added 10 µl of our shuttle to 50 µl of our competent cells and in another tube 10 µl of ddH<sub>2</sub>O to 50 µl competent cells.

TUESDAY, 25/6/2019

---

## O/N culture (Lilli)

### Aim:

Grow an O/N culture so we can peform a plasmid miniprep.

### Materials:

100 ml LB media and 25 µl Zeocin

### Methode:

The methode was proceeded according to protocol 2.

6 ml LB + Zeocin mixture was placed in each of the 6 tubes.

(This O/N culture was discarded the following day)

WEDNESDAY, 26/6/2019

---

## O/N culture (Lilli)

### Aim:

Grow an O/N culture so we can peform miniprep

### Materials:

Low salt LB meda and Zeocin

## Method:

The method was proceeded according to protocol 2.

50 ml LB media and 2,5 µl Zeocin.

6 ml of mixture was placed in 6 falcon tubes each.

---

THURSDAY, 27/6/2019

## Miniprepkit (Lilli)

### Aim:

To purify our pPICzaB shutter vector. So we can do an analytic digestion with our samples.

### Materials:

Our O/N culture, 250 µl Resuspension Solution, 250 µl Lysis solution, 350 µl Neutralization Solution, 500 µl Wash Solution and 50 µl Elution Buffer

### Method:

Method was proceeded according to Thermo Scientific GeneJET Plasmid Miniprep Kit (Quick Protocol)

Values for purified DNA was also checked. (Can be found in Lillis lab journal if needed)

## Digestion

### Aim:

To cut our shuttle vector.

### Materials:

200 ng DNA, 20 µl digeston mixture, EcoR1

### Method:

Proceeded according to protocol 3.

6 samples.

---

MONDAY, 1/7/2019

## Gel analysis on Lillis pPICZaB samples (Lilli)

### Aim:

See if the digsteion was succesful.

### Materials:

The 6 digested samples , ladder (1 µl GeneRuler, 1 µl Purple Dye, 4 µl ddH<sub>2</sub>O)

### Method:

Methoded proceeded according to protocol 4

Samples: Total volume of 6  $\mu$ l ( 5  $\mu$ l digested samples and 1  $\mu$ l Purple Dye)

### Result:

The gel was empty, no results. Probably something was wrong with the gel or the Purple Dye.

TUESDAY, 2/7/2019


---

## Gel analysis on Lillis pPICZaB samples #2 (Lilli)

Everything was proceeded like yesterday except that the Purple Dye was exanged for BLue Dye.

### Result:

The results this time were pretty OK. We can see the bands and hey are about the right size they should be. The shutter vector samples are saved.



Ladder, sample 1-6

TUESDAY, 23/7/2019

---

PCR to add overhangs, which have overlaps to pPICZaB, to AAO\_P.pastoris, HRP\_P.pastoris and MnP, LiP, GLOX genes. (Jinwen & Johan)

refer to PCR protocol (WHICH WAS WRONG ORIGINALLY)

Result: we got nothing, according to gels.

We did this again today, refer to PCR protocol.

## HRP reaction mix

	A	B
1	ddH <sub>2</sub> O	14.5µL
2	2 mM dNTPs	5 µL
3	Forward primer (5 µM)	5 µL
4	Reverse primer (5 µM)	5 µL
5	5x Phusion® HF buffer	10 µL
6	HRP ( 10ng/µL)	10µL
7	Phusion® HF DNA polymerase	0.5 µL
8	Total	50 µL

## AAO reaction mix

	A	B
1	ddH <sub>2</sub> O	14.5µL
2	2 mM dNTPs	5 µL
3	Forward primer (5 µM)	5 µL
4	Reverse primer (5 µM)	5 µL
5	5x Phusion® HF buffer	10 µL
6	AAO ( 10ng/µL)	10µL
7	Phusion® HF DNA polymerase	0.5 µL
8	Total	50 µL

## PCR program setting up for HRP

**Table3** ^

	A	B	C	D
1	Description of Step	Temp (°C)	Time: (hh:mm:ss)	Number of cycles
2	Initial denaturation	98	<u>00:00:30</u>	1x
3	Denaturation	98	<u>00:00:10</u>	2x
4	Annealing temp. + 4°C	67	00:00:30	
5	Extension	72	00:00:38	
6	Denaturation	98	<u>00:00:10</u>	6x
7	Annealing temp. + 2°C	63	00:00:30	
8	Extension	72	00:00:38	
9	Denaturation	98	<u>00:00:10</u>	25x
10	Extension	72	00:00:38	
11	Final extension	72	00:07:00	1x
12	Storage	4	∞	1x

PCR program setting up for HRP

**Table4** ^

	A	B	C	D
1	Description of Step	Temp (°C)	Time: (hh:mm:ss)	Number of cycles
2	Initial denaturation	98	<u>00:00:30</u>	1x
3	Denaturation	98	<u>00:00:10</u>	2x
4	Annealing temp. + 4°C	63	00:00:30	
5	Extension	72	00:01:00	6x
6	Denaturation	98	<u>00:00:10</u>	
7	Annealing temp. + 2°C	61	00:00:30	
8	Extension	72	00:01:00	25x
9	Denaturation	98	<u>00:00:10</u>	
10	Extension	72	00:01:00	
11	Final extension	72	00:07:00	1x
12	Storage	4	∞	1x

Gel result:

WEDNESDAY, 24/7/2019

---

Gibson assembly of AAO/HRP/MnP/GLOX/LiP to pPICZaB

prepare reaction mix for AAO



Table5			
	A	B	C
1		test group	Positive Control**
2	AAO_PCRpurification	2.4 µl	10 µl
3	pPICZaB_PCR purification	7.2 µl	
4	Gibson Assembly Master Mix (2X)	10 µl	10 µl
5	Deionized H2O	0.4 µl	0
6	Total Volume	20 µl	20 µl



TUESDAY, 30/7/2019

## PCR for pPICZaB\_HRP minipreps(1,2,3) (Jinwen)

see PCR protocol

## miniprep for pPICZaB\_AAO\_2A\_eGFP(E1,E2,E3,E4,E5) & PCR & Gel analysis

see miniprep, PCR, Gel analysis protocol

## Colony PCR for pPICZaB\_HRP colonies (Jinwen)

see colony PCR protocol

## PCR for pPICZaB\_HRP minipreps(1,2,3)(Johan)

The pPICZaB\_HRP minipreps was removed from the freezer and put on ice. Then the following PCR protocol was followed:

 PCR (Dreamtaq)

The annealing temperature was set to 50°C and the extension time was set to 1 minute and 30 seconds. A master mix of dNTPs, 2 mM; phosphorylated forward primer, 5 µM, phosphorylated reverse primer, 5 µM and 10X DreamTaq Buffer. A 1% agarose gel was made and the three samples were analysed next to a 1kb gene ruler ladder.

WEDNESDAY, 31/7/2019

MONDAY, 5/8/2019

## Sequencing of PCR products. (Johan)

### Aim:

To investigate if transformation was succesfull.

### Materials:

ddH<sub>2</sub>O, Samples that code for diffrent genes and SeqPrimer.

### Methode:

Teplates were made of 15 µl of purified DNA with a concentration of about 10 ng/µl of the DNA and then 2 µl of Seqprimer of a concentration of 10 pmol/µl the tubes were then seal and sent for sequencing. The serialnumber and the content of the test tubes can be found in the table below:

	A	B
1	EF30500622	HRP-2A-AAO-2
2	EF30500623	HRP-2A-AAO-1
3	EF30500624	PPIC-GLOX-3
4	EF30500625	PPIC-MnP-2
5	EF30500626	PPIC-MnP-3
6	EF30500612	PPIC-MnP-1
7	EF30500611	PPIC-HRP-1
8	EF30500610	PPIC-HRP-2
9	EF30500609	PPIC-HRP-3

WEDNESDAY, 7/8/2019

## Gibson assembly (Jinwen)

### Aim:

insert HRP\_2A\_eGFP into pPICZaB & insert AAO\_2A\_eGFP into pPICZaB

### Materials:

HRP\_2A\_eGFP DNA fragments

AAO\_2A\_eGFP DNA fragments

ddH<sub>2</sub>O

Gibson assembly Master Mix

### Methods:

refer to the Gibson assembly protocol

## Transformation (Jinwen)

**Aim:**

to transform the assembled products from 8/7 and also from 26/7 into DH5a for further construct amplification

**Materials:**

Gibson assembled products (HRP\_2A\_eGFP+pPIC, AAO\_2A\_eGFP+pPIC, AAO+pPIC)

DH5a super competent cells from NEB

ddH<sub>2</sub>O (negative control)

pPICZaB (positive control)

SOB LB media

LSLB(+zeocin) plates

**Methods:**

refer to the transformation protocol

## Ordered DNA fragments and primers for signal peptide mutation side project (Jinwen)

Sequence confirmed for pPICZaB-MnP, HRP, HRP\_2A\_AAO construct, and for GLOX, we had a point mutation, which may come from the distribution kit, or caused by amplification in DH5A.

THURSDAY, 8/8/2019

---

We got plenty colonies from transformation yesterday

### O/N culture from 6 colonies on each plate (Jinwen, Qian, Erik, Irina)

see O/N culture protocol

#### Plan for the next day:

1. colony PCR and Gel analysis in the morning, and glycerol stock from O/N of each sample
2. if we got any right size gene in the gel, prepare midiprep O/N for the next day after tomorrow from the relevant glycerol stock

FRIDAY, 9/8/2019

---

## Colony PCR, glycerol stock & Gel analysis (Jinwen, Qian, Johan, Erik, Jenny)

**Aim:**

verify the right construct among colonies on the transformation plates, through making colony PCR directly on cells in O/N yesterday with pPIC sequencing primers.

**Methods:**

refer to glycerol stock protocol

600 micro liters culture from each sample, mixed with 400 micro liters 50% glycerol stock


refer to colony PCR protocol

400 microliters culture from each sample were spun down (6000 rpm, 5 mins), cells were taken by dipping into the pellet with a tip and resuspended into 30 microliters water.

refer to gel analysis protocol


Results:

pPICZaB\_AAO (Jinwen)




conclusion & discussion: all samples from AAO failed, may due to the wrong Gibsonassembly products was taken, redo AAO's transformation with Qian's groups' AAO+pPIC Gibson assembly product. Samples may be mixed up with AAO-2A-eGFP, this may explain why we got band around 3kb.

pPICZaB\_AAO\_2A\_eGFP (Jenny)




#2& #3 Glycerol stocks were discarded, #6 was taken to make O/N for the nextday's midiprep.

pPICZaB\_AAO\_2A\_eGFP (Johan)




From the gel, the upward band in Lane 1&3 is correct, but Lane 1&3 also contains the band around 300bp, so we supposed that Sample 3 was contaminated with cells carrying the empty plasmid. All samples related with this gel were discarded.

pPICZaB\_HRP\_2A\_eGFP (Erik)



the amplicon we expected here for HRP-2A-eGFP is around 2100bp, so we got the right construct in the #2, #3, #5, #6 samples of Erik's. 1&4 's glycerol stock was discarded.



#2 was discarded, the others are all fine

## O/N culture (50ml) for the next day's midiprep (Jinwen & Qian)

### Aim:

to prepare DNA for transformation in yeast, which needs at least 1 micro grams DNA.

### methods:

refer to O/N protocol

## Plan for the next:

1. midi prep for the O/N culture of AAO-2A-eGFP & HRP-2A-eGFP (done by 10/8)
2. DNA preparation for transformation in yeast (DIgestion with SacI) (done by 10/8)
3. Transformation of the Gibson assembly product for AAO & LiP (done by 11/8 )

MONDAY, 12/8/2019

---

## Sequencing

Send construct of pPICZaB-HRP\_2A\_eGFP & pPICZaB\_AAO\_2A\_eGFP & AAO? for sequencing with pPIC sequencing fwd primer.

	A	B
1	EF30500574	pPICZaB-AAO-midiprep
2	EF30500575	pPICZaB-HRP-2A-eGFP midiprep
3	EF30500576	pPICZaB-AAO-2A-eGFP midiprep

We got colonies on transformation plates for pPICZaB-AAO and pPICZaB-Lip GA products, and grew O/N culture on Sunday.

## Colony PCR with pPICZ sequencing primer and gel analysis for pPICZaB-AAO and pPICZaB-Lip constructs.

### Aim:

to check if we get the right construct or not.

### method:

refer to colony PCR protocol

### result:

uv room computer , by qian

WEDNESDAY, 14/8/2019

	A	B
1	EF30500577	pPICZaB-AAO-midiprep
2	EF30500578	pPICZaB-LiP midiprep
3		

FRIDAY, 16/8/2019

PCR to add overhangs to LiP/MnP in pPICZaB construct, and to mAAO fragment ordered from IDT.

Components in PCR for adding overhangs				
	A	B	C	D
1	LiP PCR	1 ng PICZaB-LiP	25uM OvH_ProK_LiP_fwd	25uM pPIC-XXX_Rev
2	MnP PCR	1 ng PICZaB-MnP	25uM OvH_ProK_MnP_fw d	25uM pPIC-XXX_Rev
3	mAAO PCR to add overlap with proa	2 ng mAAO	25uM OvH_Proa_mAAO_f wd	25uM mAAO_Rev
4	mAAO PCR to add overlap with proK	2 ng mAAO	25uM OvH_ProK_mAAO_f wd	25uM mAAO_Rev
5	pPICZaB linerize w/o a signal peptide	1 ng PICZaB	25uM mutLinerize_fwd	25uM mutLinerize_Rev

refer to protocol

Gel analysis of adding overhangs result