

Autor: Maurice Mager

erstellt: 10.07.2019 17:04

Eintrag 1/9: promoter placeholder construct

aktualisiert: 17.09.2019 11:57

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

Autor: Maurice Mager

erstellt: 11.07.2019 16:27

Eintrag 2/9: RBS placeholder construct

aktualisiert: 17.09.2019 11:57

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

Autor: Maurice Mager

erstellt: 11.07.2019 16:36

Eintrag 3/9: Ter placeholder construct

aktualisiert: 17.09.2019 11:57

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

Autor: Marian Krämer

erstellt: 23.08.2019 16:19

Eintrag 4/9: Noch kein Eintragstitel

aktualisiert: 23.08.2019 16:19

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

Autor: Maurice Mager

erstellt: 13.09.2019 22:09

Eintrag 5/9: PRO, RBS and TER spaceholder without oriT and eYFP

aktualisiert: 20.10.2019 20:37

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

25.08

so far: transformations on spec plates resulted in huge background growth

we dont know wether its the cells or the plates but gradient agar plates showed that at 4x spec concentration in agar plates, the should be no background but colonies with the spec utex backbone part should still grow.

28.08

made agar plates with 4x spectinomycin concentration

03.09

remade the spaceholder GG with the table below and transformed into Top10

05.09

colonies were big enough to pick from

06.09

only promoter ONC were green, so those were prepped with the Macherey Nagel kit

New colonies were picked from plates

07.09

only one ONC was green again (RBS.SH.1) but i prepped all. The yield however was too low for a digest, hence they were retransformed for new ONCs.

test digesting RBS.SH.1, PRO.SH.1, PRO.SH.2 with bsmbI (see protocol)

results see gel picture no.1

lane 1: Pro SH.1

lane 2:Pro SH.2

lane 3: RBS SH.1

lane 4: generuler 1 kb ladder

Promoter spaceholder constructs were clearly correctly assembled (6 kb band for the ori part and 2 kb band for the cut out parts from 5' to 3' connector) and could be send for sequencing.

retransformation of RBS.SH.2 and TER.SH1/2 yielded green colonies, so they seem to only express the gfp on plates.

10.09

prepared ONCs of RBS.SH1/2 and TER.SH.1/2

sequencing of PRO.SH.1 didnt work but it was because johanna used the wrong primer

11.09

prepped and digested for 2h at 37°C (see table)

repeated seq of Pro.SH.1

12.09

sequencing of Pro.SH.1 looks good, rev sequencing needs to be sent away

GG spacer parts

	A	B	C	D	E	F	G	H	I
1	pro								
2	1_05	1 µl							
3	Pro SH	1 µl							
4	3_07	1 µl							
5	4_04	1 µl							
6	5_03	1 µl							
7	6_02	1 µl							
8	utex_ori 10 ng/ µl	1 µl							
9	bsal	0,5 µl							
10	t4 Puffer	1 µl							
11	ligase	1 µl							
12									
13	rbs								
14	1_05	1 µl							
15	2_20	1 µl							
16	RBS SH	1 µl							
17	4_04	1 µl							
18	5_03	1 µl							
19	6_02	1 µl							
20	utex_ori 10 ng/ µl	1 µl							
21	bsal	0,5 µl							
22	t4 Puffer	1 µl							

206										
-----	--	--	--	--	--	--	--	--	--	--

Sheet1

 [gg_spaceholder_parts_4612415_2.xlsx](#)


test digest bsmbl

	A	B	C	D	E	F	G	H	I	J
1	component	volume								
2	CutSmart	5 µl								
3	Bsmbl	1 µl								
4	Plasmid	1 ng								
5	A.dest	to 50 µl								
6										
7										
8	digest for 2h at 37°C									
9										
10										
11										
12										
13										
14										
15										
16										
17										
18										
19										
20										

Sheet1

 [test_digest_bsmbl_4612415_3.xlsx](#)

gel_no.1.jpg



digest table 11.09

	A	B	C	D	E	F	G	H	I
1	ingredients	RBSSH.1	RBSSH.2	TERS.H.1	TERS.H.2				
2	DNA	10	2	5	4				
3	cutsmart buffer	5	5	5	5				
4	Bsmbl	1	1	1	1				
5	A.dest	34	42	39	40				
6									
7									
8	digest at 37°C for 2h								
9									
10									
11									
12									
13									
14									
15									
16									
17									
18									
19									
20									

Sheet1

 [digest_table_11_4612415_5.xlsx](#)

Autor: Maurice Mager

erstellt: 17.09.2019 10:03

Eintrag 6/9: spaceholder with oriT

aktualisiert: 01.10.2019 13:03

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

17.09

since the oriT part is correct, placeholder constructs can be assembled with this part
started golden gate (program for SH constructs, leave out the final digest) with the composition in Table 1

18.09

transformed 5 µl of GG into Top10

19.09

prepared 12 ONC of Pro.SH and 3 ONC of RBS.SH. PRO.SH colonies did not fluorescent green, just as in the construct without oriT.
RBS.SH GG only yielded in 3 presumably positive colonies

20.09

prepped 6 Pro.SH and 3 Pro.SH ONCs

prepared bsmbl digests according to standard protocol

no correct bands

21.09

no correct bands for the next 6 preps either. golden gate needs to be repeated

1.10

repeated GG for both parts

transformed into Top10

sent 3 samples of PROSHoriT away for sequencing

sent 2 samples of RBSSHoriT away for sequencing

Labfolder Table

	A	B	C	D	E	F	G	H
1	Pro placeholder			RBS placeholder				
2	part	amount	actual	part	amount	actual		
3	1_05	24 ng	30	1_05	24 ng	30		
4	SH Pro Rene	37 ng	37	2_20	24 ng	28		
5	3_07	24 ng	21	RBS SH Rene	37 ng	35		
6	sfYFP	37 ng	27	sfYFP	37 ng	27		
7	5_03	24 ng	29	5_03	24 ng	29		
8	bom site	30 ng	27	bom site	30 ng	27		
9	spec UTEX ori	10 ng	12	spec UTEX ori	10 ng	12		
10	total: 7µl							
11	add 1µl of T4 Lig, 1 µl T4-Buffer, 1µl Bsal to the mix							
12								
13								
14								
15								
16								
17								
18								
19								
20								

Sheet1

Autor: Maurice Mager

erstellt: 20.09.2019 14:52

Eintrag 7/9: using PRO.SH-oriT to insert promoters for PCC7942

aktualisiert: 01.10.2019 17:14

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

23.09

PRO.SH1 without oriT is used to build promoter measurement constructs for PCC7942

PRO.SH1 is digested with BsmBI according to protocol

PRO.SH.1 is purified over an agarose gel

started 20x Lvl1,5 GG reactions with purified Pro.SH1 according to table below

used the normal Golden Gate Standard cycle protocol

24.09

Transformation into Dh5a

wrong enzyme ! needs to be repeated with Bsal

1.10.09

SPoT is done, this project is postponed

Lvi 1,5 GG

	A	B	C	D	E	F	G	H
1	23x MM							
2	ingredient	volume						
3	digested spacerholder 20 ng/ μ	23 μ l						
4	T4 Buffer	23 μ l						
5	Bsmbl	20 μ l						
6	T4Ligase	20 μ l						
7	A. dest	120 μ l						
8								
9	add 1 μ l of 60 ng Promoter part to 9 μ l of Mastermix							
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								

Sheet1

 [lvi_1,5_gg_4527531_2.xlsx](#)

Lvl 0,5 GG

	A	B	C	D	E	F	G
1	23x MM						
2	ingredient	volume					
3	digested spaceholder 20 ng/ μ l	23 μ l					
4	T4 Buffer	23 μ l					
5	Bsal	20 μ l					
6	T4Ligase	20 μ l					
7	A. dest	120 μ l					
8							
9	add 1 μ l of 60 ng Promoter part to 9 μ l of Mastermix						
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							

Sheet1

 [lvl_0,5_gg_4527531_3.xlsx](#)

Autor: Maurice Mager

erstellt: 01.10.2019 17:22

Eintrag 8/9: inserting promoters into PRO spacer + oriT

aktualisiert: 01.10.2019 22:48

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

PRO.SH+oriT will be used to insert all promoters from the Marburg Collection for measurements in UTEX2973 and PCC7942

1.10

Bsal digest according to standard protocol

Gel extraction (21,0 ng/μl)

LVL 0,5 GG according to table "LVL0,5 GG"

LVL 0,5 GG

	A	B	C	D	E	F	G
1	23x MM						
2	ingredient	volume					
3	digested SH	23 μl					
4	T4Buffer	23 μl					
5	Bsal	23 μl					
6	T4Ligase	23 μl					
7	A.dest	133 μl					
8							
9							
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							

Sheet1

Autor: Joana Esslen

erstellt: 14.10.2019 19:25

Eintrag 9/9: Spaceholder (Pro/RBS) Lvl1 GG

aktualisiert: 20.10.2019 21:12

In Projekt: Toolbox: building placeholder constructs lvl 1

Keine Tags verwendet

12.10.2019lvl 1 Golden Gate:

1. Parts on ice, buffer on ice
2. Prepared Mastermixes in PCR-tube (out of parts)
3. Added 1 μL lvl 1 UTEX bb ($c = 12.0 \text{ ng}/\mu\text{L}$)
4. Added 1 μL T4 Lig buffer
6. Enzymes in ice
7. Added 1 μL of T4 ligase
8. Added 1 μL of Bsa1
9. Pipet-mixing
10. Golden Gate program in Mastercycler (Eppendorf) (GG Spaceholder)

Mastermix 1 Pro

name of part	c (part) / ($\text{ng}/\mu\text{L}$)	V (part) / μL
1_05	26	1.00
Pro SH	40	1.00
3_07	27	1.00
4_26	28	1.00
5_03	27	1.00
6_02	22	1.00

Mastermix 2 Pro

name of part	c (part) /(ng/ μ L)	V (part) / μ L
1_05	26	1.00
Pro SH	40	1.00
3_07	27	1.00
sfYFP	31	1.00
5_03	27	1.00
6_02	22	1.00

Mastermix 3 RBS

name of part	c (part) /(ng/ μ L)	V (part) / μ L
1_05	26	1.00
2_20	23	1.00
RBS SH	40	1.00
4_26	28	1.00
5_03	27	1.00
6_02	22	1.00

Mastermix 4 RBS

name of part	c (part) /(ng/ μ L)	V (part) / μ L
1_05	26	1.00
2_20	23	1.00
RBS SH	40	1.00
sfYFP	31	1.00
5_03	27	1.00
6_02	22	1.00

13.10.2019

Transformation:

1. 5 μ L of the GG lvi1 in top10
2. incubation on ice for 30 minutes
3. heatshock (42°C, 90 sec)
4. incubation in ice for 5 minutes
5. 500 μ L of LB medium added
6. incubation (37°C, 1h)
7. centrifugation
8. remove supernatant
9. plate on agar plate with Spec
10. in incubator 37°C

14.10.2019

1. picked 3 colonies from each plate for overnight cultures (all picked colonies were green)
2. in incubator 37°C, 210rpm

15.10.2019

nothing grew in 1_1

Plasmid DNA purification (Macherey-Nagel kit):

see protocoll; warmed Elution buffer to 80°C and just added 30 μ L Elution buffer

retrafo of 1_3

new overnight cultures of 1Pro (1_4, 1_5)

Bsmbl digest of 1 2; 2 1; 2 2; 2 3; 3 1; 3 2; 3 3; 4 1; 4 2; 4 3

1. mixed components:

1_2: 4.0 μL of 1_2: (253 $\text{ng}/\mu\text{L}$) + 40.0 μL of water added in one tube

2_1: 3.0 μL of 2_1 (354 $\text{ng}/\mu\text{L}$) + 41.0 μL of water added in other tube

2_2: 1.3 μL of 2_2 (794 $\text{ng}/\mu\text{L}$) + 42.7 μL of water added in other tube

2_3: 1.4 μL of 2_3 (751 $\text{ng}/\mu\text{L}$) + 42.6 μL of water added in other tube

3_1: 1.8 μL of 3_1 (585 $\text{ng}/\mu\text{L}$) + 42.2 μL of water added in other tube

3_2: 1.9 μL of 3_2 (527 $\text{ng}/\mu\text{L}$) + 42.1 μL of water added in other tube

3_3: 1.1 μL of 3_3 (951 $\text{ng}/\mu\text{L}$) + 42.9 μL of water added in other tube

4_1: 1.4 μL of 4_1 (743 $\text{ng}/\mu\text{L}$) + 42.6 μL of water added in other tube

4_2: 2.7 μL of 4_2 (373 $\text{ng}/\mu\text{L}$) + 41.3 μL of water added in other tube

4_3 :2.9 μL of 4_3 (353 $\text{ng}/\mu\text{L}$) + 41.1 μL of water added in other tube

2. added 5.00 μL of NEBuffer and 1.00 μL of Bsmbl in each tube

3. 2h, 37°C in Mastercycler

4. Agarosegel with EtBr (100 V, 45 min), see picture below

--> sent 3_3 to sequencing

--> not looking good

16.10.2019

1. picked colonies from retransformed plates **1_3** for overnight cultures
2. in incubator 37°C, 200 rpm

Plasmid DNA purification (Macherey-Nagel kit) of **1_4** and **1_5**:

see protocol; warmed Elution buffer to 80°C and just added 30 µL Elution buffer

-> no good concentrations

18.10.2019

Transformation:

1. 5 µL of the 2nd half of the GG lvi1 in DH5α
2. incubation on ice for 30 minutes
3. heatshock (42°C, 90 sec)
4. incubation in ice for 5 minutes
5. 500 µL of LB medium added
6. incubation (37°C, 1h)
7. centrifugation
8. remove supernatant
9. plate on agar plate with Spec
10. in incubator 37°C

19.10.2019

1. picked 3 white colonies of each plate for overnight cultures
2. in incubator 37°C, 200 rpm

20.10.2019

Plasmid DNA purification (Macherey-Nagel kit):

see protocoll; warmed Elution buffer to 80°C and just added 30 µL Elution buffer

concentration too low for digest/sequencing

Gel1bearb.PNG

