

Cell extraction with autolysis combined with bead-milling

This Method is used for cell extraction via autolysis and subsequent cell wall disruption with bead-milling.

Sample preparation

For autolysis, *C. vulgaris* and *C. sorokiniana* were prepared as shown in the following table, baker's yeast served as a control.

Table 1: List of samples and their pH for autolysis

Sample	Sample ID	pH-value	cell dry weight [mg]
baker's yeast	BY1	3	300
baker's yeast	BY2	3	300
baker's yeast	BY3	12	300
baker's yeast	BY4	12	300
<i>C. vulgaris</i>	CV1	3	300
<i>C. vulgaris</i>	CV2	3	300
<i>C. vulgaris</i>	CV3	12	300
<i>C. vulgaris</i>	CV4	12	300
<i>C. sorokiniana</i>	CS1	3	300
<i>C. sorokiniana</i>	CS2	3	300
<i>C. sorokiniana</i>	CS3	12	300
<i>C. sorokiniana</i>	CS4	12	300

For each organism two samples were used (double determination). To have the same amount of biomass in each sample cell dry weight was calculated from OD-CDW correlations. For *Chlorella vulgaris* we used our own OD-CDW correlation from dilution. For baker's yeast a correlation from A. Kusterer [1] was used. For *C. sorokiniana* we used the correlation from C.Y. Lim [2].

Autolysis protocol

1. The samples were centrifuged at 3000 rpm for 20 minutes. Supernatant was discarded.
2. Samples were resuspended in DSN medium.
3. pH-value of the samples was set to pH 3 with HCl or pH 12 with NaOH with a pH-Meter (digital pH-Meter wtw)
4. Samples were heated to 50 °C for 41 h and shook at 120 rpm (AquaShaker Kühner AG Schweiz)
5. Samples were neutralized with NaOH or HCl, respectively.

Cell Disruption

1. 3 mL of every sample was used for bead-milling.
2. Prepared tubes with beads from “quick DNA fungal/bacterial Mini-Prep Kit” by Zymo were used.
3. Samples were bead-milled with “retsch RKM Typ NN2” at maximum speed for 20 minutes.
4. To get rid of the beads, filters from “quick DNA fungal/bacterial Mini-Prep Kit” by Zymo were used. The samples were put on the filter and centrifuged at 11,000 rpm for 5 minutes. The flow through was kept for further analyses with the rFAN-Assay. The Assay was performed as described in the according protocol.

Literature

- [1] A. Kusterer, „Reaktionstechnische Optimierung von Parallelreaktoren für kontrollierte Bioprozesse“, Dissertation, TUM Munich, 2007.
- [2] Lim, et al. „A strategy for urban outdoor production of high-concentration algal biomass for green biorefining“, Bioresource Technology, v.135, 2013, pp. 175-181.