

### **Q5 Site-directed mutagenesis**

To insert the tags before the sfGFP a Q5 site directed mutagenesis was performed. Therefore, non-overlapping primers were designed using the NEBase Changer tool. Annealing temperatures were calculated using the NEB Tm Calculator. The reaction composition is as follows in table 1.

Table 1: Q5 Site-directed mutagenesis reaction composition.

<b>Component</b>	<b>50 <math>\mu</math>L Reaction</b>	<b>Final concentration</b>
5x Q5 Reaction buffer	10 $\mu$ L	1x
10 mM dNTPs	1 $\mu$ L	200 $\mu$ M
10 $\mu$ M Forward Primer	2.5 $\mu$ L	0.5 $\mu$ M
10 $\mu$ M Reverse Primer	2.5 $\mu$ L	0.5 $\mu$ M
Template DNA	variable	1 pg–10 ng
Q5 High-Fidelity DNA Polymerase	0.5 $\mu$ L	0.02 U/ $\mu$ L
Nuclease-Free Water	To 50 $\mu$ L	

Table 2: Thermocycling conditions for Q5 site-directed mutagenesis

<b>Step</b>	<b>Temperature</b>	<b>Time</b>
Initial Denaturation	98 °C	30 seconds
	98°C	10 seconds
25 cycles	50-72°C	30 seconds
	72°C	20-30 seconds/kb
Hold	4 °C	

After PCR, the product was phosphorylated, then ligated and finally digested with DpnI. The reaction composition for the non-radioactive phosphorylation is as follows. Since ATP is already present in the T4 ligase buffer, the same approach was used for both reactions. The T4 polynucleotide kinase exhibits 100% activity in this buffer.

Table 3: T4 PNK and ligase reaction composition.

<b>Component</b>	<b>50 <math>\mu</math>L Reaction</b>
T4 DNA Ligase Buffer (10x)	5 $\mu$ L
T4 PNK	1 $\mu$ L (10 units)
T4 Ligase	1 $\mu$ L (10 units)
DNA	Up to 300 pmol
Nuclease-Free Water	To 50 $\mu$ L

The approach was then digested with DpnI. 1  $\mu$ L DpnI was added to the entire reaction mixture and incubated overnight at 37°C to ensure complete digestion. Afterwards 1  $\mu$ L was transformed into competent Dh5a cells. After successful transformation, the DNA was purified and the mutagenesis was verified by Sanger sequencing.