



DYS SEE

RPA.



Protocols

RPA Protocol (Kit: TwistAmp Basic from TwistDX)

- This protocol is destined for single-plex use only.

Sample Set-Up (single-plex)

1. Prepare reaction mix in 1.5ml tube:

Primer A (10 μ M)	2,4 μ l
Primer B (10 μ M)	2,4 μ l
Primer Free Rehydration buffer	29,5 μ l
Template and water to	13,2 μ l
(Total volume	47,5 μ l).

Vortex and spin briefly.

2. Add reaction mix to a TwistAmp basic reaction. Pipette to mix.
3. Add 2,5 μ l of 280mM Magnesium Acetate (MgOAc) (supplied in the kit) and mix well to start reaction.
 - Note: RPA reaction starts as soon as MgOAc is added.
4. Incubate at 39 degrees Celsius for 20 minutes. For low template copy number, remove strip after 4 minutes, vortex and spin briefly, replace in heating device.
5. After 20 minutes, clean amplicons before running on agarose gels.
 - Note: If tubes are opened after amplification, there is a great risk of contamination of work surfaces with amplicon. Ensure that appropriate avoidance measures are taken.

Kit Positive Control Set-up (single-plex)

1. Prepare reaction mix in 1,5ml tube:

Positive control primer mix	8 μ l
Primer Free Rehydration buffer	29,5 μ l
Positive control DNA template	1 μ l
Water	9 μ l
(Total volume	47,5 μ l)

Vortex and spin briefly

2. Add reaction mix to a TwistAmp Basic reaction. Pipette to mix.
3. Add 2,5mM of Magnesium Acetate (MgOAc) (supplied in the kit) and mix well to start reaction.

- Note: RPA reaction starts as soon as MgOAc is added.
4. Incubate at 39 degrees Celsius for 20 minutes. For low template copy number, remove strip after 4 minutes, vortex and spin briefly, replace in heating device.
 5. After 20 minutes, clean amplicons before running on agarose gels.
- Note: If tubes are opened after amplification, there is a great risk of contamination of work surfaces with amplicon. Ensure that appropriate avoidance measures are taken.

Modified protocol to reduce reaction volume

To create 5 μ l final reaction volumes, single 50 μ l reactions should be prepared as described above, while being store on ice to prevent amplification. Aliquots (portions) of 5 μ l are then transferred into different tubes (@paper 200 μ l PCR tubes) and then placed into the incubator. (Lillis, et al. , 2016)



Protocols