

07. (July) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

TUESDAY, 30/7/2019

Transfection :

- 96 well plate for HiBit assay
- 6 well plate for Biotin purification
- 100µL medium out of the 150µL were exchanged in the 96 well plate. 1.5mL out of the 2mL medium were exchanged in the 6 well plate.
- time of transfection : 11:00 - 11:15

96 well plate for HiBit transfection scheme

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B		1	1	1	1	1	1	7	8	9	10	
C		2	2	2	2	2	2	7	8	9	10	
D		3	3	3	3	3	3	7	8	9	10	
E		4	4	4	4	4	4	7	8	9	10	
F		5	5	5	5	5	5	7	8	9	10	
G		6	6	6	6	6	6	7	8	9	10	
H												

6 well plate transfection s...

	1	2	3
A	2	2	5
B	4	4	5

conditions								
	V8	V9	V10	V11	V14	V15	V27	V30
1	/	/	/	/	/	/	/	100
2	40	/	/	/	/	/	20	40
3	/	40	/	/	/	/	20	40
4	40	/	/	20	20	/	20	/
5	/	40	/	20	20	/	20	/
6	40	/	20	/	20	/	20	/
7	40	/	20	/	/	20	20	/
8	/	40	20	/	/	20	20	/
9	40	/	/	20	/	20	20	/
10	see below							

-> condition 10 is 40ng V25 + 60nG V30

-> for the 6.5 wells these ng respond to :

- 20ng = 1.3µL
- 40nG = 2.6µL
- 60ng = 3.9µL
- 100nG = 6.5µL

Alejandro

cell culture: Transfection

- 550 µL medium was exchanged before transfection

Table1

	condition	V8	V11	V14	V27	V28	V30
1	Mock	-	-	-	-	-	500 ng
2	His	200 ng	100 ng	100 ng	-	100 ng	-
3	scAvidin	200 ng	100 ng	100 ng	100 ng	-	-


- Transfection:

Table2

	A	B
1	DNA per well	500 ng
2	P3000 Reagent per well	1 µL
3	Lipofectamine 3000 reagent per well	0.75 µL
4	OptiMEM per well	2 x 25 µL



- Transfect cells according to the following table. Use the indicated volume of DNA and P3000™ Reagent with each of the two volumes of Lipofectamine™ 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.



Timeline		Steps	Procedure Details (Two Reaction Optimization)			
Day 1	1	Seed cells to be 70-90% confluent at transfection	Component	10-well	24-well	96-well
			Adherent cells	1-4 × 10 ⁵	0.5-2 × 10 ⁵	0.25-1 × 10 ⁵
Day 1	2	Dilute Lipofectamine™ 3000 Reagent in Opti-MEM™ Medium (2 tubes) - Mix well	Opti-MEM™ Medium	5 µL × 2	25 µL × 2	125 µL × 2
			Lipofectamine™ 3000 Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
Day 1	3	Prepare master mix of DNA by diluting DNA in Opti-MEM™ Medium, then add P3000™ Reagent (2 µL/µg DNA) - Mix well	Opti-MEM™ Medium	10 µL	50 µL	250 µL
			DNA (0.5-5 µg/µL)	0.2 µg	1 µg	5 µg
Day 1	4	Add diluted DNA to each tube of diluted Lipofectamine™ 3000 Reagent (1:1 ratio)	P3000™ Reagent (2 µL/µg DNA)	0.4 µL	2 µL	10 µL
			Diluted DNA (with P3000™ Reagent)	5 µL	25 µL	125 µL
Day 1	5	Incubate	Diluted Lipofectamine™ 3000 Reagent	5 µL	25 µL	125 µL
			Incubate for 10-15 minutes at room temperature.			
Day 1	6	Add DNA-lipid complex to cells	Component	10-well	24-well	96-well
			DNA-lipid complex	10 µL	50 µL	250 µL
Day 1	7	Visualize/analyze transfected cells	DNA amount	100 ng	500 ng	2500 ng
			P3000™ Reagent	0.2 µL	1 µL	5 µL
			Lipofectamine™ 3000 Reagent used	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
Incubate cells for 2-4 days at 37°C. Then, analyze transfected cells.						

Moritz

cell culture: transfection

- note: cell density seemed low (40 %)
- every condition was transfected in 10 wells

Transfection scheme - qPCR 96-well plate 30/07/19 (ng per well)

	Condition	V8	V10	V11	V14	V15	V17	V
1	1	-	-	-	-	-	-	100 ng
2	4	40 ng	-	20 ng	20 ng	-	20 ng	-
3	7	40 ng	20 ng	-	-	20 ng	20 ng	-

- Transfection:

Transfectionmix 96-well plate 30/0...			A	B
1	DNA per well	100 ng		
2	P3000 Reagent per well	0.2 µL		
3	Lipofectamine 3000 reagent per well	0.15 µL		
4	OptiMEM per well	2 x 5 µL		

- Transfect cells according to the following table. Use the indicated volume of DNA and P3000™ Reagent with each of the two volumes of Lipofectamine™ 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.

image.png

Timeline		Steps	Procedure Details (Two Reaction Optimization)			
Day 1	1	Seed cells in to 70-90% confluent at 4 reactions	Cell suspension	5.0 ml	2.0 ml	0.25-1 x 10 ⁶
	2	Adherent cells	Opti-MEM® Medium	5 µL x 2	25 µL x 2	125 µL x 2
Day 2	3	Dilute Lipofectamine™ 2000 Reagent in Opti-MEM® Medium (2 Labeled- No well)	Lipofectamine™ 2000 Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 µL
	4	Prepare master mix of DNA to diluting DNA in Opti-MEM® Medium, then add P3000® Reagent- No well	Opti-MEM® Medium	10 µL	50 µL	250 µL
Day 3	5	Add Diluted DNA to each tube (2 Labeled- Lipofectamine™ 2000 Reagent 15 Labeled)	DNA (5-5 µg/µL)	0.2 µg	1 µg	5 µg
	6	Incubate	P3000® Reagent (2 µL/µg DNA)	0.4 µL	2 µL	10 µL
Day 4	7	Add DNA-lipid complex to cells	Diluted DNA with P3000® Reagent	5 µL	25 µL	125 µL
		Visualize/analyze transfected cells	Diluted Lipofectamine™ 2000 Reagent	5 µL	25 µL	125 µL