Growth curve of Vibrio natriegens

To determine the growth parameters of wild type *Vibrio natriegens* in complex medium a growth curve was generated.

- 1. Select three colonies from fresh agar plate and inoculate a 5 mL starter culture of BHIN. Grow culture at 37 °C and 150 rpm overnight.
- 2. Inoculate 3x 150 mL BHIN (500 mL shaking flask) with preculture to an initial OD_{600} of 0.1
- 3. Cultivate Vibrio natriegens at 37 °C and 150 rpm
- 4. Measure OD_{600} every 15 minutes in triplicates for at least 3 hours
- 5. Plot OD_{600} values against time to determine the growth rate μ [1/h] and doubling time [h]