

## **Quantification of algae cell disruption: Anthrone assay**

### **Principle**

Carbohydrates condense with anthrone to form a green color complex after dehydration with concentrated  $\text{H}_2\text{SO}_4$  to form "Furfural". This can be measured by using colorimetrically at 620nm. Anthrone reacts with dextrans, monosaccharides, disaccharides, polysaccharides, starch, gums, and glycosides.

### **Reagents**

1. Anthrone solution: 200% Anthrone reagent dissolved in concentrated  $\text{H}_2\text{SO}_4$  (w/v), e.g. 200 mg Anthrone reagent in 100 ml concentrated  $\text{H}_2\text{SO}_4$
2. Standard Glucose solution: 10% Glucose solution (w/v) gradually diluted 1:5 with 1ml end volume per sample
3. Extracted Carbohydrates in neutral solution, possibly without high salt contents with 1ml end volume per sample

### **Procedure**

1. Add 4 ml of anthrone solution to the working standard samples, respectively. Mix thoroughly .
2. Add 4 ml of anthrone solution to the samples. Mix thoroughly.
3. Heat mixture to 100 °C for 10 minutes. Shock cool the samples in ice bath immediately.
4. Measure colorimetrically at 620nm.
5. Plot standard curve and extract the equation to calculate the amount of carbohydrates in your samples.