

# 09. (September) 2019

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Theresa Keil

FRIDAY, 6/9/2019

- continued with WB
- washed 3x10min in TBST
- 2h at RT in sec AB diluted
  - goat-anti-mouse 1:20000 in 5% BSA
  - goat-anti-rat 1:10000 in 5% BSA
 dilution to try to reduce background

new WB VLP-anti-MCP

- harvested samples like VLP harvest protocol
- lysed all samples dilution 1:1
- diluted 1:1 with Laemmli, incubated 10 min at 95 °C
- SDS-Page and WB after protocol

**Table3**

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	marker	Mock SNL	mock CL	V4 SNL	V4 CL	MCP-MS2 3 SNL	MCP-MS2 3 CL	3* SNL	3* CL	4 SNL	4 CL	5 SNL	5 CL	6 SNL	6 CL

VLP samples received

750 µl supernatant and all cells are used for lysis and DNase treatment

RNA extraction (Exosomes, VLPs and correspondin cells)

the standard trizol protocol

- 750 µl SN were mixed with 750 µl Trizol
- after addition of chloroform, the eppi was shaken heavily for 3 min and the incubated at RT for 5 min
- after centrifugation only 700 µl aquaous phase were transfered to new eppi
- 800 µl Isopropanol were added
- 30 min incubation at - 20 °C
- all 20 µl SN were discarded and the pellet + 20 µl SN were incubated at 60 °C for 20 min
- further 20 µl RNA se-free water were given to the samples and incubated again at 60 °C for 10 min
- the samples were frozen at - 80 °C via liquid N2 freeing