

# 06. (June) 2019

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Johanna Wallner

FRIDAY, 28/6/2019

## Sequencing results :

- V10 : #3,#5,#8 were positive -> continue with V10#5
- 500mL LB medium + 500µL Amp

## restriction :

- Insert DNA creation from primers 22 and 23 / primers 32 and 33 according to standard restriction protocol
- restriction digest of V10#5 (backbone) and V4 (insert)

scheme			
	A	B	C
1	DNA	0.5µL	3.3µL
2	10xCutSmartB uffer	5µL	5µL
3	NheI-HF	1µL	1µL
4	Mlu-HF	1µL	1µL
5	MilliQ	42.5µL	39.7µl

- incubate at 37°C for 15min
- both samples were mixed with 15 µL DNA loading dye and put on an electrophoresis gel for 20min
- V4 was discarded, because there were 2 bands visible and we already amplified the insert via PCR

## ligation :

- T4 Buffer 10x : 2µL
- T4 Ligase : 1µL
- Vector DNA (V10) : 50ng
- Insert DNA (V4) : 37.5ng
- nuclease free water : 17.04µL

## restriction digest :

- DNA (Nhe-MCP-Mlu1) : 7.3µL
- 10x Cut Smart Buffer : 5µL
- Nhe1-HF : 1µL
- Mlu1-Hf : 1µL
- MilliQ : 35.7µL

## purification :

- cleanup of restriction digest using NEB Monarch PCR DNA cleanup kit
- ligation on digested V10 backbone and MC-fragment to form V11 (see above)
  - incubate for 10min at room temperature
  - heat inactivation at 65°C for 10min

transformation :

- according to standard transformation protocol started with overnight culture

cell culture :

- splitting of motherplates
- 6 wells of a 6 well plate a 500.000 cells/well
- 30 wells of a 96 well plate a 20.000 cells per well
- new DMEM medium prepared