



UPDATED 9.10.2020 16:26



by Franka Butzbach (created by Philip Schulz)



Introduction

PCR product purification with Omega Cycle Pure Kit

Materials

- › centrifuge tubes 1.5μL
- › cycle pure kit
- › MQ

Procedure

Procedure

1. calculate how much volume your sample has: eg from PCR amplification microstrip 25 μ L - 5(that we run a GEL with) = 20 μ L sample volume
2. transfer this volume of PCR product in a new 1.5 μ L centrifuge tube
3. add 4-5 volumes of CP Buffer to each tube (in our case= 100 μ L)
4. vortex roughly and centrifuge so that everything is mixed well and at the bottom
5. insert a HliBind minicolumn into the provided 2 μ L collection tubes and transfer everything in the HliBind column-2 μ L tube
6. centrifuge at max. speed of 13'000g for 60 sec. Discard the filtrate and reuse the collection tube.
7. add 700 μ L DNA wash buffer that is diluted with 100 % ethanol (green bottle: if crossed et already added) centrifuge at max. speed for 60 sec. Discard the filtrate and reuse the collection tube
8. repeat this step 7. for a second DNA wash buffer step!
9. centrifuge the empty Hibind coloumn at max speed for 2 min. to dry the column
10. transfer the empty Hibind column to a new clean 1.5 μ L tube
11. Add 30-50 μ L Elution buffer or MQ directly to the centre of the coloumn. let it sit at room temp for 2 min. centrifuge at max speed for 2 min.
12. take out the column, store DNA at -20°C