





DH5α Mix & Go Cell Transformation





Introduction

This protocol transforms plasmid DNA into DH5α cells.

Reagents

-  Plasmid DNA
-  New England Biolabs Mix & Go Cells
-  SOB media
-  Luria Broth/Kanamycin plates

Equipment

-  Parafilm
-  Pipette
-  Centrifuge
-  Ice and bucket
-  Incubation cabinet

Procedure

1. Add 1-5 μL of plasmid DNA to a tube of thawed Mix & Go cells on ice, mixing gently for a few seconds.
 - a. *Try to keep the added volume of DNA less than 5% of the total.*
2. After the transformation mixture has incubated on ice for 5-10 minutes, add 4 volumes of SOC (400 μL of SOB:100 μL of transformation mixture) and incubate for 1 hour at 37°C with gentle shaking at 200-300 rpm.
3. Spread 50-100 μL of the mixture onto a pre-warmed (37°C) culture plate containing LB/Kanamycin.
4. Parafilm the edges of the plate and incubate at 37°C for the colonies to grow.