



GFP_KanR (823 bp)





Reception and resuspension

WEDNESDAY, 7/29/2020

IDT sequence has been received like 6 other sequences.
Resuspension finale à une concentration de 10ng/ul.

GG#3

WEDNESDAY, 7/29/2020

Done by Florent and Mathieu

Note 1 : You need to extract and quantify with a nanodrop the plasmid before to begin the GG (here we used pBAD-Moclo serie 1 with a concentration of 77 ng/µL)

Note 2 : Note : this manipulation was done in the same time with kanR_GFP

Material :

- IDT sequence insert ilvA_ccdB resuspended at 10 ng/µL (1 000 ng of sequence lyophylized = powder was given by IDT)
- thermocycleur
- H₂O mQ
- 75 ng of plasmid pBAD-MOCLO
- kit Golden Gate NEB

Protocol for cloning golden gate for pBADMoclo:

Add in this order :

1. 15 µL of inserts with ratio 2:1 (insert:plasmid) (10 ng/µL) (H₂O for control condition)
2. 2 µL of T4 DNA ligase buffer
3. 1 µL of plasmids pBADMoclo (75 ng/µL)
4. 1 µL of QSP H₂O
5. 1µL of NEB Golden Gate

NB: DONT FORGET TO PUT THE TUBES IN THE ICE AND DONT FORGET TO HAVE A CONTROL CONDITION

Table1

	A	B	C	D	E	F
1		plasmide de destination :	Insert with amplicon form (150 ng) :	T4 DNA ligase Buffer	NEB Golden Gate Enzyme mix	H ₂ O
2	Volume (µL) for Golden Gate cloning	1 µL	15 µL (C = 10 ng/µL)	2 1-2		qsp 20 µL

We prepared 3 GG reactionnal tubes :

- pBAD-MOCLO + knt_ccdB_gene272_Bsal
- pBAD-MOCLO + knt_GFP
- pBAD-MOCLO only

(Heat at 37°C for 1 min and then 16°C for 1min) for 30 cycles and 60°C for 5 min

After the GG cycle, we transformed bacteria DH5alpha.

Table2



	A	B	C	D	E	F
1	GG#X	Plasmids	Number of colonies			Middles
2			50 µL	100 µL	Concentrate	
3	GG#3	pBad-MC_sfGFP_knt	0	0	0	LB + AMP + glucose
4		pBAD-MC_ccdB_knt	0	0	0	
5		control H2O (sans insert)	0	0	0	

Transformation 2 GG#3

THURSDAY, 7/30/2020

Seconde transformation based on GG#3 simultaneously to PCR #5 to determine if transformation is effective if PCR results are those expected (it's not the case)

Done by Florent.

We started by transforming chemocompetent DH5alpha.

Protocol of transformation [Inaccessible Protocol](#)

We than transformed the golden gate products :

- GG#3 reaction tube pBAD-MC knt-GFP
- GG#3 reaction tube pBAD-MC knt-ccdB
- GG#3 reaction tube pBAD digéré sans insert

[Inaccessible Entry](#)

PCR#5 from GG#3 and GG#4

THURSDAY, 7/30/2020

Done by Mathieu

We had no transformants with GG#3 and GG#4 so we did a PCR and an electrophoresis migration to verify if it's our Golden Gate or our transformation method that don't function.

Material:

- 2019GO-5-Seq-pBAD-F
- 2019GO-6-Seq-pBAD-R
- oligo EC632 (forward)
- oligo EC633 (reverse)
 - -> NB : 1 mutation between oligos and plasmid sequence for F and R.

Table1

	A	B	C	D	E
1	Nom de l'oligo	Cible	Séquence (5'-3')	Tm (°C, IDT-oligo analyzer)	hairpin Tm (°C, IDT)
2	2019GO-5-Seq-pBAD-F	pBAD24	5'- TGGCTATGCCATAGCATTAT CC	55,6	50,0
3	2019GO-6-Seq-pBAD-R	pBAD24	5'- GGT CTG ATT TAA TCT GTA TCA GGC TG	55,3	38,7
4	oligo EC632	pZE12	5' gcccttcgttccacctc		
5	oligo EC633	pZE12	tcaaaggcgtaatacggtt		

- results of GG#3 and GG#4
- PCR tube
- H2O milli Q
- Thermocycler
- toothpick
- ice

Protocol

1. Add composants in this order (cf below). Attention, all this tubes are in ice

We use the kit " Q5 High-Fidelity 2 X Master MIX" from New England Biolabs.

Table2

	A	B	C	D
1	Order of introduction	Components	Volume (μ L) for 50 μ L reaction	Final concentration
2	1)	Template DNA (pBAD-MOCLO +/- insert)	1 (0 for the negative control)	<1.000 ng
3	2)	Q5-High_Fidelity 2X Master Mix	25	0.5 μ M
4	3)	Oligo Forward	2,5	0.5 μ M
5	4)	Oligo Reverse	2,5	1 X
6	5)	Nuclease free water	19 (20 for the negative control)	/

We prepare 7 tubes with the mix PCR

- GG#3 reaction tube pBAD-MC knt-GFP
- GG#3 reaction tube pBAD-MC knt-ccdB
- GG#3 reaction tube pBAD digéré sans insert
- GG#4 reaction tube pZE12 + Luc_consensus_Nter_Bsal
- GG#4 reaction tube Luc_GFP_gene 28
- GG#4 reaction tube Luc_knt_gene 36
- GG#4 reaction tube pZE12 + GFP_Luc1
- GG#4 reaction tube pZE12 + GFP_Luc2
- pBAD natif
- pZE12 natif
- pZE12 digéré par DPN1
- GG#4 reaction tube pZE12 digéré sans insert

PCR Program:

Determination of the Temperature of elongation :

size of pBAD-MOCLO : 45 110 pb

size of ilvA_ccdB : 1 582 pb

size of stable_YFP : 367 pb

$$4\ 510 + 1\ 582 = 6\ 092 \text{ pb}$$

Processivité de la polymérase : 20-30 seconds / kb

6kb X 30 second = 180 seconds

1. 98°C - 2min
2. 98°C - 10 s
3. 52°C - 30 s
4. 72°C - **60** s
5. Repeat from step 2 (5 times)
6. 98°C - 10 s
7. 52°C - 30 s
8. 72°C - **60** s
9. Repeat from step 6 (**20** times)
10. 72°C - 2 min
11. 4°C - Infinity

N-B : We changed the time for the temperature of elongation : in the protocol of PCR #1, the time was 3 min and in this protocol it's 1 min.

The products are analysed by electrophoresis

Preparation of electrophoresis gel (cf file "Electrophoresis pBAD- MOCLO")

Electrophoresis of PCR#5 products

THURSDAY, 7/30/2020

Done by Mathieu

Protocol :

NB: Don't forget your gloves for the heat AND the Ethidium Bromide (BET) !!!

For the gel

- Weigh 4,8g of agarose (0,8% --> 0,8g pour 100 ml)
- Measure 600 ml of TAE buffer (0.5X) in a test tube
- Put the agarose and then add the buffert in a 1L (glass) bottle.
- Keep the gel boiling in the microwave
- Cast the gel and put the combs

For 1 well :

- on 1 PCR tube, put : 10µL of mix PCR + 2µL of blue charge. Then put 12µL of this mix in 1 well
- for the ladder : just put 5 µL of ladder

For the migration :

- Set the machine on 100 Volts during **25 minutes**.

Electrophoresis on :

We prepare 7 tubes with the mix PCR


- GG#3 reaction tube pBAD-MC knt-GFP
- GG#3 reaction tube pBAD-MC knt-ccdB
- GG#3 reaction tube pBAD digéré sans insert
- GG#4 reaction tube pZE12 + Luc_consensus_Nter_Bsal
- GG#4 reaction tube Luc_GFP_gene 28
- GG#4 reaction tube Luc_knt_gene 36
- GG#4 reaction tube pZE12 + GFP_Luc1
- GG#4 reaction tube pZE12 + GFP_Luc2
- pBAD natif
- pZE12 natif
- pZE12 digéré par DPN1
- GG#4 reaction tube pZE12 digéré sans insert

Depot map's + Results

1. Ladder
2. GG#3 reaction tube pBAD-MC knt-GFP
3. GG#3 reaction tube pBAD-MC knt-ccdB
4. GG#3 reaction tube pBAD digéré sans insert
5. Ladder
6. GG#4 reaction tube pZE12 + Luc_consensus_Nter_Bsal
7. GG#4 reaction tube Luc_GFP_gene 28
8. GG#4 reaction tube Luc_knt_gene 36
9. Ladder
10. GG#4 reaction tube pZE12 + GFP_Luc1
11. GG#4 reaction tube pZE12 + GFP_Luc2
12. pBAD natif
13. Ladder


14. pZE12 natif
15. pZE12 digéré par DPN1
16. GG#4 reaction tube pZE12 digérée sans insert
17. Ladder

Results



Le Golden Gate n'a pas fonctionné, seuls les plasmides natifs semblent avoir été amplifiés par PCR.

 image.png



GG#7

WEDNESDAY, 8/5/2020

Done by Goddefroy

New protocol !!!!!!!!

Material :

- IDT sequenc inserts resuspended at 10 ng/ μ L (1 000 ng of sequence lyophylized = powder was given by IDT)
- thermocycleur
- H₂O mQ
- plasmids (pBAD-Mc or pZ)
- kit Golden Gate NEB

Protocol modifications :

- T4 DNA ligase : 1000U become 500U
- Nb of cycles is variable (15min 37°C, 30 cycles, 45 cycles, OV (overnight))
- insert : 100 ng become 150 ng
- Bsal : 30U become 15U

Protocol for cloning golden gate for pBADMoclo:

Add in this order :

1. 15 μ L of inserts with ratio 1:1 (insert:plasmid) (10 ng/ μ L) (H₂O for control condition)
2. 2,5 μ L of T4 DNA ligase buffer
3. 4uL of plasmid
4. 0,25 μ L of T4 DNA ligase
5. Bsal-HF-V2 (20U/uL) : 0,75uL
6. H₂O MQ : 3uL

NB: DONT FORGET TO PUT THE TUBES IN THE ICE AND DONT FORGET TO HAVE A CONTROL CONDITION

Table1								
	A	B	C	D	E	F	G	H
1		plasmide de destination :	Insert with amplicon form (150 ng) :	T4 DNA ligase Buffer	Bsal-HF-V2 (20U/uL)	T4 DNA ligase	H ₂ O	
2	Volume (μ L) for Golden Gate cloning	75ng	15 μ L (C = 10 ng/ μ L)	2,5 uL	0,75uL	0.25uL	3 uL	

We prepared 10 GG reactionnal tubes :

- pBAD-MC + YFP stable 15min 37°C
- pBAD-MC + YFP stable 30 cycles
- pBAD-MC + knt_GFP 15min 37°C
- pBAD-MC + knt_GFP 30 cycles
- pBAD-MC + knt_GFP45 cycles
- pBAD-MC + knt_GFP OV
- pBAD-MC digested without insert 15min 37°C
- pBAD-MC digested without insert 30 cycles
- pBAD-MC digested without insert 45 cycles
- pBAD-MC digested without insert OV

After the GG cycle, we transformed bacteria DH5alpha for 15min and 30 cycles conditions, no for the others.

 Inaccessible Entry

Table2

	A	B	C	D	E	F
1	GG#X	Plasmids	Number of colonies			Middles
2			50 µL	100 µL	Concentrate	
3	GG#7	pBAD-MC_YFP stable 15'	3	20	98	LB + AMP + glucose
4		pBAD-MC_YFP stable 1h30	0	5	38	
5		pBad-MC_sfgfp_knt 15'	3	22	80	
6		pBad-MC_sfgfp_knt 1h30	3	8	37	
7						
8						
9		pBad-MC digested alone 15'	4	5	88	
10		pBad-MC digested alone 1h30	0	0	11	
11						
12						
13						

PCR colony #7 from GG#7

WEDNESDAY, 8/5/2020

Done by Godeffroy

Material:

- 2019GO-5-Seq-pBAD-F
- 2019GO-6-Seq-pBAD-R

Table1					
	A	B	C	D	E
1	Nom de l'oligo	Cible	Séquence (5'-3')	Tm (°C, IDT-oligo analyzer)	hairpin Tm (°C, IDT)
2	2019GO-5-Seq-pBAD-F	pBAD24	5'- TGGCTATGCCATAGCATTAT CC	55,6	50,0
3	2019GO-6-Seq-pBAD-R	pBAD24	5'- GGT CTG ATT TAA TCT GTA TCA GGC TG	55,3	38,7

- colony GG#7 knt_GFP 15'
- colony GG#7 knt_GFP 1h
- colony GG#7 YFP 154
- colony GG#7 YFP 15'
- colony GG#7 H2O
- colony GG#7 native colony

Protocol

1. With a toothpick, stripe 1 colony on 1 cm; then, put it on an eppendorf containing 100 µL of H2O. Homogenise the tube. You have the 'Template DNA' ready. NOTE : in PCR#1 and 2, I pipette 1µL of this 100µL H2O + plasmid as DNA Template. I was wrong You need to pipette 19+1 µL H2O + plasmid.
2. Add composants in this order (cf below). Attention, all this tubes are in ice

We use the kit " Q5 High-Fidelity 2 X Master MIX" from New England Biolabs.

Table2

	A	B	C	D
1	Order of introduction	Components	Volume (μ L) for 50 μ L reaction	Final concentration
2	1)	Template DNA (pBAD-MOCLO +/- insert)	1 (0 for the negative control)	<1.000 ng
3	2)	Q5-High_Fidelity 2X Master Mix	25	0.5 μ M
4	3)	2019GO-5-Seq-pBAD-F	2,5	0.5 μ M
5	4)	2019GO-6-Seq-pBAD-R	2,5	1 X
6	5)	Nuclease free water	19 (20 for the negative control)	/

Plan de dépôt voir photos

PCR Program:

98° 10 sec

25 times

1. 55° 30 sec
2. 72° 1 min

78° min

Sub-cloning transformation GG#7

THURSDAY, 8/6/2020

Stries with transformants from [Inaccessible Entry](#) (transformants YFP and knt_GFP) on middle LB + Amp to obtain isolated transformants which will be used for selection characterisation , sequencing and stocking.

Culture of colonies of sub-cloning GG#7 on selective middle + stocking + preparation tubes for PCR#9

FRIDAY, 8/7/2020

Done by Florent

Voir  Inaccessible Entry

- Sub-clones knt-GFP were put on LB+AMP, LB +AMP+SPE 125uL et LB+AMP+SPE 250 uL.
- Sub-clones YFP were put on LB+AMP

PCR colony #10

MONDAY, 8/10/2020

PCR colony done with knt_GFP sub-clones from GG#7 grew on LB+amp+spe middle.

Done by Godeffroy

Material:

- 2019GO-5-Seq-pBAD-F
- 2019GO-6-Seq-pBAD-R

Table1					
	A	B	C	D	E
1	Nom de l'oligo	Cible	Séquence (5'-3')	Tm (°C, IDT-oligo analyzer)	hairpin Tm (°C, IDT)
2	2019GO-5-Seq-pBAD-F	pBAD24	5'- TGGCTATGCCATAGCATTAT CC	55,6	50,0
3	2019GO-6-Seq-pBAD-R	pBAD24	5'- GGT CTG ATT TAA TCT GTA TCA GGC TG	55,3	38,7

- colony sub-clone knt_GFP GG#7 on box LB+amp (noir)
- colony sub-clone knt_GFP GG#7 on box LB+amp+spe (rouge)

Protocol

1. With a toothpick, stripe 1 colony on 1 cm; then, put it on an eppendorf containing 100 µL of H2O. Homogenise the tube. You have the 'Template DNA' ready. NOTE : in PCR#1 and 2, I pipette 1µL of this 100µL H2O + plasmid as DNA Template. I was wrong You need to pipette 19+1 µL H2O + plasmid.
2. Add composants in this order (cf below). Attention, all this tubes are in ice

PCR Program:

98° 10 sec

25 times

1. 55° 30 sec
2. 72° 1 min

78° min


We use the kit " Q5 High-Fidelity 2 X Master MIX" from New England Biolabs.

Table2

	A	B	C	D
1	Order of introduction	Components	Volume (μ L) for 50 μ L reaction	Final concentration
2	1)	Template DNA (pBAD-MOCLO +/- insert)	1 (0 for the negative control)	<1.000 ng
3	2)	Dream taq Buffer	0,5	0.5 μ M
4	3)	2019GO-5-Seq-pBAD-F	0,5	0.5 μ M
5	4)	2019GO-6-Seq-pBAD-R	0,5	1 X
6	5)	Nuclease free water	19 (20 for the negative control)	/
7	6)	Dream taq polymerase	0,125	

Plan de dépôt voir photos


image.png



We can see a bande at the good length so it's possibly good
we have an insert with the good length.

pBAD+KNT-GFP sous - clone GG#⁷
pBAD+KNT-GFP (spectino) —————

 image.png



Electrophoresis PCR colony #10 from sub-clones knt-GFP from GG#7

MONDAY, 8/10/2020

Done by Godeffroy

Protocol :

For the gel

- Weigh 4,8g of agarose (0,8% --> 0,8g pour 100 ml)
- Measure 600 ml of TAE buffer (0.5X) in a test tube
- Put the agarose and then add the buffert in a 1L (glass) bottle.
- Keep the gel boiling in the microwave
- Cast the gel and put the combs

For 1 well :

- on 1 PCR tube, put : 10µL of mix PCR + 2µL of blue charge. Then put 12µL of this mix in 1 well
- for the ladder : just put 5 µL of ladder

For the migration :

- Set the machine on 100 Volts during **25 minutes**.

Electrophoresis on :


- colony sub-clone knt_GFP GG#7 sur boite LB+amp (black)
- colony sub-clone knt_GFP GG#7 sur boite LB+amp+spe (red)

Depot map's + Results

See photos


Results

image.png



The bands are at the good lenght (eviron 1kb) so it 's probably our insert but bands very weak, probably a TM (temperature melting) problem for the PCR. **But there were no arabinose on the middle so finally it's not our insert.**

 image.png



Electrophoresis PCR colony #10 from sous clones knt-GFP from GG#7 (2)

TUESDAY, 8/11/2020

Done by Florent

Protocol :

For the gel

- Weigh 4,8g of agarose (0,8% --> 0,8g pour 100 ml)
- Measure 600 ml of TAE buffer (0.5X) in a test tube
- Put the agarose and then add the buffert in a 1L (glass) bottle.
- Keep the gel boiling in the microwave
- Cast the gel and put the combs

For 1 well :

- on 1 PCR tube, put : 10µL of mix PCR + 2µL of blue charge. Then put 12µL of this mix in 1 well
- for the ladder : just put 5 µL of ladder

For the migration :

- Set the machine on 100 Volts during 25 minutes.


Electrophoresis on :

- colony sub-clone knt_GFP GG#7 sur boite LB+amp (black)
- colony sub-clone knt_GFP GG#7 sur boite LB+amp+spe (red)

Depot map's + Results


See photos

Results




The bands are at the good lenght (environ 1kb) so it 's probably our insert but bands very weak, probably a TM (temperature melting) problem for the PCR. **But there were no arabinose on the middle so finally it's not our insert.**

📎 image.png



📎 image.png



Observation of knt-GFP and YFP strains with UV table

TUESDAY, 8/11/2020

Done by Mathieu

knt-GFP GG#7 amp (1, 2, 3, 4) : no fluorescence

knt-GFP GG#7 sous clone amp + spectino (1, 4) : no fluorescence

YFP GG#7 (1, 2, 3, 4) : no fluorescence

Miniprep pBAD knt_GFP sub-clone GG#7

TUESDAY, 8/11/2020

Done by Baptiste



To have enough DNA to sequence it.

Material :

- Monarch^R Plasmid Miniprep Kit (#T10105 lot: 10040564)

Protocol :

With 5 minutes of incubation and H2O for eluant.



Résultat :

Obtention of a pBAD-MC tube of concentration of 547 ng/uL

Measurement of DO of sous-clones from GG#7 and GG#8 grew in liquid middle by night

WEDNESDAY, 8/12/2020

Done by Baptiste

Mesurement of DO of sous-clones from GG#7 and GG#8 grew in liquid middle by night before observe them with microscope to mesure their fluorescence.

Results :

- pBAD-MC H2O 1: 0.91 1
- pBAD-MC H2O 2 : 0.85 2
- pBAD-MC H2O 3 : 0.75 3
- sub-clone pBAD-MC knt_GFP from GG#8 1: 0.50 4
- sub-clone pBAD-MC knt_GFP from GG#8 2 : 0.48 5
- sub-clone pBAD-MC knt_GFP from GG#8 3 : 0.49 6
- sub-clone pBAD-MC YFP from GG#7 1: 0.50 7
- sub-clone pBAD-MC YFP from GG#7 2 : 0.50 8
- sub-clone pBAD-MC YFP from GG#7 3 : 0.52 9

Culture by night of colonies knt-GFP and YFP for microscopique fluorescence observation and Clariostar observation + results observations

WEDNESDAY, 8/12/2020

Done by Baptiste

Initial OD of our culture was 0,1.

On middle liquid LB + AMP + Ara

-pBAD-MC H2O :

- col 1 : 549 uL
- col 2 : 588 uL
- col 3 : 666 uL

-pBAD-MC knt_GFP :

- col 4 : 1mL
- col 5: 1041 uL
- col 6: 1020 uL

-pBAD-MC YFP :

- col 7 : 1mL
- col 8 : 1041 uL
- col 9 961 uL

On middle liquid LB + chloro + Ara

- Witness which express GFP constitutively

Results :

 Inaccessible Entry

Cultures of knt-GFP and YFP for fluorescence observation on UV table in LB

WEDNESDAY, 8/12/2020

Incubation overnight

- pBAD-MC H2O (col 1, 2, 3)
- pBAD-MC knt_GFP (col 4, 5, 6)
- pBAD-MC YFP (col 7, 8, 9)
- Witness which express GFP constitutively 10)
- Witness YFP (11)


Results



All the colonies seem to express fluorescence even negative witness

Sending DNA from miniprep colony knt-GFP 4 from GG#8 for sequencing

WEDNESDAY, 8/12/2020



Fluorescence microscope observation 1

WEDNESDAY, 8/12/2020

Culture of colonies knt-GFP on selective middle

THURSDAY, 8/13/2020

Objective: Culture of pBAD-MC_knt_GFP colonies on solid middles to show spectinomycin resistance.

Protocol:

Culture of 4 colonies pBAD-MC_knt_GFP and pBAD-MC with no insert on three middles :

- LB + Amp + Ara + Spe (50 ng/mL)
- LB + Amp + Ara + Spe (100 ng/mL)
- LB + Amp + Ara


NB : Arabinose (Ara) is necessary because MC is an inducible vector

Results:



No colony grew on the middles with spectinomycin so there is no resistance given by our entangled gene.

image.png



Culture by night of colonies knt-GFP and YFP for microscopique fluorescence observation

THURSDAY, 8/13/2020

Same protocol and same colonies as Inaccessible Entry

Done by Baptiste

Initial OD of our culture was 0,1.

On middle liquid LB + AMP + Ara

-pBAD-MC H2O :

- col 1 : 549 uL
- col 2 : 588 uL
- col 3 : 666 uL

-pBAD-MC knt_GFP :

- col 4 : 1mL
- col 5: 1041 uL
- col 6: 1020 uL

-pBAD-MC YFP :

- col 7 : 1mL
- col 8 : 1041 uL
- col 9 961 uL

On middle liquid LB + chloro + Ara

- col 11 : Witness which express GFP constitutively
- col 10 : Witness which express YFP constitutively (not stable YFP)

FRIDAY, 8/14/2020

Results :

1 : 0,34

2 : 0,35

3 : 0,34

4 : 0,34

5 : 0,33

6 : 0,37

7 : 0,24

8 : 0,27

9 : 0,28

10 : 0,09

11 : 0,45

NB : the cultures have been diluted by factor 10 because DO was superior to 1,5

Inaccessible Entry

Results microscope

Culture of knt-GFP colonies one spectinomycin gradient

THURSDAY, 8/13/2020

Middle = LB + Arabinose + amp + spe gradient

Résults


Transformants and witnesses seem to present no resistance to spectinomycine.



Gel 2 UV



Gel 2



Clariostar knt-GFP and YFP

FRIDAY, 8/14/2020

Goal

Use Clariostar to measure fluorescence of our transformed bacteria.

Plate plan :

Table1								
	A	B	C	D	E	F	G	H
1	blanc 1	blanc 2	1	2	3	4	5	6
2	7	8	9	10	11			

blanc 1 : LB liquid + Amp + Arabinose

blanc 2 : LB liquid + Cmp + Arabinose

Numbers 1 to 9 refer to colonies of Inaccessible Entry

10 : YFP witness

11 : GFP witness

NB : our YFP witnesses do not express the same YFP as our transformed colonies which express our stable YFP.

Results

 [Clariostar_YFP_14082020_3_valeurs.xlsx](#)

 [Clariostar_OD600_14082020_3_valeurs.xlsx](#)

 [Clariostar_GFP_14082020_3_valeurs.xlsx](#)


The results seem to show that there is no fluorescence of our bacteria due to our biobricks knt-GFP and YFP.

Fluorescent microscope observation 2

FRIDAY, 8/14/2020

It seems that there is was no clear fluorescence of our bacteria expressing knt-GFP compared to the negative and positive controls.

Negative witness :




Positive witness :



Knt-GFP :

 DH5alpha + pBAD-KntGFP c11 image3_w2GFP.png



Cytometry YFP and knt-GFP

WEDNESDAY, 9/2/2020

Cytometry test to determine if our transformants of stable YFP and knt-GFP express fluorescence.


Results :

 Résultats cytométrie.pptx

We observe a clear fluorescence due to YFP but results for knt-GFP seem not clear, another cytometry experiment is required.


Negative control
(empty vector)

Mean 0.63




Stable YFP - LOV

Mean 1.83




Negative control
(empty vector)

Mean 0.63




Knt-GFP

Mean 0.87




Negative control
(empty vector)

Mean 0.63




Stable YFP - LOV

Mean 1.83




Negative control
(empty vector)

Mean 0.63




Positive control
YFP (pTH3)


Mean 50.54




Témoin *E. coli* vecteur vide
+ arabinose




Témoin *E. coli* YFPstable-LOV
+ arabinose



Témoin *E. coli* knt-GFP
+ arabinose




Témoin *E. coli* GFP-uv
+ arabinose



Témoin *E. coli* pTH3 (YFP)

+ arabinose



Cytometry YFP and knt-GFP 2

THURSDAY, 9/3/2020

The results of the first cytometry were not sufficient to determine if knt-GFP permit fluorescence so we did a second cytometry test.

Results :

 iGEM-Cytometry 2_data_analyses.docx

It seems that transformants of stable YFP express fluorescence compared to negative control but less fluorescence than positive control.

There is no clear proof fluorescence due to knt-GFP.

iGEM2020 Cytometry data


Fluorescence detection:

Laser blue (excitation): 488 nm (20 mW)

Emission: 536 nm (/40 nm)

First experiment (Partec CyFlow ML)


StableYFP(LOV)-versus negative control



Green fluorescence (arbitrary units)

Filename	# of Events	Median	Geometric Mean
Negative control	856154	46	0
YFPstable(LOV)	99147865536	111	0


kntGFP-versus negative control



Green fluorescence (arbitrary units)

Filename	Median	Geometric Mean
Negative control	46	0
coli-LGFPKNTbisgain FL6-350_021.FCS	41	0

Positive YFP control condition-versus negative control




Green fluorescence (arbitrary units)

Filename	Median	Geometric Mean
Negative control	46	0
YFP positive control	611	4.94

Second experiment (performed with PARTEC CUBE1 cytometer)


StableYFP(LOV)-versus negative control



Green fluorescence (arbitrary units)

Filename	Median	Geometric Mean
Negative control	4.49	3.55
StableYFP(LOV)	13.49	9.51


kntGFP-versus negative control



Green fluorescence (arbitrary units)

Filename	Median	Geometric Mean
Negative control	4.49	3.55
Knt-GFP	6.63	4.94

Positive YFP control condition-versus negative control



Green fluorescence (arbitrary units)

Filename	Median	Geometric Mean
Negative control	4.49	3.55
YFP positive control	101	8.70

Culture of YFP and knt-GFP transformants on middles with different antibiotics concentration gradients

THURSDAY, 9/3/2020


Composition of media is written on the photos. We tried the following aminoglycoside antibiotics: spectinomycin, streptomycin and kanamycin.

1 : positive control DH5alpha (plasmid with SpecR)


2 : DH5alpha (pBAD24-MoClo-knt-GFP)

3 : negative control DH5alpha (pBAD24-MoClo)


118774506_3482563185140124_145799137588464
8640_n.jpg




No difference between 2 and 3 for spectinomycin resistance



118780229_3562645253754703_4969674500013
421743_n.jpg




No difference between 2 and 3 for streptomycin resistance



118771223_254389875635796_7414920641372441
991_n.jpg



No difference between 2 and 3 for kanamycin resistance



It seems that knt-GFP did not confer resistance against aminoglycoside antibiotics to our bacteria.

Comparison sequences knt-GFP designed and sequenced

WEDNESDAY, 9/30/2020

We send DNA inserted in our transformants bacterias Stable-YFP to IDT for sequencing with the amores :

5'- TGGCTATGCCATAGCATTATCC (forward)

5'- GGT CTG ATT TAA TCT GTA TCA GGC TG (riverse)

When we received the results of the sequencing we compared it to the sequence we designed with CAMEOS and Benchling with BLAST.

The result is that the query cover is 100% with an E-value of 0,00 and an identity percentage of 99,87% for forward and reverse primer.

For the sequencig with the forward primer it seems that there is a single mutation of type substitution at the nucleotide 63 of our gene where a C has been replaced by a A. However, the sequencing shows that at this position (753) sequencing's quality isn't optimal, so there is a possiblity there is in fact no mutation.

For the sequencig with the reverse primer it seems that there is a single mutation of type substitution at the nucleotide 708 of our gene where a G has been replaced by a T.


In conclusion because two difrent mutation have been detected with the two primers, we can suppose there is no mutation of our gene knt-GFP.

Screen captures :

Forward sequences :

Blast :


Description	Max Score	Total Score	Coverage	E-value	Identity	Accession
GST094_1039847_1039847	1403	1403	100%	0.0	99.87%	Query_14039



Sequencing :



Reverse sequences :



Sequencing :

