
Plasmid Extraction


Nanjing_NFLS

质粒 (Plasmids)

- Plasmids are loops of DNA **outside the genome**.
- Plasmid size ranged from 1 kb-200 kb, much smaller than the E. coli genome (4.6 Mb).
- Plasmids are communication tools for microbes.


monoclonal

- Replication start site: determines the number of copies of the plasmid

- Polyclonal sites.
- Easy insertion of DNA fragments

Plasmid Extraction: Base Cleavage, Over-Column Purification

Practice using a pipette

- 1 step for sampling, 2 stop for adding samples
- Tune a gun.
- Put a gun on the head.
- Liquid extraction.
- Hold the tip of the gun against the wall of the barrel to sample, or stick the tip below the liquid surface to sample.
- Immediately after the addition of the sample, leave the nozzle from the liquid surface, do no suction back.

Practice adding liquid to the walls of the tube: 300 µL, 30 µL, 3 µL.

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Plasmid Extraction

9. The liquid in the collection tube was discarded, and after loading back into the adsorption column, the column was idled at 13,000 rpm for 2 minutes to remove ethanol impurities.
10. 取一只干净的1.5ml离心管，标记pV和你的序号，将步骤8后的吸附柱装进离心管，打开盖子，室温放置5分钟，挥发酒精。
11. 向吸附柱内添加50ul Elution buffer，室温静置5分钟，12000rpm离心2分钟。

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