

M9 Minimal Medium Preparation

1. 5X M9 salt solution preparation:

To 400mL of ddH₂O add:

- a. 64g ($\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$)
- b. 15g (KH_2PO_4)
- c. 2.5g (NaCl)
- d. 5.0g (NH_4Cl)

Make to 1 L with ddH₂O

Sterilize by autoclaving

2. Preparation of 1 M MgSO₄ solution:

To 100 mL ddH₂O add:

- a. 24.65 g ($\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$)

3. Preparation of 40% glucose (w/v) solution:

To 100 mL distilled/deionized water add:

- a. 40 g (glucose)

Caution: Add glucose to stirring water in beaker; Do not attempt to add water to GLUCOSE!

4. Preparation of 1M CaCl₂ solution:

To 100 mL ddH₂O add:

- a. 147.014g ($\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$)

Make to 1 L with dH₂O

5. Method for preparation of (1L of media) minimal medium:

I. 200 mL 5X M9 salts solution

II. 800 mL ddH₂O

III. 2 mL of 1M MgSO₄ solution

IV. 0.1 mL of 1M CaCl₂ solution

(Our Lab put I-IV solutions together in order to make mixture buffer)

V. 10 mL of 40% glucose