

## Protocol for PCR & Agarose Gel Electrophoresis - System B

### *·Material*

NovoProtein 2×Taq Master Mix

Template P16-1 Plasmid

Template P16-2 Plasmid

Template P18 Plasmid

Primer

Double Distillation Water (ddH<sub>2</sub>O)

DNA Marker

### *·Step*

Group 1\2\3\4

- ① Adding 10μl of Takara Prime STAR HS, 1μl of tablet, 4μl of primer, 5μl of ddH<sub>2</sub>O into a PCR tube. Different group use specific tube with distinctive sign. Those four group are arranged as table follow:

Group	Primer	Template
1	Seq-P16-F1 Seq-P16-R1	P16-1
2	Seq-P16-F1 Seq-P16-R1	P16-2
3	Seq-P16-F2 Seq-P16-R2	P16-1 or P16-2
4	Seq-P16-F2 Seq-P18-R2	P18

- ② Place those PCR tubes into Peltier thermal cyclor.
- ③ Set the protocol as follow: begin at 98°C for 5 mins, then keep 98°C for 10 secs for denaturation, decrease to 54°C for 15 secs, 72°C for 1 min and repeat that cycle 30 times, finally maintain 4°C infinity.
- ④ Adding 5μl different samples for each into different wells(Adding DNA marker in the first well), then run the gel at 120V for about 30 minutes. Check the result under UV rays.

Group 5

- ① Adding 10μl of Takara Prime STAR HS, 1μl of tablet, 4μl of primer, 5μl of ddH<sub>2</sub>O into a PCR tube. Different group use specific tube with distinctive sign. Those four group are arranged as table follow:

Group	Primer	Template
5	Seq-P16-F1 Seq-P18-R1	P18

- ② Place the PCR tube into Peltier thermal cycler.
- ③ Set the protocol as follow: begin at 98°C for 5 mins, then keep 98°C for 10 secs for denaturation, decrease to 52°C for 15 secs, 72°C for 1 min and repeat that cycle 30 times, finally maintain 4°C infinity.
- ④ Adding 5 $\mu$ l different samples for each into different wells(Adding DNA marker in the first well), then run the gel at 120V for about 30 minutes. Check the result under UV rays.