

# Site Directed Mutagenesis

---

## Introduction

[LucidChart Overviews](#) | [DNA](#) | [Combined](#)

[Procedures List](#)

[Q5® Site-Directed Mutagenesis Kit Protocol](#)

[Product Manual](#)

## Materials

- › [Q5 Site-Directed Mutagenesis Kit](#)
  - › [NEB® 5-alpha Competent \*E. coli\* \(High Efficiency\)](#)
  - › [Q5®-Hot Start High-Fidelity 2X Master Mix](#)
  - › [SOC Outgrowth Medium](#)
  - › [Control SDM Plasmid](#)
  - › [Control SDM Primer Mix](#)
  - › [KLD Enzyme Mix](#)
  - › [KLD Reaction Buffer](#)
  - › [pUC19 Vector](#)
- › [Forward and Reverse Primers](#)
- › [Template DNA](#)
- › [Nuclease-free water](#)
- › [Ice](#)
- › [Selection Plates](#)

## Procedure

### Step I: Exponential Amplification (PCR)

1. Assemble the following reagents in a thin-walled PCR tube.

	A	B	C
1		<b>25 µl RXN</b>	<b>Final Concentration</b>
2	<b>Q5 Hot Start High-Fidelity 2X Master Mix</b>	12.5 µl	1X
3	<b>10 µM Forward Primer</b>	1.25 µl	0.5 µM
4	<b>10 µM Reverse Primer</b>	1.25 µl	0.5 µM
5	<b>Template DNA (1–25 ng/µl)</b>	1 µl	1-25 ng
6	<b>Nuclease-free water</b>	9.0 µl	

2. Mix reagents completely, then transfer to a thermocycler.

3. Perform the following cycling conditions:

Thermocycling Conditions for a Routine PCR:

*\* For a Q5-optimized annealing temperature of mutagenic primers, please use [NEBaseChanger™](#), the online NEB primer design software. For pre-designed, back-to-back primer sets, a  $T_a = T_m + 3$  rule can be applied, but optimization may be necessary.*

	A	B	C
1	<b>Step</b>	<b>Temperature °C</b>	<b>Time</b>
2	<b>Initial Denaturation</b>	98	30 seconds
3		98	10 seconds
4	<b>25 Cycles</b>	50-72	10-30 seconds
5		72	20-30 seconds/kb
6	<b>Final Extension</b>	72	2 minutes
7	<b>Hold</b>	4-10	

## Step II: Kinase, Ligase & DpnI (KLD) Treatment

4. Assemble the following reagents:

	A	B	C
1		<b>Volume</b>	<b>Final Concentration</b>
2	<b>PCR Product</b>	1 µl	
3	<b>2X KLD Reaction Buffer</b>	5 µl	1X
4	<b>10X KLD Enzyme Mix</b>	1 µl	1X
5	<b>Nuclease-free Water</b>	3 µl	

5. Mix well by pipetting up and down and incubate at room temperature for 5 minutes.

## Step III: Transformation

6. Thaw a tube of NEB 5-alpha Competent *E. coli* cells on ice.

7. Add 5  $\mu$ l of the KLD mix from Step II to the tube of thawed cells. Carefully flick the tube 4-5 times to mix. Do not vortex.

8. Place the mixture on ice for 30 minutes.

9. Heat shock at 42°C for 30 seconds.

10. Place on ice for 5 minutes.

11. Pipette 950  $\mu$ l of room temperature SOC into the mixture.

12. Incubate at 37°C for 60 minutes with shaking (250 rpm).

13. Mix the cells thoroughly by flicking the tube and inverting, then spread 50-100  $\mu$ l onto a selection plate and incubate overnight at 37°C. It may be necessary (particularly for simple substitution and deletion experiments) to make a 10- to 40-fold dilution of the transformation mix in SOC prior to plating, to avoid a lawn of colonies