

Reagent

apt-1 (100 μ M)

apt-2 (100 μ M)

activator (100 μ M)

Cas12a (6nM in buffer A)

crRNA (3nM in RNase free buffer)

ssDNA-FQ (50nM, in RNase free buffer)

Kanamycin

Buffer A pH7.4 40 mM HEPES, 100 mM NaCl, 20 mM MgCl₂

Preparation of solutions

Solution 1-probe solution 80 μ L

Reagent	Volume/ μ L
apt-1 4 \times	20
apt-2 4 \times	20
Activator 8 \times	10
Buffer A	30

The concentrations of the three oligonucleotides in the solution are 25nM,25nM and 12.5nM, respectively.

1.Incubate at 90 $^{\circ}$ C for 10min

2.slow cooling on the benchtop

Solution 2-reporter solution 120 μ L

Reagent	Volume/ μ L
Cas12a (6nM in buffer A)	40
crRNA (3nM in RNase free buffer)	40
ssDNA-FQ	40

1.Mix 40 μ L Cas12a (6nM in buffer A) and 40 μ L crRNA (3nM in RNase free buffer) and wait for 30 minutes at 25 $^{\circ}$ C

2.Add 40 μ L ssDNA-FQ (50nM) to the mixture

Solution 3-kanamycin

120 μ M kanamycin solution

Detection

Solution/ μ L Group	1 Control	2 (60 μ M kana)	3 (120 μ M kana)
S1	20		
S2	20		
S3	0	20	40
Buffer	40	20	0
Total volume	100	100	100

Each tube:

- 1.Mix Buffer with Solution 3
- 2.Mix 20 μ L solution 1 (probe) and the mixture got from step 1. Incubate for 5 min.
- 3.Add 60 μ L solution 2 (reporter) and incubate at 25 $^{\circ}$ C for 35 min.
- 4.Incubate at 60 $^{\circ}$ C for 5 min (stop enzymatic reaction)
- 5.Detect the fluorescence intensity (the excitation wavelength is 485nm and the emission was monitored at 520nm)