

## Electrocompetent BL21 Cell Transformation

### Introduction

This protocol yields transformed electrocompetent cells.

### Reagents

- 🔗 BL21 electrocompetent cells
- 🔗 Assembly product DNA
- 🔗 SOB liquid media
- 🔗 LB/Kanamycin plates

### Equipment

- 🔗 Electroporator
- 🔗 Electroporation cuvette
- 🔗 Eppendorf tubes
- 🔗 Ice in a bucket
- 🔗 Pipette and tips
- 🔗 Incubation cabinet

### Procedure

1. Thaw BL21 electrocompetent cells on ice. Place the electroporation cuvette on ice to chill.
2. Add 1  $\mu\text{L}$  of the desired assembly production to 50  $\mu\text{L}$  of electrocompetent cells in a sterile Eppendorf tube. Mix gently by pipetting up and down.
3. Transfer the contents of the Eppendorf tube to the chilled electroporation cuvette.
  - a. *Cut a larger opening on a pipette tip so as to not damage the cells.*
4. Insert the electroporation cuvette into the electroporator and carry out the procedure at 1.8 kV.
5. Remove transformed cells, place in a second Eppendorf tube, and add 950  $\mu\text{L}$  of SOB media.
6. Incubate the tube at 37°C for 60 minutes. Shake vigorously. Warm selection plates at 37°C.
7. Spread 100  $\mu\text{L}$  of the transformed cells onto the plates and incubate overnight at 37°C.