

Heat shock competent cells

Aim of the experiment

The purpose of this experiment is to transfer the plasmid into E.coli competent cells by heat shock.

Materials

- E.coli Rosetta 2 (DE3) pLysS、E.coli BL21、E.coli DH5alpha Competent bacteria;
- pC013-Twinstrep-SUMO-huLwCas13a、
pC0061 PsmCas13 (B05) His6-TwinStrep-SUMO-Bsal、
p2CT-His-MBP-Lba_Cas13a_WT;

- LB medium;

- Chloramphenicol (Cm, 35ug/ml) :

Dissolve 0.35g of Amp by adding it to 8 mL of absolute ethanol and then vortexing. Add absolute ethanol to bring the volume to 10 mL. Store at -20 °C;

- Ampicillin (Amp, 100ug/ml) :

Dissolve 0.1g of Amp by adding it to 8 mL of deionized water and then vortexing. Add deionized water to bring the volume to 10 mL and filter-sterilize with a 0.22-μm syringe filter. Store at -20 °C.

Procedure

1. Thaw one vial of E.coli competent cells on ice for 30 min, and then add 2 μL of plasmid. Incubate on ice for 5 min.
 2. Heat-shock the cells by placing the vial into a 42 °C pre-heated water bath for 90 s, and then cold-shock the cells on ice for 2 min.
-

3. Add 200 μL of LB medium to the cells and plate 100 μL of cell suspension on a pre-warmed LB plate containing 100 $\mu\text{g/mL}$ ampicillin and 35 $\mu\text{g/mL}$ Chloramphenicol. Incubate the plate overnight in a 37 $^{\circ}\text{C}$ incubator.