

In vitro Transcription

Aim of the experiment

We used in vitro transcription (IVT) to synthesize RNA from DNA templates. The first protocol describes IVT from complete double stranded templates. For generation of crRNAs we developed a protocol using incomplete templates to facilitate variation of target sequences. crRNA templates were completed using DNA polymerase and transcribed in a one-batch reaction.

Materials

Single tube reaction system (20 μ L) :

	Positive control group	Experimental groupA	Experimental groupB	Experimental groupC
T7 Control template	2 μ L			
DNA Template (100 ng/ μ L)		1 μ L	2 μ L	3 μ L
10X Transcription Buffer	2 μ L	2 μ L	2 μ L	2 μ L
rNTP (100 mM ench)	5 μ L	5 μ L	5 μ L	5 μ L
T7 RNA Polymerase	2 μ L	2 μ L	2 μ L	2 μ L
RNase Free H ₂ O	9 μ L	10 μ L	9 μ L	8 μ L
Total	20 μ L	20 μ L	20 μ L	20 μ L

Procedure

1. According to the experimental materials prepare the transcription mix.
2. The incubation:

(1) Gently absorb and mix the reaction system with a pipette (do not brush strong

oscillation). Instantaneous centrifugation ensures that all the solution is collected at the bottom of the tube.

(2) The reaction tube was placed in the PCR instrument (or other thermostatic equipment) and incubated with 2-4h at 37°C (partially inactivated at 70 °C).

Take out the reaction products and place them on ice to measure the concentration.

3. Take out the reaction products and place them on ice to measure the concentration.
 4. Next test or -80C store as soon as possible.
-