

## Subcloning: digestion & ligation reaction protocol

### Introduction

Subcloning is the procedure of moving a DNA fragment in a certain plasmid to another. This process requires digestion of a “parent” plasmid to obtain the DNA fragment and/or insertion plasmid with complementary sticky ends to the sequence it will be ligated with in the final construct. For that, a ligation reaction is also required. This protocol describes a single digestion and ligation reaction of a DNA fragment and plasmid with PstI and EcoRI restriction sites.

This protocol is adapted from iGEM's 3A Assembly protocol:

### Materials

- DNA plasmids
- ddH<sub>2</sub>O
- Digestion cutSMART 10x buffer
- PstI restriction endonuclease
- EcoRI restriction endonuclease
- T4 DNA ligase
- Ligase 10x buffer

### Procedure

The following procedure is for one digestion and ligation reaction, respectively.

#### Digestion

1. Mix 0.5  $\mu$ l Digestion cutSMART 10x buffer with 4.3  $\mu$ l ddH<sub>2</sub>O and add 0.1  $\mu$ l of each enzyme. This is your master mix.
2. Transfer 4  $\mu$ l of the master mix into a PCR tube.
3. Add 100 ng of DNA plasmid to the same PCR tube by adding 4  $\mu$ l of DNA with concentration 25 ng/  $\mu$ l.
4. Tap on the tube to mix or centrifuge for a few seconds to spin down the liquid.
5. Put the PCR tube into a thermocycler with the following settings for incubation and heat-inactivation of the enzymes.

Step	Temperature (°C)	Time (min)
1	37	30
2	80	20

3	4	Infinite
---	---	----------

The sample may be stored in the -20 freezer.

### Ligation (10 µl reaction mixture, 2 DNA parts)

1. Mix 1 µl 10X T4 DNA ligase buffer with 4.5 µl ddH<sub>2</sub>O in a PCR tube.
2. Add 2 µl of each part you want to ligate from the above digestion.
3. Add 0.5 µl of T4 DNA ligase.
4. Put the PCR tube into a thermocycler with the following settings for incubation and heat-inactivation of the enzyme.

Step	Temperature (°C)	Time (min)
1	16	30
2	80	20
3	4	Infinite

The sample may be stored in the -20°C freezer.

## References

Parts.igem.org. 2020. *3A Assembly*. [online] Available at: <[http://parts.igem.org/Help:Protocols/3A\\_Assembly](http://parts.igem.org/Help:Protocols/3A_Assembly)> [Accessed 5 December 2020].