

Fluorescence intensity quantification protocol

To measure fluorescence intensity in whole bacterial cells.

Adapted from iGEM Wageningen 2019

(<https://www.protocols.io/view/gfp-rfp-plate-reader-assay-784hryw>)

Materials

- Bacteria
 - From plates or glycerol stocks
- Luria-Bertani (LB) medium
- Antibiotics (if required)
- Phosphate buffered saline (PBS) 1X, pH 7.4
- Eppendorf tubes
- 96-well plate, clear bottom
- CLARIOstar Plus Microplate Reader

Procedure

1. Calibrate the microplate reader to obtain a standard curve:
 - a. For OD600 measurement, calibrate with silica beads.
 - b. For GFP fluorescence intensity measurement, calibrate with fluorescein.
 - c. For RFP fluorescence intensity measurement, calibrate with Texas Red.
2. Set the following protocols for CLARIOstar Plus Microplate Reader:
 - a. To measure OD600:

Measurement type: absorbance, endpoint

Excitation wavelength: 600 nm

Top optic

Temperature: 37°C

- b. To measure GFP fluorescence intensity:

Measurement type: fluorescence (FI), endpoint

Excitation wavelength: 483 nm

Emission wavelength: 530 nm

Gain: 680

Dichroic filter: 502.5

Focal height: 10 mm

Top optic

Temperature: 37°C

c. To measure RFP fluorescence intensity:

Excitation wavelength: 560 nm

Emission wavelength: 615 nm

Gain: 910

Dichroic filter: 582.5

Focal height: 10 mm

Top optic

Temperature: 37°C

3. Prepare bacteria cultures (from plates or glycerol stocks) in 10 mL LB media supplied with antibiotics (if required).
4. Incubate cultures overnight at 37°C, 150 rpm.
5. Transfer 1 mL of overnight culture into an Eppendorf tube.
6. Centrifuge at 13000 rpm for 5 minutes to pellet cells.
7. Resuspend the cell pellet in 1 mL PBS.
8. Repeat centrifuging and washing steps.
9. In a clear bottom 96-well plate, load 160 μ L PBS and 40 μ L bacteria suspension.
 - a. As negative control, load 160 μ L PBS and 40 μ L non-fluorescent bacteria suspension.
 - b. As blank, load 200 μ L PBS.
10. Measure OD600 and fluorescence intensity of the samples using the protocols above.
11. Adjust fluorescence intensity measurements for OD600 values.