

Mini-prep kit – plasmid extraction from bacteria cells

You can extract high/low copy number plasmid out of bacteria cells with this kit.

If you need to keep the starter sterile use sterile conditions.

High Copy Number Protocol

Harvesting -

1. Place **2 mL** of bacterial culture into a 2 mL microcentrifuge tube.
2. Centrifuge for 1 min at full speed (table microcentrifuge) and discard supernatant.
3. Repeat steps 1-2 (total of 4 mL).

Resuspension -

4. Add 200 μ L of PD1 Buffer (RNase A added) and pipette (avoid bubbles).

Lysis -

5. Add 200 μ L of PD2 Buffer and mix gently by inverting the tube 10 times.
Do not vortex, avoid shearing genomic DNA.
6. Allow mixture to stand for 2 minutes at room temperature until lysate clears.

Neutralization -

7. Add 300 μ L of PD3 Buffer and mix immediately by inverting the tube 10 times.
Do not vortex.
8. Centrifuge for 2 min at 13,000 rpm (table microcentrifuge).

DNA Binding -

9. Place a PD Column in a 2 mL Collection Tube (provided in the kit), mark with plasmid name.
10. Pour the clear lysate (supernatant) from Step 7 to the PD Column.
11. Centrifuge at 13,000 rpm for 30 seconds.
12. Discard the flow-through and return the PD Column to the 2 mL Collection Tube.

Wash -

13. Add 400 μ L of W1 Buffer in the PD Column.
14. Centrifuge at 13,000 rpm for 30 seconds.
15. Discard the flow-through and return the PD Column to the 2 mL Collection Tube.
16. Add 600 μ L of Wash Buffer (ethanol added) to PD Column.
17. Centrifuge at 13,000 rpm for 30 Seconds.

18. Discard the flow-through and return the PD Column to the 2 mL Collection Tube.
19. Centrifuge again for 3 min at full speed to dry the column matrix. If the column matrix is not dry, repeat centrifuge.

DNA Elution -

20. Transfer the dried PD Column to a clean 1.5ml microcentrifuge tube (not provided), mark it properly.
21. Add 50 μ L DDW directly onto the centre of the membrane. Avoid residual buffer adhering to the wall of the column.
22. Allow to stand for 2 min until the liquid is absorbed.
23. Centrifuge for 2 min at 13,000 rpm to elute plasmid DNA.
24. After plasmid purification is complete measure the concentration in the Nano-Drop.
25. Plasmid can be stored at -20°C.