

Author: Allyssa Hinkle  
 Entry 4/4: Assay Protocol  
 In Project: iGEM TUK 2020 Labjournal  
 No tags associated

created: 22.10.2020 06:32  
 updated: 27.10.2020 18:29

## Lab Journal

### Table of Contents

## 5/25/2020-Made Buffers

Stock Solutions				Notes
	<b>Buffers</b>			
	<b>Basic Solution</b>			
	Company	Carl Roth GmbH+Co.	$\text{Na}_2\text{HPO}_4$	Di-natriumhydrogen phosphat-dihydrat
		Formular Weight: 177.99g/mol		Calculation
	Calculated	44.4975g per 250ml		$177.99\text{g/mol} \cdot 0.25 = 44.4975\text{g}$ per 250ml
	Actual	44.5g	Added powder to 250ml bottle	
		Making buffer with ultra filtrated $\text{H}_2\text{O}$		
		Fill 250ml bottle to 200ml with UF $\text{H}_2\text{O}$		
		Moved this mixture to 1 L container		

		Added H <sub>2</sub> O to 1L creating 0.25 molar concentration		
	Acidic Solution			
	Company	Sigma	Citric Acid	trisodium salt dihydrate
		Formular Weight: 294.1g/mol		Calculation
	Calculated	For 250ml solution: 73.525g		294.1g/mol*0.25=73.525g in 250ml
	Calculated	For 100ml solution: 29.41g		73.525g(100ml/250ml)=29.41g
	Actual	29.42g		
		Filled 250ml bottle with 29.42g of powder and willed to 50ml with UF H <sub>2</sub> O		
		Stirred/shook till clear then filled to 100ml		
		Created 1 molar solution		
		Citric acid monohydrate		
		Formular weight: 210,14g/mol		210,14g/mol*0.25=52,535g in 250ml

After recalculating to avoid second dilution, this is corrected stock solution amounts:

Na <sub>2</sub> HPO <sub>4</sub>	NEED		
177,99 g/mol	0,2 M	200 ml	
CALCULATION: $177,99 \text{ g/mol} * 0,2 \text{ M} * 0,2 \text{ l} = 7,119 \text{ g}$			
Citric acid	NEED		
210,14 g/mol	0,1 M	200 ml	
CALCULATION: $210,14 \text{ g/mol} * 0,1 \text{ M} * 0,2 \text{ l} = 4,203 \text{ g}$			
pH 7: 164,7 ml Na <sub>2</sub> HPO <sub>4</sub> + 35,3 ml Citric acid			
	For 200 ml of Buffer		
pH 4: 77,1 ml Na <sub>2</sub> HPO <sub>4</sub> + 122,9 ml Citric acid			

## 5/26/2020-Buffer Dilutions

Calculations					
Enzyme	0.5u/mg in bottle		1u=1μmol/min	1μmol=0.1μ ABTS	10m incubation
	10u/ml	=	20mg/ml		
	0.5u/mg		c=1mM ABTS	mM=mmol/L	

				$c=n/v$	$v=1\text{ml}=0.001\text{L}$	
				1 $\mu\text{mol}$ ABTS (n) in 1 ml (v)		
				1 $\mu\text{mol}/1\text{ml}$		
				1 $\mu/10\text{m}=0.$ 1 $\mu$ =need 1/5mg		
<b>Citric Acid Buffer</b>						
	0.1 M want	(100ml)	10ml buffer	or for 200ml	20ml buffer	
	1M have		90ml H <sub>2</sub> O		180ml H <sub>2</sub> O	
<b>Basic Buffer</b>	0.2 M Want	(100ml)	80ml buffer	or for 200ml	160ml buffer	
	0.25 M Have		20ml H <sub>2</sub> O		40ml H <sub>2</sub> O	
Prepared	Basic Buffer	0.2M				
	and					
	Citric Acid Buffer	0.1M				
	with above buffers we made 3 different pH buffers as follows					
	pH want	0.2M Na <sub>2</sub> HPO <sub>4</sub>	Citric Acid	Actual pH	Final pH*	
	4	7.71ml	12.29ml	8,51	3,95	
	5	10.3ml	9.7ml	9	5,04	
	7	16.47ml	3.53ml	9,28	7,04	
	*This is the final pH after adjustment using conc. Hydrochloric acid at approximately 32% and 10M NaOH					

ABTS		$C_{18}H_{16}N_4O_6S_4^-$ ( $NH_4$ ) <sub>2</sub> -Mr548.7				
	c=n/v	54,87	=	0.05487g/ml		
		1000				
		548.7g=1mol				
		54.87 g=100mmol				
Actual	enzyme	0.278mg in 1.39 ml 7.0 pH solution			v=m/c	
	Claculate d ABTS:	actual weighted /c				
		0.0584g/0. 0547g	=	1.068ml UF H <sub>2</sub> O		


### 5/27/20-ABTS abs spectra with *T. veriscolor* at different pH's

ABTS:	1mU	(1ml)	=	0.01 ml	=	1 0 uL
	100mU					
Enzy me:	0.2mg/ml	=	1ml	=	0.01 ml	= 10uL
	20mg/ml					
Buffer:	total volume	=	980uL			
	Filled 9 eppendorf tubes with:					
	pH 4 buffer 980uL	10uL UF H <sub>2</sub> O	10uL ABTS			
	pH 4 buffer 980uL	10uL UF H <sub>2</sub> O	10uL Enzyme			

	pH 4 buffer 980uL	10uL ABTS	10uL Enzyme				
	pH 5 buffer 980uL	10uL UF H <sub>2</sub> O	10uL ABTS				
	pH 5 buffer 980uL	10uL UF H <sub>2</sub> O	10uL Enzyme				
	pH 5 buffer 980uL	10uL ABTS	10uL Enzyme				
	pH 7 buffer 980uL	10uL UF H <sub>2</sub> O	10uL ABTS				
	pH 7 buffer 980uL	10uL UF H <sub>2</sub> O	10uL Enzyme				
	pH 7 buffer 980uL	10uL ABTS	10uL Enzyme				
	UV-Visible System Spectrometer-Chem Station						
	Using pH 7 buffer to blank for enzyme						
	See folders /documents on Allyssa's flashdrive for data						
Time point	1	pH 4 buffer with ABTS					
	2	pH 4 buffer with ABTS and enzyme					
	3	pH 4 buffer with ABTS and enzyme		2m			
	4	pH 4 buffer with ABTS and enzyme		4m			

	5	pH 4 buffer with ABTS and enzyme	8m		
	6	pH 4 buffer with ABTS and enzyme	10m		
	7	pH 4 buffer with ABTS and enzyme	12m		
	8	pH 4 buffer with ABTS and enzyme	15m		
Time point	1	pH 5 buffer with ABTS			
	2	pH 5 buffer with ABTS and enzyme			
	3	pH 5 buffer with ABTS and enzyme	2m		
	4	pH 5 buffer with ABTS and enzyme	4m		
	5	pH 5 buffer with ABTS and enzyme	6m		
	6	pH 5 buffer with ABTS and enzyme	8m		
	7	pH 5 buffer with ABTS and enzyme	10m		
	8	pH 5 buffer with ABTS and enzyme	12m		
Time point	1	pH 7 buffer with ABTS			
	2	pH 7 buffer with ABTS and enzyme			
	3	pH 7 buffer with ABTS and enzyme	2m		

4	pH 7 buffer with ABTS and enzyme	4m		
5	pH 7 buffer with ABTS and enzyme	6m		
6	pH 7 buffer with ABTS and enzyme	8m		
7	pH 7 buffer with ABTS and enzyme	10m		
8	pH 7 buffer with ABTS and enzyme	12m		
9	pH 7 buffer with ABTS and enzyme	13m		
	*For future experiments maybe lower enzyme concentration			

 [5-27 pH4 data.xlsx](#)

 [5-27 pH5.xlsx](#)

 [5-27 pH7.xlsx](#)

### 6/2/2020-*T. versicolor* Positive Control ABTS Assay

pH	Concentrations							ABTS	pH4	200µl	20µl	180µl
	Enzyme ug per ml											
4	0	2	10	20	50	100	200	ABTS	ABTS	200µl	20µl	180µl
4	0	2	10	20	50	100	200	ABTS	Enzyme	200µl	20µl	180µl



control	0	2	10	20	50	100	200	ABTS pH4 control		pH5			
5	0	2	10	20	50	100	200	ABTS		ABTS	200µl	20µl ABTS	180µl pH5 buffer
5	0	2	10	20	50	100	200	ABTS		Enzyme	200µl	20µl Enzyme	180µl buffer
control	0	2	10	20	50	100	200	ABTS pH5 control		pH7			
7	0	2	10	20	50	100	200	ABTS		ABTS	200µl	20µl ABTS	180µl pH7 buffer
7	0	2	10	20	50	100	200	ABTS		Enzyme	200µl	20µl Enzyme	180µl buffer
control= enzyme without ABTS								14 controls without ABTS *pH4 and 5					
								8 controls without Enzyme in all pH					
0.5µl Conc value µg/ml enzyme+0. 5µlABTS+199µl buffer								Can't pipette 0.5microl, need to adjust dilutions with buffer and increase to 5microl each					
5µl Conc value µg/ml enzyme+5µlABTS+190µl buffer													

## Concentration calculations

Dilution amounts		
Concentration enzyme	$\mu\text{l}$ enzyme/previous dilution	Notes
2	100 total volume=20microl 20 enzyme+80microl buffer	
10	100 total volume=50microl 20 enzyme+50microl buffer	
20	200 total volume=80microl 50 enzyme+120microl buffer	*final volume 130microl after r
50	200 total volume=100microl 100 enzyme+100microl buffer	*final volume 120microl after r
100	200 total volume=100microl 200 enzyme+100microl buffer	*final volume 100microl after r
200	200 total volume=20microl stock enzyme+180microl buffer	*final volume 100microl after r

6/3/2020-*T. versicolor* Positive Control ABTS Assay*T. versicolor* prep

Enzyme dilution prep			
pH 4		pH 5	
2	100 total volume=20microl 20 enzyme(pH4)+80microl buffer(pH4)	2	100 total volume=20microl 20 enzyme(pH5)+80microl buffer(pH5)
10	100 total volume=50microl 20 enzyme(pH4)+50microl buffer(pH4)	10	100 total volume=50microl 20 enzyme(pH5)+50microl buffer(pH5)
20	200 total volume=80microl 50 enzyme(pH4)+120microl buffer(pH4)	20	200 total volume=80microl 50 enzyme(pH5)+120microl buffer(pH5)
50	200 total volume=100microl 100 enzyme(pH4)+100microl buffer (pH4)	50	200 total volume=100microl 100 enzyme(pH5)+100microl buffer (pH5)
100	200 total volume=100microl 200 enzyme(pH4)+100microl buffer (pH4)	100	200 total volume=100microl 200 enzyme(pH5)+100microl buffer (pH5)
200	200 total volume=20microl stock enzyme (pH7)+180microl buffer (pH4)	200	200 total volume=20microl stock enzyme (pH7)+180microl buffer (pH5)

ABTS Dilution prep			
ABTS	200µl	20µl ABTS	180µl pH4 buffer
ABTS	200µl	20µl ABTS	180µl pH5 buffer
ABTS	200µl	20µl ABTS	180µl pH7 buffer

96 well plate set up:

Well	3	4	5	6	7	8	9	
Enzyme Concentration	0	2	10	20	50	100	200	
pH	Concentrations Enzyme µg per ml							
4	5µl ABTS + 195µl buffer	5µl 2enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 10enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 20enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 50enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 100enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 200enzyme µg/ml + 5µl ABTS + 190µl buffer	A B T S
4	5µl ABTS + 195µl buffer	5µl 2enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 10enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 20enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 50enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 100enzyme µg/ml + 5µl ABTS + 190µl buffer	5µl 200enzyme µg/ml + 5µl ABTS + 190µl buffer	A B T S

								5µl 200enzy µg/ml +195µl buffer	A B T S p H 4c on tr ol
cont rol	195µl buffer+5u l buffer	5µl 2enzy µg /ml +195µl buffer	5µl 10enzy µg/ml +195µl buffer	5µl 20enzy µg/ml +195µl buffer	5µl 50enzy µg/ml +195µl buffer	5µl 100enzy µg/ml +195µl buffer			
5	5µlABTS +195µl buffer	5µl 2enzy µg /ml +5µlABTS+ 190µl buffer	5µl 10enzy µg/ml +5µlABTS+ 190µl buffer	5µl 20enzy µg/ml +5µlABTS+ 190µl buffer	5µl 50enzy µg/ml +5µlABTS+ 190µl buffer	5µl 100enzy µg/ml +5µlABTS+ 190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	A B T S	
5	5µlABTS +195µl buffer	5µl 2enzy µg /ml +5µlABTS+ 190µl buffer	5µl 10enzy µg/ml +5µlABTS+ 190µl buffer	5µl 20enzy µg/ml +5µlABTS+ 190µl buffer	5µl 50enzy µg/ml +5µlABTS+ 190µl buffer	5µl 100enzy µg/ml +5µlABTS+ 190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	A B T S	
cont rol	195ul buffer+5u l buffer	5µl 2enzy µg /ml +195µl buffer	5µl 10enzy µg/ml +195µl buffer	5µl 20enzy µg/ml +195µl buffer	5µl 50enzy µg/ml +195µl buffer	5µl 100enzy µg/ml +195µl buffer	5µl 200enzy µg/ml +195µl buffer	A B T S p H 5c on tr ol	
7	5µlABTS +195µl buffer	5µl 2enzy µg /ml +5µlABTS+ 190µl buffer	5µl 10enzy µg/ml +5µlABTS+ 190µl buffer	5µl 20enzy µg/ml +5µlABTS+ 190µl buffer	5µl 50enzy µg/ml +5µlABTS+ 190µl buffer	5µl 100enzy µg/ml +5µlABTS+ 190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	A B T S	

7	5µl ABTS +195µl buffer	5µl 2enzy µg /ml +5µl ABTS+ 190µl buffer	5µl 10enzy µg/ml +5µl ABTS+ 190µl buffer	5µl 20enzy µg/ml +5µl ABTS+ 190µl buffer	5µl 50enzy µg/ml +5µl ABTS+ 190µl buffer	5µl 100enzy µg/ml +5µl ABTS+ 190µl buffer	5µl 200enzy µg/ml +5µl ABT S+190µl buffer	A B TS
5µl Con c valu e µg /ml enzy me+ 5µl AB TS+ 190 µl buff er		*pH5 abts buffer mix for b4						
					Premix ABTS and buffer in future assays			



 [6-3\\_TV\\_assay.xlsx](#)

Plate 6/3



### 6/8/2020-*T. versicolor* Positive Control ABTS Assay

Enzyme dilution prep			
pH 4		pH 5	
2	100 total volume=20microl 20 enzyme(pH4)+80microl buffer (pH4)	2	100 total volume=20microl 20 en
10	100 total volume=50microl 20 enzyme(pH4)+50microl buffer (pH4)	10	100 total volume=50microl 20 en
20	200 total volume=80microl 50 enzyme(pH4)+120microl buffer(pH4)	20	200 total volume=80microl 50 en
50	200 total volume=100microl 100 enzyme(pH4)+100microl buffer (pH4)	50	200 total volume=100microl 100
100	200 total volume=100microl 200 enzyme(pH4)+100microl buffer (pH4)	100	200 total volume=100microl 200
200	200 total volume=20microl stock enzyme (pH7)+180microl buffer (pH4)	200	200 total volume=20microl stock

ABTS Dilution prep			
ABTS	200µl	20µl ABTS	180µlpH4buffer
ABTS	200µl	20µl ABTS	180µlpH5buffer
ABTS	200µl	20µl ABTS	180µlpH7buffer


Calculations for mastermix ABTS/Buffer (need one for each pH)		
32x190 Buffer (each pH)	pH4	pH5
32x5 ABTS (each buffer dilution)	160 microl ABTS buffer + 6080 microl buffer	160 microl ABTS buffer + 6080 microl buffer

	Well	3	4	5	6	7
	Enzyme Concentration	0	2	10	20	50
	pH	Concentrations Enzyme ug per ml				
a	4	5µl ABTS+195µl buffer	5µl 2enzyme µg/ml +5µl ABTS+190µl buffer	5µl 10enzyme µg/ml +5µl ABTS+190µl buffer	5µl 20enzyme µg/ml +5µl ABTS+190µl buffer	5µl 50enzyme µg/ml +5µl ABTS+190µl buffer
b	4	5µl ABTS+195µl buffer	5µl 2enzyme µg/ml +5µl ABTS+190µl buffer	5µl 10enzyme µg/ml +5µl ABTS+190µl buffer	5µl 20enzyme µg/ml +5µl ABTS+190µl buffer	5µl 50enzyme µg/ml +5µl ABTS+190µl buffer
c	control	195µl buffer+5ul buffer	5µl 2enzyme µg/ml +195µl buffer	5µl 10enzyme µg/ml +195µl buffer	5µl 20enzyme µg/ml +195µl buffer	5µl 50enzyme µg/ml +195µl buffer
d	5	5µl ABTS+195µl buffer	5µl 2enzyme µg/ml +5µl ABTS+190µl buffer	5µl 10enzyme µg/ml +5µl ABTS+190µl buffer	5µl 20enzyme µg/ml +5µl ABTS+190µl buffer	5µl 50enzyme µg/ml +5µl ABTS+190µl buffer
e	5	5µl ABTS+195µl buffer	5µl 2enzyme µg/ml +5µl ABTS+190µl buffer	5µl 10enzyme µg/ml +5µl ABTS+190µl buffer	5µl 20enzyme µg/ml +5µl ABTS+190µl buffer	5µl 50enzyme µg/ml +5µl ABTS+190µl buffer
f	control	195ul buffer+5ul buffer	5µl 2enzyme µg/ml +195µl buffer	5µl 10enzyme µg/ml +195µl buffer	5µl 20enzyme µg/ml +195µl buffer	5µl 50enzyme µg/ml +195µl buffer
g	7	5µl ABTS+195µl buffer	5µl 2enzyme µg/ml +5µl ABTS+190µl buffer	5µl 10enzyme µg/ml +5µl ABTS+190µl buffer	5µl 20enzyme µg/ml +5µl ABTS+190µl buffer	5µl 50enzyme µg/ml +5µl ABTS+190µl buffer
h	7	5µl ABTS+195µl buffer	5µl 2enzyme µg/ml +5µl ABTS+190µl buffer	5µl 10enzyme µg/ml +5µl ABTS+190µl buffer	5µl 20enzyme µg/ml +5µl ABTS+190µl buffer	5µl 50enzyme µg/ml +5µl ABTS+190µl buffer



[6-8\\_TV\\_assay.xlsx](#)

6/8 plate



6/15/2020-Preparation of Solutions

Preparing more 4/5/7 buffers			40 ml total	
Prepared	Basic Buffer	0.2M		
	and			
	Citric Acid Buffer	0.1M		
with above buffers we made 3 different pH buffers as follows				
	pH want	0.2M Na <sub>2</sub> HPO <sub>4</sub>	Citric Acid	Actual
	4	15.42ml	24.58ml	8..960
	5	20.6ml	19.4ml	9..100
	7	32.94ml	7.06ml	9..220
*This is the final pH after adjustment using conc. Hydrochloric acid at approximately 32% and 10M NaOH				
Prepared				
BLUE tubes for ABTS/Buf mix	Calculations for mastermix ABTS/Buffer (need one for each pH)			
	32x95 Buffer (each pH)	pH4	pH5	pH7
	32x5 ABTS (each buffer dilution)	160 microl ABTS buffer + 3040 microl buffer	160 microl ABTS buffer + 3040 microl buffer	160 ml
Prepared				
Actual				
	ABTS	0.0548g/0.0548g	=	1.00ml
		0.548g/0.0547g	=	ml UF

6/16/2020-Chlamy ABTS Assay cancelled due to miscommunication

6/18/2020-Chlamy Lysate ABTS Assay

Go from 20-200microg/ml, increase time for test to 2 hours

Enzyme dilution prep		
pH 4		pH 5
20	200 total volume=80microl 50 enzyme(pH5)+120microl buffer(pH4)	20 200 total volume=80m
50	400 total volume=200microl 100 enzyme(pH4)+200microl buffer (pH4)	50 400 total volume=200
100	400 total volume=200microl 200 enzyme(pH4)+200microl buffer (pH4)	100 400 total volume=200
200	400 total volume=40microl stock enzyme (pH7)+360microl buffer (pH4)	200 400 total volume=40m

ABTS Dilution prep			
ABTS	200µl	20µl ABTS	180µlpH4buffer
ABTS	200µl	20µl ABTS	180µlpH5buffer
ABTS	200µl	20µl ABTS	180µlpH7buffer

	Add 100 microl of mastermix to each well, add 100 of each dilution to each well.	
BLUE tubes for ABTS/Buf mix	Calculations for mastermix ABTS buffer/Buffer pH (need one for each pH)	
	32 wellsx95 microl Buffer (each pH)	pH4
	32 wellsx5 microl ABTS (each buffer dilution)	160 microl ABTS buffer + 3040 r

		Repli cate 1					Repli cate 2			
We ll	3	4	5	6	7	8	9	10	11	

	Enzyme Concentration	20	50	100	200	0	20	50	100	200
	pH					Concentrations Enzyme $\mu\text{g}$ per ml				
a	4	5 $\mu\text{l}$ 20enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 50enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 100enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 200enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	100 $\mu\text{l}$ Buffer +5 $\mu\text{l}$ ABT S+95 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 20enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 50enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 100enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer	5 $\mu\text{l}$ 200enzyme $\mu\text{g}$ /ml +5 $\mu\text{l}$ ABT S+19 0 $\mu\text{l}$ buffer







h	7	5µl	5µl	5µl	5µl	5µl	5µl	5µl	5µl	5µl
		20en	50enz	100en	200e	100	5µl	5µl	100e	200e
		zy µg	y µg	zy µg	µg	Buffer	y µg	y µg	µg	µg
		/ml	/ml	/ml	/ml	r	/ml	/ml	/ml	/ml
		+5µl	+5µl	+5µl	+5µl	+5µl	+5µl	+5µl	+5µl	+5µl
		ABT	ABT	ABT	ABT	ABT	ABT	ABT	ABT	ABT
		S+19	S+19	S+19	S+19	S+95	S+19	S+19	S+19	S+19
		0µl	0µl	0µl	0µl	µl	0µl	0µl	0µl	0µl
		buffer	buffer	buffer	buffer	buffer	buffer	buffer	buffer	buffer

b8-b11 mixed in just 5microl abts,  
not abts buffer

e8-e11 mixed in just 5microl abts,  
not abts buffer

h8-h11 mixed in just 5microl abts,  
not abts buffer

 [6-18 Chlamy assay.xlsx](#)

## 6/22/2020-*T. versicolor* Positive Control ABTS Assay

Go from 20-200microg/ml, increase time for test to 4 hours



Step 1				
	ABTS Dilution prep			
	ABTS	200µl	20µl ABTS	180µlpH4buffer
	ABTS	200µl	20µl ABTS	180µlpH5buffer
	ABTS	200µl	20µl ABTS	180µlpH7buffer

OPTIONAL STEP 2	>>>Otherwise just mix ABTS and Buffer in wells separately	
BLACK mix		
Calculations for mastermix ABTS/Buffer (need one for each pH)		
32x190 Buffer (each pH)	pH4	pH5
32x5 ABTS (each buffer dilution)	160 microl ABTS buffer + 6080 microl buffer	160 microl ABTS buffer + 6080 microl buffer

	<p><b>Step 3:</b> Fill in main wells with ABTS /Buffer for each</p>					
	<p><b>Replicate 1</b></p>					<p><b>Replicate 2</b></p>
Well	3	4	5	6	7	8

Enzyme Concentration	20	50	100	200	0	20	
pH						Concentrations Enzyme ug per ml	
4	5µl 20enzyme µg/ml +5µlABTS S+190µl buffer	5µl 50enzyme µg/ml +5µlABTS +190µl buffer	5µl 100enzyme µg/ml +5µlABTS +190µl buffer	5µl 200enzyme µg/ml +5µlABTS S+190µl buffer	5µlABTS +195µl buffer	5µl 20enzyme µg/ml +5µlABTS +190µl buffer	
4	5µl 20enzyme µg/ml +5µlABTS S+190µl buffer	5µl 50enzyme µg/ml +5µlABTS +190µl buffer	5µl 100enzyme µg/ml +5µlABTS +190µl buffer	5µl 200enzyme µg/ml +5µlABTS S+190µl buffer	5µlABTS +195µl buffer	5µl 20enzyme µg/ml +5µlABTS +190µl buffer	

control	5µl 20enzy µg /ml +195µl buffer	5µl 50enzy µg /ml +195µl buffer	5µl 100enzy µg/ml +195µl buffer	5µl 200enzy µg/ml +195µl buffer	195µl buffer+5u l buffer	5µl 20enz /ml +195µl buffer
5	5µl 20enzy µg /ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABTS +190µl buffer	5µl 100enzy µg/ml +5µlABTS +190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20enz /ml +5µlA +190µl buffer
5	5µl 20enzy µg /ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABTS +190µl buffer	5µl 100enzy µg/ml +5µlABTS +190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20enz /ml +5µlA +190µl buffer
control	5µl 20enzy µg /ml +195µl buffer	5µl 50enzy µg /ml +195µl buffer	5µl 100enzy µg/ml +195µl buffer	5µl 200enzy µg/ml +195µl buffer	195ul buffer+5u l buffer	5µl 20enz /ml +195µl buffer

7	5µl 20enzy µg /ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABTS +190µl buffer	5µl 100enzy µg/ml +5µlABTS +190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20enz /ml +5µlA +190µ buffer
7	5µl 20enzy µg /ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABTS +190µl buffer	5µl 100enzy µg/ml +5µlABTS +190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20enz /ml +5µlA +190µ buffer

Step 4: Mix Enzyme dilutions	Enzyme dilution prep		
pH 4			pH 5
20	200 total volume=80microl 50 enzyme(pH5)+120microl buffer(pH4)	20	200 total volume=80m
50	400 total volume=200microl 100 enzyme(pH4)+200microl buffer (pH4)	50	400 total volume=200l
100	400 total volume=200microl 200 enzyme(pH4)+200microl buffer (pH4)	100	400 total volume=200l
200	400 total volume=40microl stock enzyme (pH7)+360microl buffer (pH4)	200	400 total volume=40m

Step 5: In second well plate, add around 10microl of each enzyme dilution, then when assay begins, add enzyme to plate and insert into machine.

[6-22\\_TV\\_assay.xlsx](#)

## 6/23/2020-Prepare new stocks, *T. versicolor* Positive Control ABTS Assay

MONOHYDRATE!						
Citric Acid Buffer	210.14g/mol molecular weight		210.14g/mol x 0.25 L volume		>> 52.535g /250ml	
	0.1 M want	(100ml)	10ml buffer	or for 200ml	20ml buffer	
	1M have		90ml H <sub>2</sub> O		180ml H <sub>2</sub> O	
					created 1M solution	
Prepared	Basic Buffer	0.2M				
	and					
	Citric Acid Buffer	0.1M				
with above buffers we made 3 different pH buffers as follows						
	pH want	0.2M Na <sub>2</sub> HPO <sub>4</sub>	0.1M Citric Acid	Actual pH	Final pH*	
	4	7.71ml	12.29ml	3..90		
	5	10.3ml	9.7ml	4..92		
	7	16.47ml	3.53ml	7..04		
*This is the final pH after adjustment using conc. Hydrochloric acid at approximately 32% and 10M NaOH						
ABTS Stock						
		0.0546g/0. 0547g	=	0.998ml UF H <sub>2</sub> O		

Preparation of Assay

Step 1				
	ABTS Dilution prep			
	ABTS	200µl	20µl ABTS	180µlpH4buffer
	ABTS	200µl	20µl ABTS	180µlpH5buffer
	ABTS	200µl	20µl ABTS	180µlpH7buffer

OPTIONAL STEP 2	>>>Otherwise just mix ABTS and Buffer in wells separately	
BLACK mix		
Calculations for mastermix ABTS/Buffer (need one for each pH)		
32x190 Buffer (each pH)	pH4	pH5
32x5 ABTS (each buffer dilution)	160 microl ABTS buffer + 6080 microl buffer	160 microl ABTS buffer + 6080 microl buffer

		<p><b>Step 3:</b>                  Fill in main wells with ABTS /Buffer for each</p>			<p><b>***ROW SIX SKIPPED FOR ENZYME, ROW 7 USED!!!</b></p>		
		<p><b>Replicate 1</b></p>					<p><b>Replicate 2</b></p>
Well	3	4	5	6	7	8	

	Enzyme Concentration	20	50	100	200	0	20
	pH					Concentrations Enzyme ug per ml	
a	4	5µl 20enzyme µg/ml +5µlABTS S+190µl buffer	5µl 50enzyme µg /ml +5µlABTS S+190µl buffer	5µl 100enzyme µg/ml +5µlABTS S+190µl buffer	5µl 200enzyme µg/ml +5µlABTS S+190µl buffer	5µlABTS +195µl buffer	5µl 20enzyme /ml +5µl S+1 buff
b	4	5µl 20enzyme µg/ml +5µlABTS S+190µl buffer	5µl 50enzyme µg /ml +5µlABTS S+190µl buffer	5µl 100enzyme µg/ml +5µlABTS S+190µl buffer	5µl 200enzyme µg/ml +5µlABTS S+190µl buffer	5µlABTS +195µl buffer	5µl 20enzyme /ml +5µl S+1 buff

c	control	5µl 20enzy µg/ml +195µl buffer	5µl 50enzy µg /ml +195µl buffer	5µl 100enzy µg/ml +195µl buffer	5µl 200enzy µg/ml +195µl buffer	195µl buffer+5 ul buffer	5µl 20e /ml +19 buff
d	5	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABT S+190µl buffer	5µl 100enzy µg/ml +5µlABT S+190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20e /ml +5µ S+1 buff
e	5	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABT S+190µl buffer	5µl 100enzy µg/ml +5µlABT S+190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20e /ml +5µ S+1 buff
f	control	5µl 20enzy µg/ml +195µl buffer	5µl 50enzy µg /ml +195µl buffer	5µl 100enzy µg/ml +195µl buffer	5µl 200enzy µg/ml +195µl buffer	195ul buffer+5 ul buffer	5µl 20e /ml +19 buff




g	7	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABT S+190µl buffer	5µl 100enzy µg/ml +5µlABT S+190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20e /ml +5µ S+1 buff
h	7	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µl 50enzy µg /ml +5µlABT S+190µl buffer	5µl 100enzy µg/ml +5µlABT S+190µl buffer	5µl 200enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl 20e /ml +5µ S+1 buff

Enzyme dilution prep			Step 4: Mix Enzyme dilutions		
pH 4		pH 5		pH 7	Double the amount
20	200 total volume=80microl 50 enzyme(pH5)+120microl buffer (pH4)	20	200 total volume=80microl 50 enzyme(pH5)+120microl buffer (pH4)	20	200 total volume=80microl 50 enzyme(pH5)+120microl buffer (pH4)
50	400 total volume=200microl 100 enzyme(pH4)+200microl buffer (pH4)	50	400 total volume=200microl 100 enzyme(pH4)+200microl buffer (pH4)	50	400 total volume=200microl 100 enzyme(pH4)+200microl buffer (pH4)
100	400 total volume=200microl 200 enzyme(pH4)+200microl buffer (pH4)	100	400 total volume=200microl 200 enzyme(pH4)+200microl buffer (pH4)	100	400 total volume=200microl 200 enzyme(pH4)+200microl buffer (pH4)
200	400 total volume=40microl stock enzyme (pH7)+360microl buffer (pH4)	200	400 total volume=40microl stock enzyme (pH7)+360microl buffer (pH4)	200	400 total volume=40microl stock enzyme (pH7)+360microl buffer (pH4)

Step 5: In second well plate, add around 10microl of each enzyme dilution, then when assay begins, add enzyme to plate and insert into machine.

6-23 plate



### 6/29/2020-*T. versicolor* Positive Control ABTS Assay

ABTS Stock calculation			
	0.0547g/0.0548g	=	0.998ml UF H2O

Go from 2-20mg/ml, increase time for test to 4 hours

Step 1	ABTS Dilution prep			
	ABTS	200µl	20µl ABTS	180µlpH4buffer
	ABTS	200µl	20µl ABTS	180µlpH5buffer
	ABTS	200µl	20µl ABTS	180µlpH7buffer

OPTIONAL STEP 2	>>>Otherwise just mix ABTS and Buffer in wells separately	
BLACK mix		
Calculations for mastermix ABTS/Buffer (need one for each pH)		
32x190 Buffer (each pH)	pH4	pH5
32x5 ABTS (each buffer dilution)	160 microl ABTS buffer + 6080 microl buffer	160 microl ABTS buffer + 6080 microl buffer

		<p><b>Step 3:</b>                  Fill in main wells with ABTS /Buffer for each</p>					
		<p><b>Replicate 1</b></p>					<p><b>Replicate 2</b></p>
Well	3	4	5	6	7	8	

	Enzyme Concentration	2	10	15	20	0	2
	pH					Concentrations Enzyme ug per ml	
a	4	5µl 2enzyme µg/ml +5µlABTS S+190µl buffer	5µl 10enzyme µg/ml +5µlABTS S+190µl buffer	5µl 15enzyme µg/ml +5µlABTS S+190µl buffer	5µl 20enzyme µg/ml +5µlABTS S+190µl buffer	5µlABTS +195µl buffer	5µl µg/ml +5µl S+1 buff
b	4	5µl 2enzyme µg/ml +5µlABTS S+190µl buffer	5µl 10enzyme µg/ml +5µlABTS S+190µl buffer	5µl 15enzyme µg/ml +5µlABTS S+190µl buffer	5µl 20enzyme µg/ml +5µlABTS S+190µl buffer	5µlABTS +195µl buffer	5µl µg/ml +5µl S+1 buff

c	control	5µl 2enzy µg /ml +195µl buffer	5µl 10enzy µg /ml +195µl buffer	5µl 15enzy µg /ml +195µl buffer	5µl 20enzy µg/ml +195µl buffer	195µl buffer+5 ul buffer	5µl µg/ml +19 buff
d	5	5µl 2enzy µg /ml +5µlABT S+190µl buffer	5µl 10enzy µg /ml +5µlABT S+190µl buffer	5µl 15enzy µg /ml +5µlABT S+190µl buffer	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl µg/ml +5µ S+1 buff
e	5	5µl 2enzy µg /ml +5µlABT S+190µl buffer	5µl 10enzy µg /ml +5µlABT S+190µl buffer	5µl 15enzy µg /ml +5µlABT S+190µl buffer	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl µg/ml +5µ S+1 buff
f	control	5µl 2enzy µg /ml +195µl buffer	5µl 10enzy µg /ml +195µl buffer	5µl 15enzy µg /ml +195µl buffer	5µl 20enzy µg/ml +195µl buffer	195ul buffer+5 ul buffer	5µl µg/ml +19 buff


g	7	5µl 2enzy µg /ml +5µlABT S+190µl buffer	5µl 10enzy µg /ml +5µlABT S+190µl buffer	5µl 15enzy µg /ml +5µlABT S+190µl buffer	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl µg/ml +5µ S+1 buff
h	7	5µl 2enzy µg /ml +5µlABT S+190µl buffer	5µl 10enzy µg /ml +5µlABT S+190µl buffer	5µl 15enzy µg /ml +5µlABT S+190µl buffer	5µl 20enzy µg/ml +5µlABT S+190µl buffer	5µlABTS +195µl buffer	5µl µg/ml +5µ S+1 buff

Enzyme dilution prep			Step 4: Mix Enzyme dilu
pH 4			pH 5
2	100 total volume=20microl 20 enzyme(pH4)+80microl buffer(pH4)	2	100 total volume=20mic
10	100 total volume=10microl 100 enzyme(pH4)+90microl buffer(pH4)	10	100 total volume=10mic
15	100 total volume=15microl 100 enzyme(pH4)+85microl buffer(pH4)	15	100 total volume=15mic
20	200 total volume=80microl 50 enzyme(pH4)+120microl buffer(pH4)	20	200 total volume=80mic

Step 5: In second well plate, add around 10microl of each enzyme dilution, then when assay begins, add enzyme to plate and insert into machine.

 [6-28\\_TV\\_assay.xlsx](#)

6-28 plate.png




## 7/1/2020-Drop Assay, Lysate ABTS Assay, Bradford

Frozen cells
pH 5 buffer
Open cells with sonification, centerfuge
Supernatent is used for lysate
More Details in "Lysate Protocoll" Document
Extract Lysate
Perform Bradford
Calculate concentration
Perform Assay
Perform Drop Assay


BSA Stock solution				
2mg/ml				
Dilutions chart				
				For Bradfordassay 20microl of each Lysate dilution needed
BSA-Stock	pH7 Buffer	Conc. ( $\mu\text{g}/\mu\text{l}$ )	$\mu\text{g}$ Protein in Test	Lysate dilutions
0 $\mu\text{L}$	2000 $\mu\text{L}$ for 2ml	0	0	0
50 $\mu\text{L}$	1950 $\mu\text{L}$ for 2ml	0.05	1	0.05
100 $\mu\text{L}$	1900 $\mu\text{L}$ for 2ml	0.1	2	0.1
200 $\mu\text{L}$	1800 $\mu\text{L}$ for 2ml	0.2	4	0.25
400 $\mu\text{L}$	1600 $\mu\text{L}$ for 2ml	0.4	8	0.5
800 $\mu\text{L}$	1200 $\mu\text{l}$ for 2mL	0.8	16	Undiluted



7/1 Bradford1.png



7/1 Bradford2.png



7/1/2020 Drop Test continued:

Drop Test					
Plate set up					
Positive control		UD1	0.5	0.25	0.1
Negative control		UD2	0.5	0.25	0.1
100µL of ABTS distributed on a LB plate					
Drops of 5µL Volume on each point from the corresponding dilution					


Pos control=Just Laccase at 25 dilution

Neg control=pH7 Buffer

0.1-Undiluted

10microl on each spot

Droptest.png



7/1/2020 Assay continued:

BarLac Ecoli lysate assay

Buffer		Actual after mix	
pH7		7.04	
16.47ml	3.53ml		

ABTS Dilution				
	ABTS	200µl	20µl ABTS	180µlpH7buffer

		Lysate 1				L
Well	3	4	5	6	7	
Enzyme Concentration	0.25 UnDi	0.5 UnDi	Undiluted	0	0.	
pH						
a						
b						

c						
d	7	100 µl 0.25 UNDI µg/ml +5µlABTS+ 95µl buffer	100 µl 0.5 UNDI µg/ml +5µlABTS+9 5µl buffer	100 µl UNDI µg/ml +5µlABTS+9 5µl buffer	100 µl Buffer +5µlABTS +95µl buffer	10 50 /n + 95
e	7	100 µl 0.25 UNDI µg/ml +5µlABTS+ 95µl buffer	100 µl 0.5 UNDI µg/ml +5µlABTS+9 5µl buffer	100 µl UNDI µg/ml +5µlABTS+9 5µl buffer	100 µl Buffer +5µlABTS +95µl buffer	10 50 /n + 95
f	control	100µl 0.25 UNDIµg/ml +100µl buffer	100µl 0.5 UNDIµg/ml +100µl buffer	100µl UNDI µg/ml +100µl buffer	200 µl Buffer	10 50 /n bu
g						
h						

 [7-1\\_TV\\_assay.xlsx](#)

7/8/2020-*T. versicolor* Positive Control ABTS Assay

Prepared	Basic Buffer	.2M		
	and			
	Citric Acid Buffer	.1M		
	with above buffers we made 3 different pH buffers as follows			
	pH want	.2M Na <sub>2</sub> HPO <sub>4</sub>	.1M Citric Acid	Actual pH
		4 15.42ml	24.58ml	3.93
		5 20.6ml	19.4ml	4.97
		7 32.94ml	7.06ml	7.03

Step 1	ABTS Dilution prep		
ABTS	200µl	20µl ABTS	180µl pH4 buffer
ABTS	200µl	20µl ABTS	180µl pH5 buffer
ABTS	200µl	20µl ABTS	180µl pH7 buffer

OPTIONAL STEP 2	>>> Otherwise just mix ABTS and Buffer in wells separately	
BLACK mix		
Calculations for mastermix ABTS/Buffer (need one for each pH)		
32x190 Buffer (each pH)	pH4	pH5
32x5 ABTS (each buffer dilution)	160 microl ABTS buffer + 6080 microl buffer	160 microl ABTS buffer + 6080 microl buffer

		Step 3: Fill in main wells with ABTS /Buffer for each							
		pH4			pH5				
Well	1	2	3	4	5	6	7	8	
Enzyme Concentration	R1	R2	R3	Control	R1	R2	R3	C	
Lac case conc					Concentrations Enzyme ug per ml				

a	0	5µlAB TS+19 5µl buffer	5µlAB TS+19 5µl buffer	5µlAB TS+195 µl buffer	200µl buffer	5µlAB TS+19 5µl buffer	5µlAB TS+19 5µl buffer	5µlAB TS+19 5µl buffer	2 b
b	2	5µl 2enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 2enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 2enzy µg/ml +5µlAB TS+190 µl buffer	5µl 2enzy µg/ml +195µl buffer	5µl 2enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 2enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 2enzy µg/ml +5µlA BTS+1 90µl buffer	5 2 µ -
c	10	5µl 10enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 10enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 10enzy µg/ml +5µlAB TS+190 µl buffer	5µl 10enzy µg/ml +195µl buffer	5µl 10enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 10enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 10enzy µg/ml +5µlA BTS+1 90µl buffer	5 1 µ -
d	15	5µl 15enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 15enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 15enzy µg/ml +5µlAB TS+190 µl buffer	5µl 15enzy µg/ml +195µl buffer	5µl 15enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 15enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 15enzy µg/ml +5µlA BTS+1 90µl buffer	5 1 µ -



e	20	5µl 20enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 20enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 20enzy µg/ml +5µlAB TS+190 µl buffer	5µl 20enzy µg/ml +195µl buffer	5µl 20enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 20enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 20enzy µg/ml +5µlA BTS+1 90µl buffer	5 2 µ -
f	50	5µl 50enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 50enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 50enzy µg/ml +5µlAB TS+190 µl buffer	5µl 50enzy µg/ml +195µl buffer	5µl 50enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 50enzy µg/ml +5µlA BTS+1 90µl buffer	5µl 50enzy µg/ml +5µlA BTS+1 90µl buffer	5 5 µ -
g	100	5µl 100enz y µg /ml +5µlA BTS+1 90µl buffer	5µl 100enz y µg /ml +5µlA BTS+1 90µl buffer	5µl 100enz y µg /ml +5µlAB TS+190 µl buffer	5µl 100enz y µg /ml +195µl buffer	5µl 100enz y µg /ml +5µlA BTS+1 90µl buffer	5µl 100enz y µg /ml +5µlA BTS+1 90µl buffer	5µl 100enz y µg /ml +5µlA BTS+1 90µl buffer	5 1 y /

h	200	5µl	5µl	5µl		5µl	5µl	5µl	
		200enz	200enz	200enz		200enz	200enz	200enz	
		y µg	y µg	y µg	5µl	y µg	y µg	y µg	5
		/ml	/ml	/ml	200enz	/ml	/ml	/ml	2
		+5µlA	+5µlA	+5µlAB	y µg	+5µlA	+5µlA	+5µlA	y
		BTS+1	BTS+1	TS+190	/ml	BTS+1	BTS+1	BTS+1	/
		90µl	90µl	µl	+195µl	90µl	90µl	90µl	-
		buffer	buffer	buffer	buffer	buffer	buffer	buffer	t

Enzyme dilution prep			Step 4: Mix Enzyme dilu
pH 4			pH 5
2	100 total volume=20microl 20 enzyme(pH4)+80microl buffer(pH4)	2	100 total volume=20mic
10	100 total volume=50microl 20 enzyme(pH4)+50microl buffer(pH4)	10	100 total volume=50mic
15	100 total volume=15microl 100 enzyme(pH4)+85microl buffer(pH4)	15	100 total volume=15mic
20	200 total volume=80microl 50 enzyme(pH4)+120microl buffer(pH4)	20	200 total volume=80mic
50	200 total volume=100microl 100 enzyme(pH4)+100microl buffer (pH4)	50	200 total volume=100mi
100	200 total volume=100microl 200 enzyme(pH4)+100microl buffer (pH4)	100	200 total volume=100mi
200	200 total volume=20microl stock enzyme (pH7)+180microl buffer (pH4)	200	200 total volume=20mic

 [7-8\\_TV\\_assay.xlsx](#)

### 7/16/2020-*E. coli* marLac and BaLac BL21 ABTS Assay

**Dilutions depend on initial volume. Need at least 300microl of each	
4 hour duration	
20microl as "undiluted" sample for T. vesticolor	

Enzyme dilution prep

2	200 total volume=2microl 200dilution+198microl buffer
200	200 total volume=20microl stock enzyme(pH7)+180microl buffer(pH7)

ABTS Dilution prep			
ABTS	300µl	30µl ABTS	270µlpH4buffer
ABTS	300µl	30µl ABTS	270µlpH7buffer

Calculations for mastermix ABTS/Buffer (need one for each pH)		
40x190 Buffer (each pH)	pH4	pH7
40x5 ABTS (each buffer dilution)	200 microl ABTS buffer + 7600 microl buffer	200 microl ABTS buffer +

	<b>MarLac obtained:</b>		
	65.35 microM for 750microl	<b>Control Enzyme: T. versicolor</b>	<b>Purified Enzyme: MarLac</b>
	Well	5	6
			7

	Enzyme Concentration	0.04mM	0.06535mM	0 enzyme
	pH			
a	7			
b	7	5µl 2enzyme dilution +5µlABTS+190µl buffer	buffer	0 µl enzyme 95µl
c	7	5µl 2enzyme dilution +5µlABTS+190µl buffer	buffer	0 µl enzyme 95µl
d	7	5µl 2enzyme dilution +5µlABTS+190µl buffer	195 µl 1enzyme+5µlABTS+0 µl buffer	0 µl enzyme 95µl
e	7	5µl 2enzyme dilution +5µlABTS+190µl buffer	195 µl 1enzyme+5µlABTS+0 µl buffer	0 µl enzyme 95µl
f	Control 7	5µl 2enzyme dilution +5µlABTS+190µl buffer	195 µl 1enzyme+5µlABTS+0 µl buffer	0 µl enzyme 95µl
g		5µl 2enzyme dilution +0µlABTS+195µl buffer	buffer	0 µl enzyme 95µl

h		Buffer	195 $\mu$ l 1 enzyme+0 $\mu$ l ABTS+5 $\mu$ l buffer	Buffer
			**Buffer will be substituted with incubation buffer (pH7) in this assay	

		Had airbubbles, had to remove some before putting into machine.			p o b
BaLac obtained:		Control Enzyme: T. versicolor	Purified Enzyme: BaLac	Enzyme control	
30microM for 500microl	Well	5	6	7	8
	Enzyme Concentration	0.03mM	0.03mM	0 enzyme	P c
	pH				

a	4				
b	4				
c	4				
d	4				
e	4	5µl Benzy dilution +5µlABTS +190µl buffer	195 µl 1enzyme+ 5µlABTS+ 0µl buffer	0 µl enzyme+5µlA BTS+195µl buffer	
f	Control 4	5µl Benzy dilution +0µlABTS +195µl buffer	195 µl 1enzyme+ 0µlABTS+ 5µl buffer	0 µl enzyme+0µlA BTS+200µl buffer	
g					

h			**Buffer will be substituted with incubation buffer (pH4) in this assay	
---	--	--	---	--


 [7-16 MarLac assay.xlsx](#)

 [7-16 BaLac assay.xlsx](#)

7/16 marLac plate.png



7/16 BaLac plate.png



7/28/2020-*E. coli* marLac ABTS Assay

Salt in process	
Salt added to increase solubility of MarLac	1M stock=584.4mg in 10mL dialysis buffer
NaCl 99.5% p.a., ACS, ISO	Adjusted to 1M stock= 578.4mg in 9.897ml dialysis buffer
M 58.44 g/mol	200microl per well, 100mM salt buffer>>20microl salt buffer in each well
Dichte 2,17	Optimal mM of NaCl according to MarLac paper=100mM
Sdp 1461C	
Carl Roth GmbH+ co	

2microg level chosen for T. ver as concentration should be around 0.04mM for made enzyme, doesn't need to be exact, just acting

4 hour duration		Tversicolor
	2	200 total volume=4microl 200dilution+196microl buffer
	200	200 total volume=20microl stock enzyme(pH7)+180microl buffer(pH7)
		Same dilution will be used for both salt and normal buffer, salt will just be added in place of buffer

ABTS Dilution prep			
ABTS	300µl	30µl ABTS	270µlpH7buffer


MarLac  
obtained:



	28.45 microM for ~800microl	Control Enzyme:T. versicolor	
	Well	5	6
	Enzyme Concentrati on	0.04mM	3.57mM
	pH		
a	7+Salt		
b	7+Salt		
c	7+Salt		
d	7+Salt	5µl 2enzy dilution +5µlABTS+170µl buffer+20µl salt	5µl 200enzy dilution +5µlABTS+170µl b
e	7+Salt- ABTS	5µl 2enzy dilution +0µlABTS+175µl buffer+20µl salt	5µl 200enzy dilution +0µlABTS+175µl b

f	7+Buffer	5 $\mu$ l 2enzy dilution +5 $\mu$ l ABTS+190 $\mu$ l buffer	5 $\mu$ l 200enzy dilution +5 $\mu$ l ABTS+190 $\mu$ l b
g	7+Buffer- ABTS	5 $\mu$ l 2enzy dilution +0 $\mu$ l ABTS+195 $\mu$ l buffer	5 $\mu$ l 200enzy dilution +0 $\mu$ l ABTS+195 $\mu$ l b
h	7-Enzy	0 $\mu$ l enzyme +5 $\mu$ l ABTS+175 $\mu$ l buffer+20 $\mu$ l salt	0 $\mu$ l enzyme+5 $\mu$ l AB buffer

7/28 marLac.png



[7-28 marLac assay.xlsx](#)

7/29/2020-Diclofenac Preparation, PreHPLC Spectra

	Tversicolor dilution
	Needs to be 200microg/ml, combining with buffer/diclo dilution at 1:1 so needs 400microg/ml
400	200microl total volume=40microl stock enzyme(pH7)+160microl buffer (pH4/5/6)



Final probe prep		
	Control	Incubation
pH4	100microl pH4 buffer+ 100microl pH4 diclo	100microl pH4 laccase+ 100microl pH4 diclo
pH5	100microl pH5 buffer+ 100microl pH5 diclo	100microl pH5 laccase+ 100microl pH5 diclo
pH7	100microl pH7 buffer+ 100microl pH7 diclo	100microl pH7 laccase+ 100microl pH7 diclo
	>>Take 100microl of each and measure in spectro after 10 min incubation	
	EACH FINAL WELL WILL HAVE 20MICROG/ML (~0.0625mM?) DICLO AND 200MICROG/ML (3.57mM) LACCASE CONCENTRATIONS	
	***We should calculate diclo's concentration from same we did for ABTS	
Dicl test with HPLC prep to verify it won't be affected by the 0.2microm filter		
Heat to 95C for 10 mins to complete reaction>>centerfuge for 5 mins at 13000		
Filter with needle, 0.2microm pore PVDF? into new epi		
Measure spectra after blanking with buffer ph5		

[pH 5 Day1.xlsx](#)

[pH 5 Day2.xlsx](#)

[pH4&pH7.xlsx](#)

[pH5 HPLC Heat.xlsx](#)

[pH5 HPLC Heat customized.xlsx](#)


## 8/4/2020-*E. coli* BaLac BL21 ABTS Assay

T. versicolor							
2	200 total volume=4microl 200dilution+196microl buffer					<0. 04mM	
200	200 total volume=20microl stock enzyme (pH7)+180microl buffer(dialysis)						

ABTS Dilution prep			
ABTS	200µl	20µl ABTS	180µlpH4buffer

	<b>BaLac obtained:</b>			
	6.37 microM for 500microl	<b>Control Enzyme: T. versicolor</b>	<b>Purified Enzyme:BaLac</b>	<b>Enzyme control</b>
	Well	5	6	7

	Enzyme Concentra tion	0.04mM	0.006mM	0 enzyme
	pH			
a	4			
b	4			
c	4			
d	4			
e	4	5µl TVenzy dilution +5µlABTS+190µl buffer	195 µl 1enzyme+5µlA BTS+0µl buffer	0 µl enzyme+5µ BTS+195µ buffer
f	Control 4	5µl TVenzy dilution +0µlABTS+195µl buffer	195 µl 1enzyme+0µlA BTS+5µl buffer	0 µl enzyme+0µ BTS+200µ buffer
g	4			
h	4			
			**Buffer will be substituted with dialysis buffer (pH4) in this assay	

 [8-4\\_balac\\_assay.xlsx](#)

8/4 BarLac plate.png

