

DNA extraction with AxyPrep Kit

Aim of the experiment

The purpose of this experiment is to extract target DNA.

Materials

- RNase A: 50mg/ml, can be stored at room temperature for 6 months, long-term storage at -20°C
- Buffer S1: suspension, add RNase A and mix well. Store at 4°C
- Buffer S2: Lysis solution (including SDS/NaOH), sealed and stored at room temperature
- Buffer S3: Neutralizing solution, sealed at room temperature
- Buffer W1: washing solution, sealed at room temperature
- Buffer W2 concentrate: Desalt solution. Before use, add absolute ethanol according to the specified volume on the reagent bottle, mix well, and store in airtight at room temperature.
- Eluent: eluent, airtight storage at room temperature

Procedure

1. Before using the kit add all the RNase A into BufferS1 and save it at 4°C.
 2. Suck 4ml microbial into the EP tube and centrifuge it at 12000 rpm for 1 minute, through the supernatant.
 3. Add 250ul Buffer S2 and gently turn it upside and down for 4-6 times.
 4. Add 350ul BufferS3 and mildly turn it upside and down for 6-8 times. Centrifuge it at 12000rpm for 10 minutes.
 5. Suck the supernatant and transfer it into the centrifuge tube. Centrifuge it at 12000rpm for 1minute and through the solution.
-

6. Add 500ul BufferW1, centrifuge it at 12000rpm for 1 minute and through the solution.
 7. Add 700ul Buffer W2, centrifuge it at 12000rpm for 1 minute and through the solution.
 8. Add 700ul Buffer W2 and through the solution.
 9. Centrifuge it at 12000rpm for 1 minute.
 10. Add 60-80 ul Eluent and move the microcentrifuge tube into new 1.5ml centrifuge tube. Let stand at room temperature for 1 minute and centrifuge it at 12000rpm for 1 minute.
-