

Tissue PCR

Introduction

The purpose of this protocol is to confirm correct insertion/deletion of DNA fragments after transformation in *A.niger*.

Info about controls:

Negative control:

Everything but the template (which is the plasmid with bioblocks).

If we have a band then we have contamination.

Positive control:

The genomic DNA for insertion and the insertion plasmid for a knockout.

The bands should be smaller than the ones we get from the bioblocks (the size of a protospacer) because there is only one protospacer/sgRNA in this one.

Materials

- MyTaq plant PCR kit
- Forward Primer (10 μ M)
- Reverse Primer (10 μ M)
- MyTaq DNA Polymerase
- MilliQ water

Procedure

All reactions are set up on ice!!

1. Make a master mix solution of Mytaq plant PCR kit and each Forward and Reverse primer.

Reactant	Per reaction (10 μ L) [μ L]
MyTaq plant PCR kit	5
10 μ M Forward Primer	0.5
10 μ M Reverse Primer	0.5

2. In another PCR tube add 15 μ L MilliQ water. Using an inoculation loop or toothpick poke the fungi from the plate and transfer the loop to the PCR tube with water.
3. Transfer 4 μ L of the solution from PCR with water to the PCR tube with the master mix. Use the rest of the solution for making a stock.
4. Run reaction in a thermocycler using the following settings.

Step	Temperature	Duration	Number of Cycles
Initial denaturation	95 °C	3 minutes	1 cycle
Amplification	95 °C	15 seconds	25-30 cycles
	45-68°C	15 seconds	
	72 °C	around 45 seconds/kb*	
Final extension	72 °C	10 minutes	1 cycle
Hold	12 °C	-	1 cycle

*Not recommended to have more than 3 min extension.

5. The PCR products can then be stored at -20 °C, used directly, or purified using PCR purification or gel extraction.