

Protocol for Preparation of competent cells

Material

Shaker

Ice Maker

Low Temperature Refrigerator

Centrifuge

LB Medium

CaCl₂

Glycerol

Steps

- ① Pick a single colony from the newly activated E. coli plate and inoculate in a 5 ml LB liquid medium.
- ② Put the medium in the shaker and shake with 220rpm at 37 °C for 12h.
- ③ Draw 2ml bacterial liquid transferred to a flask containing 100ml LB liquid medium. Shake it violently at 37 °C for 3-4 h.
- ④ Draw 40ml bacterial liquid to a precooled 1.5 ml centrifuge tube and put in in ice bath for 30min.
- ⑤ Centrifugate the liquid at 4 °C, 4000 rpm for 10 min. Discarded the supernatant.
- ⑥ Add 700µl cold 0.1M CaCl₂ into the tube and put it in ice bath for 30min.
- ⑦ Centrifugate the liquid at 4 °C, 4000 rpm for 10 min. Discarded the supernatant.
- ⑧ Add 300µl cold 0.1M CaCl₂ and 100µl 10%glycerol into the tube and mix them well, then put it in -80 °C refrigerator.

Note

- ① The whole process is operated under sterile condition.