

Protocol for PCR & Agarose Gel Electrophoresis - System A

·Material

Takara Prime STAR HS (Premix)

Template P16-1 plasmid

Template P16-2 plasmid

Template P18 plasmid

Primer

Double Distillation Water (ddH₂O)

Loading Buffer

DNA Marker

·Step

Group 1\2\3\4

- ① Adding 10µl of Takara Prime STAR HS, 1µl of tablet, 4µl of primer, 5µl of ddH₂O into a PCR tube. Different group use specific tube with distinctive sign. Those four group are arranged as table follow:

Group	Primer	Template
1	Seq-P16-F1 Seq-P16-R1	P16-1
2	Seq-P16-F1 Seq-P16-R1	P16-2
3	Seq-P16-F2 Seq-P16-R2	P16-1 or P16-2
4	Seq-P16-F2 Seq-P18-R2	P18

- ② Place those PCR tubes into Peltier thermal cyclor.
- ③ Set the protocol as follow: begin at 98°C for 5 mins, then keep 98°C for 10 secs for denaturation, decrease to 54°C for 15 secs, 72°C for 1 min and repeat that cycle 30 times, finally maintain 4°C infinity.
- ④ Adding 2µl different samples for each, 1µl loading buffer, 7µl ddH₂O into different wells(Adding DNA marker in the first well), then run the gel at 120V for about 30 minutes. Check the result under the UV rays.

Group 5

- ① Adding 10µl of Takara Prime STAR HS, 1µl of tablet, 4µl of primer, 5µl of ddH₂O into a PCR tube. Different group use specific tube with distinctive sign. Those four group are arranged as table follow:

Group	Primer	Template
5	Seq-P16-F1 Seq-P18-R1	P18

- ② Place the PCR tube into Peltier thermal cycler.
- ③ Set the protocol as follow: begin at 98°C for 5 mins, then keep 98°C for 10 secs for denaturation, decrease to 52°C for 15 secs, 72°C for 1 min and repeat that cycle 30 times, finally maintain 4°C infinity.
- ④ Adding 2μl different samples for each, 1μl loading buffer, 7μl ddH₂O into different wells(Adding DNA marker in the first well), then run the gel at 120V for about 30 minutes. Check the result under UV rays.